

(19) World Intellectual Property Organization
International Bureau



(43) International Publication Date
26 July 2001 (26.07.2001)

PCT

(10) International Publication Number
WO 01/53316 A1

(51) International Patent Classification⁷: C07H 21/04, 2800 Warwick Lane, Modesto, CA 95350 (US). IRELAN,
C12Q 1/70 Nancy, A.; 3108 Tuxford Lane, Modesto, CA 95350 (US).

(21) International Application Number: PCT/US01/01458 (74) Agent: SHARP, Jeffrey, S.; Marshall, O'Toole, Gerstein,
Murray & Borun, 6300 Sears Tower, 233 South Wacker
(22) International Filing Date: 17 January 2001 (17.01.2001) Drive, Chicago, IL 60606 (US).

(25) Filing Language: English (81) Designated States (national): AU, BR, CA, NZ, ZA.

(26) Publication Language: English (84) Designated States (regional): European patent (AT, BE,
CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC,
(30) Priority Data: 09/488,295 20 January 2000 (20.01.2000) US NL, PT, SE, TR).

(71) Applicant: E. & J. GALLO WINERY [US/US]; 600
Yosemite Boulevard, P.O. Box 130, Modesto, CA 95353
(US).
Published:
— with international search report
— with amended claims

(72) Inventors: ENGEL, Stacia, R.; 3005 Merle Avenue,
Modesto, CA 95355 (US). MORENZONI, Richard, A.;
For two-letter codes and other abbreviations, refer to the "Guid-
ance Notes on Codes and Abbreviations" appearing at the begin-
ning of each regular issue of the PCT Gazette.

(54) Title: DETECTION OF FERMENTATION-RELATED MICROORGANISMS

(57) Abstract: Unique DNA sequences characteristic of fermentation-related microorganisms are provided which are useful as primers in PCR-based analysis for identification of fermentation-related microorganisms.

WO 01/53316 A1

- 1 -

DETECTION OF FERMENTATION-RELATED MICROORGANISMS

Field of the Invention

The present invention relates to assays to detect fermentation-related
5 microorganisms.

Description of the Related Art

In traditional winemaking, the indigenous yeasts ferment the grape must. Most modern winemakers, however, inoculate with a pure culture of a selected yeast
10 strain to ensure a rapid, reliable and predictable fermentation. It is thought that indigenous yeasts are suppressed by the competitive effect of addition of a high-density monoculture, but some evidence suggests that indigenous yeast can still participate in the fermentation. A range of commercial yeasts with different winemaking characteristics is available, and a number of those yeast strains may be used in a single winery.
15 Furthermore, unwanted microorganisms may be present which lead to spoilage. Therefore, there is a need for a rapid, simple and accurate method for identifying microorganisms in starter cultures and fermentations.

Past techniques for detecting and identifying fermentation-related microorganisms, especially yeast, include colony morphology, fermentation performance,
20 sugar fermentation tests, tolerance to various stresses (e.g., ethanol tolerance) phenotypes with functional relevance (e.g., flocculation) nutritional requirements (e.g., oxygen), and resistance and sensitivity levels of cycloheximide. These methods, however, have numerous disadvantages, including lengthy analysis periods, inability to differentiate, e.g., different strains of yeast, and lack of reproducibility.

Recent developments in molecular biology and protein chemistry have provided new methods for identifying microorganisms, including DNA restriction fragment length polymorphisms, protein electrophoresis patterns and chromosome
25 fingerprinting. Such techniques have been used for identifying fermentation-related microorganisms. See, for example, Casey et al., *Journal of the American Society of Brewing Chemists*, 48(3):100-106, 1990; Degre et al., *American Journal of Enology and Viticulture*, 40(4):309-315, 1989; Guillamon et al., *Systematic and Applied Microbiology*,
30

19:122-132, 1992; Hoeben et al., *Current Genetics*, 10:371-379, 1986; Mozina et al., *Letters in Applied Microbiology*, 24(4):311-315, 1997; Paffetti et al., *Research Microbiology*, 146:587-594, 1995; Panchal et al., *Journal of the Institute of Brewing*, 93:325-327, 1987; Querol et al., *Systematic and Applied Microbiology*, 15:439-446, 1992; Vezinhet et al., *Applied Microbiology and Biotechnology*, 32:568-571, 1990; and Vezinhet et al., *American Journal of Enology and Viticulture*, 43(1):83-86, 1992.

Polymerase chain reaction (PCR)-based techniques have also been used to detect fermentation-related microorganisms. See, for example, DeBarros Lopes et al., *Applied and Environmental Microbiology*, 62(12):4514-4520, 1996; Fell, *Molecular Marine Biology and Biotechnology*, 2(3):174-180, 1993; Fell, *Journal of Industrial Microbiology*, 14(6):475-477, 1995; Ibeas et al., *Applied and Environmental Microbiology*, 62(3):998-1003, 1996; Lavallee et al., *American Journal of Enology and Viticulture*, 45(1):86-91, 1994; Lieckfeldt et al., *Journal of Basic Microbiology*, 33(6):413-425, 1993; and Ness et al., *J. Sci. Food Agric.*, 62:89-94, 1993.

Ribosomal genes are suitable for use as molecular probe targets because of their high copy number. Non-transcribed and transcribed spacer sequences associated with ribosomal genes are usually poorly conserved and, thus, are advantageously used as target sequences for the detection of recent evolutionary divergence. Fungal rRNA genes are organized in units. Each unit encodes mature subunits of 18S, 5.8S, and 28S rRNA. The internal transcribed spacer (ITS) region lies between the 18S and 28S rRNA genes and contains two variable non-coding spacers (referred to as ITS1 and ITS2) and the 5.8S rRNA gene (White et al., 1990; In: *PCR Protocols*; Eds.: Innes et al.; pages 315-322). In addition, the transcriptional units are separated by non-transcribed spacer sequences (NTSs). The ITS and NTS sequences are particularly suitable for the detection of different fungal pathogens.

Kumeda et al. (*Applied and Environmental Microbiology*, 62(8):2947-2952, 1996) describes use of PCR to amplify ribosomal DNA internal transcribed spacers in order to differentiate species of *Aspergillus* Section *Flavi*. The ITS1-5.8S-ITS2 region was amplified using universal primers, and the PCR product analyzed by the principle of single-strand conformation polymorphism. In addition, Gardes et al. (In: *Methods in Molecular Biology, Vol. 50:Species Diagnostics Protocols*:

- 3 -

PCR and Other Nucleic Acid Methods, Ed. J.P. Clapp, Humana Press, Totowa, NJ, (1996) pp. 177-186) describes restriction fragment length polymorphism (RFLP) analysis of fungal ITS regions amplified by PCR.

5 The PCR amplification of fungal ITS has also been described using other than universal primers. These methods allow for more specificity in identifying classes of fungi, or particular species of fungi. Thus, Gardes and Bruns (*Molecular Ecology*, 2:113-118, 1993) identified ITS primers which allow differentiation of DNA from basidiomycetes against ascomycete DNA. Identification of specific species has been observed using PCR primers directed to unique sequences in the ITS1 and/or ITS2
10 regions of fungal pathogens. See, for example, Hamelin et al., *Applied and Environmental Microbiology*, 62(11):4026-4031, 1996; Mazzola et al., *Phytopathology*, 86(4):354-360, 1996; O'Gorman et al., *Canadian Journal of Botany*, 72:342-346, 1994; and U.S. Patent No. 5,585,238 to Ligon et al.

Of interest to the present application is the disclosure of PCT International
15 Application US99/04251 based on U.S. Application Serial No. 09/037,990 which relates to oligonucleotide primers for the ITS of ribosomal RNA gene regions of fermentation-related microorganisms such as *Botrytis cinerea*, *Penicillium*, *Brettanomyces/Dekkera*, *Saccharomyces*, *Hanseniaspora/Kloeckera*, *Candida krusei/Issatchenkia orientalis*, *Pichia kluyveri*, *Pichia anomala/Candida pelliculosa*, *Debaryomyces carsonii*, and
20 *Saccharomyces ludwigii*. Also of interest to the present application is the disclosure of PCT International Application US98/25219 based on U.S. Application Serial No. 08/986,727 which relates to oligonucleotide primers for the ITS of ribosomal RNA gene regions of fungal pathogens such as those that infect grape plants.

The present invention addresses the problem of detecting and identifying
25 fermentation-related microorganisms by PCR-based techniques.

Summary of the Invention

The present invention is directed to the identification of different fermentation-related microorganisms, particularly those involved in the production of
30 wine. The present invention provides DNA sequences which exhibit variability between different fermentation-related microorganisms. In particular, the present invention

identifies regions of DNA sequence located in the internal transcribed spacer (ITS) of the ribosomal RNA gene regions of various fermentation-related microorganisms. Primers derived from the ITS can be used in polymerase chain reaction (PCR)-based and other diagnostic assays to determine the presence or absence of specific fermentation-related microorganisms, including those involved in the production of wine. The primers can also be used as molecular probes to detect the presence of target DNA.

Thus, in one aspect, the present invention provides an isolated double stranded nucleic acid of the full length ITS1 or ITS2 region of a fermentation-related microorganism. More particularly, the DNA sequence is selected from among Sequence ID NOS: 6 to 15 and their complementary sequences.

In another aspect, the present invention provides an oligonucleotide primer for identification of a fermentation-related microorganism, wherein the primer is a divergent portion of the ITS1 or ITS2 region of a fermentation-related microorganism. More particularly, the oligonucleotide primer is selected from among Sequence ID NOS: 16 to 24. Furthermore, the oligonucleotide primers may be selected from among sequences which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24, primers which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15, and primers of 10 nucleotide bases or longer which contain at least 5 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with from 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15. A pair of the foregoing oligonucleotide primers for use in the amplification-based detection of an ITS of a fermentation-related microorganism is also provided.

In yet another aspect, a method is provided for detection of a fermentation-related microorganism which comprises: (a) obtaining DNA from a fungal culture or colony isolated from a fermentation, or from an organism present in a fermentation beverage; (b) amplifying a part of the ITS of the fermentation-related microorganism using the DNA as a template in a polymerase chain reaction with the aforementioned oligonucleotide primers; and (c) detecting the amplified part of the ITS sequence to determine whether the fermentation-related microorganism is present.

In still another aspect, kits are provided which are useful in detecting fermentation-related microorganisms comprising one or more containers, at least one of which comprise an oligonucleotide primer or a pair of oligonucleotide primers according to the invention.

5

Brief Description of the Sequences in the Sequence Listing:

	SEQ ID NO: 1	DNA sequence for the internal transcribed spacer of <i>Metschnikowia pulcherrima</i> .
	SEQ ID NO: 2	DNA sequence for the internal transcribed spacer of <i>Candida stellata</i> .
10	SEQ ID NO: 3	DNA sequence for the internal transcribed spacer of <i>Zygosaccharomyces bailii</i> .
	SEQ ID NO: 4	DNA sequence for the internal transcribed spacer of <i>Kluyveromyces thermotolerans</i> .
15	SEQ ID NO: 5	DNA sequence for the internal transcribed spacer of <i>Torulaspora delbrueckii</i> .
	SEQ ID NO: 6	DNA sequence for the ITS1 of <i>Metschnikowia pulcherrima</i> .
	SEQ ID NO: 7	DNA sequence for the ITS2 of <i>Metschnikowia pulcherrima</i> .
	SEQ ID NO: 8	DNA sequence for the ITS1 of <i>Candida stellata</i> .
20	SEQ ID NO: 9	DNA sequence for the ITS2 of <i>Candida stellata</i> .
	SEQ ID NO: 10	DNA sequence for the ITS1 of <i>Zygosaccharomyces bailii</i> .
	SEQ ID NO: 11	DNA sequence for the ITS2 of <i>Zygosaccharomyces bailii</i> .
	SEQ ID NO: 12	DNA sequence for the ITS1 of <i>Kluyveromyces thermotolerans</i> .
	SEQ ID NO: 13	DNA sequence for the ITS2 of <i>Kluyveromyces thermotolerans</i> .
25	SEQ ID NO: 14	DNA sequence for the ITS1 of <i>Torulaspora delbrueckii</i> .
	SEQ ID NO: 15	DNA sequence for the ITS2 of <i>Torulaspora delbrueckii</i> .
	SEQ ID NO: 16	Oligonucleotide Sequence MXL258.
	SEQ ID NO: 17	Oligonucleotide Sequence MXL268.
	SEQ ID NO: 18	Oligonucleotide Sequence CsF.
30	SEQ ID NO: 19	Oligonucleotide Sequence CsR.
	SEQ ID NO: 20	Oligonucleotide Sequence QAZ165.

- 6 -

SEQ ID NO: 21	Oligonucleotide Sequence QAZ622.
SEQ ID NO: 22	Oligonucleotide Sequence GSZ450
SEQ ID NO: 23	Oligonucleotide Sequence TdF1
SEQ ID NO: 24	Oligonucleotide Sequence TdR1
5 SEQ ID NO: 25	Oligonucleotide Sequence ITS5.
SEQ ID NO: 26	Oligonucleotide Sequence ITS4.

Detailed Description Of The Invention

10 The present invention provides unique DNA sequences which are useful in identifying fermentation-related microorganisms. These unique DNA sequences can be used as primers in PCR-based analysis for the identification of fermentation-related microorganisms, or as molecular probes to detect the presence of DNA from fermentation-related microorganisms. The DNA sequences of the present invention include the internal transcribed spacer (ITS) of the ribosomal RNA gene regions of
15 specific fermentation-related microorganisms, as well as primers that are derived from these regions which are capable of identifying the particular microorganism.

The DNA sequences of the invention are from the ITS of the ribosomal RNA gene region of fermentation-related microorganisms. However, the present invention is not limited to detecting the presence of the microorganisms in fermentation
20 operations, i.e., the invention can be used to detect the presence of such microorganisms from any source. There is variability in the ITS DNA sequences from different microorganisms. The ITS sequences can be aligned and compared. Primers can be designed based on regions within the ITS regions that contain the greatest differences in sequence among the fermentation-related microorganisms. The sequences and primers
25 based on these sequences can be used to identify specific microorganisms.

DNA sequences of particular interest include ITS DNA sequences from *Metschnikowia* sp., especially *Metschnikowia pulcherrima* (anamorph *Candida pulcherrima*); *Zygosaccharomyces* sp., especially *Zygosaccharomyces bailii*; *Kluyveromyces* sp., especially *Kluyveromyces thermotolerans*; *Candida* sp., especially
30 *Candida stellata* and *Torulaspora* sp., especially *Torulaspora delbrueckii* (anamorph *Candida colliculosa*). The ITS DNA sequences, as well as primers of interest, are set

forth in SEQUENCE ID NOS: 1-24. The sequences are useful in PCR-based identification of fermentation-related microorganisms.

Methods for use of the primer sequences of the invention in PCR analysis are well known in the art. For example, see U.S. Patent Nos. 4,683,195; 4,683,202 and
5 5,585,238, the contents of all of which are hereby incorporated by reference.

The primer sequences of the invention can also be used as molecular probes to detect the presence of target DNA. The T_m for the primers ranges from about 48-58 °C at 50 mM salt. The hybridization temperature is approximately 5-10 °C below the melting temperature. Thus, the primers are hybridized to target DNA typically at a
10 temperature ranging from about 43-55°C. Final wash conditions generally range from about 45-55°C at about 36 mM salt concentration. Specific hybridization as used herein means the use of a final high stringency wash in about 0.2X SSPE (salt concentration of about 36 mM) at a temperature appropriate for the particular primer. IX SSPE contains 10 mM NaH_2PO_4 , 180 mM NaCl, and 1 mM EDTA, at pH 7.4.

15 The ITS DNA sequences of the present invention can be cloned from fermentation-related microorganisms by methods known in the art. In general, the methods for the isolation of DNA from microorganism isolates are known. See, Raeder et al., *Letters in Applied Microbiology*, 2:17-20, 1985; Lee et al., *Fungal Genetics Newsletter*, 35:23-24, 1990; and Lee et al., In: *PCR Protocols: A Guide to Methods and Applications*, Innes et al. (Eds.); pages 282-287, 1990; the contents of all of which are
20 hereby incorporated by reference.

Alternatively, the ITS regions of interest can be identified by PCR amplification. Primers to amplify the entire ITS region can be synthesized according to White et al. (1990; In *PCR Protocols*; Eds.: Innes et al., pages 315-322, the contents of
25 which are hereby incorporated by reference).

The ITS sequences were determined and the sequences were compared to locate divergences which might be useful to test in PCR to distinguish the different fermentation-related microorganisms. The sequences of the ITS regions which were determined are shown as Sequence ID NOS: 1 to 5. The DNA sequences for the ITS1
30 and ITS2 regions are shown as Sequence ID NOS: 6 to 15. From the identification of divergences, numerous primers were synthesized and tested in PCR-amplification.

Purified microorganism DNA and DNA isolated from fermentation cultures and colonies were used as templates for PCR-amplification. Thus, pairs of diagnostic primers were identified, i.e., those which identified one particular fermentation-related microorganism species. Preferred primer combinations are able to distinguish between the different microorganisms in, for example, fermentation cultures. Primer sequences are set forth in Sequence ID NOS: 16 to 24. Thus, while oligonucleotide primers selected from among Sequence ID NOS: 16 to 24 are preferred, primers may also be used which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24. Additionally, primers may be used which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of corresponding SEQ ID NOS: 6 to 15, and primers of 10 nucleotide bases or longer which contain at least 5 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with from 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of corresponding SEQ ID NOS: 6 to 15.

The present invention provides numerous diagnostic primer combinations. The primers of the invention are designed based on sequence differences among the microorganism ITS regions. A minimum of one base pair difference between sequences can permit design of a discriminatory primer. In general, primers should have a theoretical melting temperature between about 55 °C to about 65 °C to achieve good sensitivity, and should be void of significant secondary structure and 3' overlaps between primer combinations. Primers are generally at least about 10 nucleotide bases, more preferably at least about 15 to about 20 nucleotide bases.

The oligonucleotide primers of the present invention are particularly useful in detecting microorganisms involved in fermentations, in particular, microorganisms selected from among *Metschnikowia* sp., especially *Metschnikowia pulcherrima* (anamorph *Candida pulcherrima*); *Zygosaccharomyces* sp., especially *Zygosaccharomyces bailii*; *Kluyveromyces* sp., especially *Kluyveromyces thermotolerans*; *Candida* sp., especially *Candida stellata* and *Torulaspora* sp., especially *Torulaspora delbrueckii* (anamorph *Candida colliculosa*). However, the primers of the present invention can also be used to detect the presence of the foregoing microorganisms from any source.

- 9 -

The present invention also relates to the preparation of "kits" containing elements for detecting fermentation-related microorganisms. Such a kit may comprise a carrier to receive therein one or more containers, such as tubes or vials. Unlabeled or detectably labeled oligonucleotide primers may be contained in one or more of the
5 containers. The oligonucleotide primers may be present in lyophilized form, or in an appropriate buffer. One or more enzymes or reagents for use in PCR reactions may be contained in one or more of the containers. The enzymes or reagents may be present alone or in admixture, and in lyophilized form or in appropriate buffers. The kit may also contain any other component necessary for carrying out the present invention, such as
10 buffers, extraction agents, enzymes, pipettes, plates, nucleic acids, nucleoside triphosphates, filter paper, gel materials, transfer materials, and autoradiography supplies.

The examples below illustrate typical experimental protocols which can be used in the isolation of ITS sequences, the selection of suitable primer sequences, the testing of primers for selective and diagnostic efficacy, and the use of such primers to
15 detect the presence of a fermentation-related microorganism. Such examples are provided by way of illustration and not by way of limitation.

Example 1

Culture of Yeast and Fungal Isolates

20

Viable isolates of *Candida stellata*, *Kluyveromyces thermotolerans*, *Metschnikowia pulcherrima*, *Torulaspora delbrueckii*, *Zygosaccharomyces bailii*, *Brettanomyces bruxellensis*, *Candida parapsilosis*, *Candida tropicalis*, *Hanseniaspora uvarum*, *Pichia anomala*, *Saccharomyces cerevisiae*, *Zygosaccharomyces bisporus*, and
25 *Zygosaccharomyces cidri* were obtained from the E & J Gallo Winery, the Gallo of Sonoma Winery, the American Type Culture Collection (ATCC) or the Centraalbureau voor Schimmelcultures, The Netherlands (see Table 1). Yeasts were grown on any of several media of choice.

Example 2

Amplification and Sequencing of the Internal Transcribed Spacer (ITS) Regions

The internal transcribed spacer region was amplified from the different isolates directly from the yeast colony using ITS5 (5'-GGAAGTAAAAGTCGTAACAAGG-3'; SEQ ID NO: 25) and ITS4 (5'-TCCTCCGCTTATTGATATGC-3'; SEQ ID NO: 26). A sterile pipette tip was used to scrape a small amount of colony off of the plate and deposited into a 200- μ l microcentrifuge tube containing 5 μ l each of GeneAmp[®] 10X PCR Buffer (PE Applied Biosystems, Foster City, CA; part no. N808-0160), 0.2 mM each of dATP, dCTP, dGTP, and dTTP (GeneAmp[®] dNTPs; PE Applied Biosystems, Foster City, CA; part no. N808-0007), approximately 25 pmole/ μ l each of ITS5 and ITS4, and 2.5 Units AmpliTaq[®] DNA polymerase (PE Applied Biosystems; part no. N808-0160). Reactions were run for 35 cycles of 30 s at 94°C, 40 s at 58°C, and 2 min at 72°C, followed by a final elongation step at 72°C for 10 min, on a Perkin Elmer GeneAmp[®] PCR System 9700 (PE Applied Biosystems). PCR products were purified using QIAquick[®] PCR Purification Kits (Qiagen Inc., Santa Clarita, CA) to remove any excess primers, nucleotides, and polymerases. Five microliters of the purified PCR products were run on a 1.2% agarose gel with 5 μ l of pGEM-3Zf(+) double-stranded DNA Control Template (0.2 g/L, PE Applied Biosystems) to approximate concentrations. All products were sequenced using the primers ITS5 and ITS4 (see sequences above; White et al., 1990; In: *PCR Protocols*; Eds.: Innes et al. pp. 315-322). Sequencing was performed on an PE Applied Biosystems 377 Automated DNA Sequencer[®] using ABI PRISM[™] dRhodamine Terminator Cycle Sequencing Ready Reaction Kits[®] (PE Applied Biosystems; part no. 403044). Cycle sequencing products were run over Centri-Sep[®] spin columns (Princeton Separations, Inc., Adelphia, NJ) to remove excess primers, dye-labeled terminators, nucleotides, and polymerases before being run on the automated sequencer.

- 11 -

Example 3

Selection of Species-Specific Primers

The ITS sequences of the *Metschnikowia* sp., especially *Metschnikowia pulcherrima* (anamorph *Candida pulcherrima*); *Zygosaccharomyces* sp., especially *Aygosaccharomyces bailii*; *Kluyveromyces* sp., especially *Kluyveromyces thermotolerans*; *Candida* sp., especially *Candida stellata* and *Torulaspora* sp., especially *Torulaspora delbrueckii* (anamorph *Candida colliculosa*) isolates were aligned and primers were designed using Oligo 5.0 (National Biosciences, Inc., Plymouth, MN) in regions of maximum sequence difference between the target species. (See Table 2)

Example 4

Primer Synthesis

Primers were synthesized on an Applied Biosystems 394 DNA/RNA Synthesizer[®] using phosphoramidite chemistry.

Example 5

Verification of Primer Specificity to Target Species

Different annealing temperatures were tested to determine the optimal temperature for PCR for individual primers. In cases with multiple species-specific primers, different primer combinations were used to determine the best primer combination and annealing temperature to amplify a single species-specific DNA fragment. Species-specific amplification products were produced from primers designed from the ITS region between the 18S and 28S ribosomal DNA subunits of each yeast species of interest.

Species-specific primers were tested against other yeast species to confirm their failure to amplify targets in those species. Specifically, the *Metschnikowia pulcherrima*-specific primers (SEQ ID NOs: 16 and 17) were tested in various combinations with and without primers ITS5 (SEQ ID NO: 25) and ITS4 (SEQ ID NO:

26) against the following species (*Candida stellata*, *Kluyveromyces thermotolerans*,
Torulaspora delbrueckii, *Zygosaccharomyces bailii*, *Brettanomyces bruxellensis*,
Candida parapsilosis, *Candida tropicalis*, *Hanseniaspora uvarum*, *Pichia anomala*,
Saccharomyces cerevisiae, *Zygosaccharomyces bisporus*, and *Zygosaccharomyces cidri*)
5 and did not amplify targets in those species.

Further, the *Candida stellata*-specific primers (SEQ ID NOs: 18 and 19)
were tested in various combinations with and without primers ITS5 (SEQ ID NO: 25) and
ITS4 (SEQ ID NO: 26) against the following species (*Kluyveromyces thermotolerans*,
Metschnikowia pulcherrima, *Torulaspora delbrueckii*, *Zygosaccharomyces bailii*,
10 *Brettanomyces bruxellensis*, *Candida parapsilosis*, *Candida tropicalis*, *Hanseniaspora*
uvarum, *Pichia anomala*, *Saccharomyces cerevisiae*, *Zygosaccharomyces bisporus*, and
Zygosaccharomyces cidri) and did not amplify targets in those species.

Further, the *Zygosaccharomyces bailii*-specific primers (SEQ ID NOs: 20
and 21) were tested in various combinations with and without primers ITS4 (SEQ ID NO:
15 26) and ITS5 (SEQ ID NO: 25) against the following species (*Candida stellata*,
Kluyveromyces thermotolerans, *Metschnikowia pulcherrima*, *Torulaspora delbrueckii*,
Brettanomyces bruxellensis, *Candida parapsilosis*, *Candida tropicalis*, *Hanseniaspora*
uvarum, *Pichia anomala*, *Saccharomyces cerevisiae*, *Zygosaccharomyces bisporus*, and
Zygosaccharomyces cidri) and did not amplify targets in those species.

Further, the *Kluyveromyces thermotolerans*-specific primer (SEQ ID
20 NO:22) was tested with ITS5 (SEQ ID NO: 26) against the following species (*Candida*
stellata, *Metschnikowia pulcherrima*, *Torulaspora delbrueckii*, *Zygosaccharomyces*
bailii, *Brettanomyces bruxellensis*, *Candida parapsilosis*, *Candida tropicalis*,
Hanseniaspora uvarum, *Pichia anomala*, *Saccharomyces cerevisiae*, *Zygosaccharomyces*
25 *bisporus*, and *Zygosaccharomyces cidri*) and did not amplify targets in those species.

Finally, the *Torulaspora delbrueckii*-specific primers (SEQ ID NOs: 23
and 24) were tested in various combinations with and without primers ITS4 (SEQ ID NO:
26) and ITS5 (SEQ ID NO: 24) against the following species (*Candida stellata*,
Kluyveromyces thermotolerans, *Metschnikowia pulcherrima*, *Zygosaccharomyces bailii*,
30 *Brettanomyces bruxellensis*, *Candida parapsilosis*, *Candida tropicalis*, *Hanseniaspora*

- 13 -

uvarum, *Pichia anomala*, *Saccharomyces cerevisiae*, *Zygosaccharomyces bisporus*, and *Zygosaccharomyces cidri*) and did not amplify targets in those species.

Example 6

5 Utilization of ITS sequences as diagnostic probes to hybridize with target DNA

1. Put chosen concentration of DNA sample in 100 ul of TE, pH 7.0.
2. Add 0.1 volume of 3.0 M NaOH, vortex to mix and incubate at 65°C for 20 min to destroy the RNA and denature the DNA.
- 10 3. Spin down condensation. Allow samples to cool to room temp. Neutralize by adding 1.0 volume [110 μ l] of 2M ammonium acetate, pH 7.0, vortex to mix. Spin down to remove solution off of cap. Refrigerate until slot blot apparatus is ready.
- 15 4. Apply to 220 μ l slot-blot apparatus according to manufacturer's protocol.
5. Label ITS sequence probe according to kit manufacturer's recommendation.
- 20 6. Prehybridize blots in 1.0% BSA, 1mM EDTA, 0.5 M NaHPO₄, pH 7.2, 7.0% sodium dodecyl sulfate for a minimum of 2 hr prior to adding the probe, and then hybridize for 16 hr at 45°C. Initial washes consist of two 30-min washes in 1X SSPE/0.1% SDS at 50°C. Transfer blots to a plastic tray and wash in 1X SSPE for 1 hr, at 50°C with shaking. The final wash should consist of 15 min at 50°C in 0.2X SSPE.

25 While the invention has been described in detail and with reference to specific embodiments thereof, it will be apparent to one skilled in the art that various changes and modifications can be made therein without departing from the spirit and scope thereof.

Table 1		
Sources of test isolates		
5	Species name	ID number
	Source	
	<i>Candida stellata</i>	CBS157
		CBS ¹
	<i>Candida stellata</i>	CBS2649
		CBS
	<i>Candida stellata</i>	GS056
		Gallo of Sonoma ²
	<i>Candida stellata</i>	GS111
		Gallo of Sonoma
10	<i>Candida stellata</i>	GS128
		Gallo of Sonoma
	<i>Candida stellata</i>	GS174
		Gallo of Sonoma
	<i>Kluyveromyces thermotolerans</i>	GS003
		Gallo of Sonoma
	<i>Kluyveromyces thermotolerans</i>	GS104
		Gallo of Sonoma
	<i>Kluyveromyces thermotolerans</i>	GS119
		Gallo of Sonoma
15	<i>Kluyveromyces thermotolerans</i>	GS127
		Gallo of Sonoma
	<i>Kluyveromyces thermotolerans</i>	GS132
		Gallo of Sonoma
	<i>Kluyveromyces thermotolerans</i>	GS137
		Gallo of Sonoma
	<i>Kluyveromyces thermotolerans</i>	GS138
		Gallo of Sonoma
	<i>Kluyveromyces thermotolerans</i>	GS166
		Gallo of Sonoma
20	<i>Metschnikowia pulcherrima</i>	GS002
		Gallo of Sonoma
	<i>Metschnikowia pulcherrima</i>	GS011
		Gallo of Sonoma
	<i>Metschnikowia pulcherrima</i>	GS018
		Gallo of Sonoma
	<i>Metschnikowia pulcherrima</i>	GS024
		Gallo of Sonoma
	<i>Metschnikowia pulcherrima</i>	GS030
		Gallo of Sonoma
25	<i>Metschnikowia pulcherrima</i>	GS043
		Gallo of Sonoma
	<i>Metschnikowia pulcherrima</i>	GS055
		Gallo of Sonoma
	<i>Metschnikowia pulcherrima</i>	GS066
		Gallo of Sonoma
	<i>Torulaspora delbrueckii</i>	GS038
		Gallo of Sonoma
	<i>Torulaspora delbrueckii</i>	NS1-3
		Gallo of Sonoma

Table 1		
Sources of test isolates		
Species name	ID number	Source
<i>Torulaspora delbrueckii</i>	NS1-5	Gallo of Sonoma
<i>Torulaspora delbrueckii</i>	NS1-9	Gallo of Sonoma
<i>Torulaspora delbrueckii</i>	NS1-11	Gallo of Sonoma
<i>Torulaspora delbrueckii</i>	NS1-16	Gallo of Sonoma
<i>Torulaspora delbrueckii</i>	NS1-19	Gallo of Sonoma
<i>Zygosaccharomyces bailii</i>	QA17	E & J Gallo Winery ³
<i>Zygosaccharomyces bailii</i>	QA30	E & J Gallo Winery
<i>Zygosaccharomyces bailii</i>	QA31	E & J Gallo Winery
<i>Zygosaccharomyces bailii</i>	QA48	E & J Gallo Winery
<i>Zygosaccharomyces bailii</i>	QA57	E & J Gallo Winery
<i>Brettanomyces bruxellensis</i>	Y153	E & J Gallo Winery
<i>Candida parapsilosis</i>	QA45	E & J Gallo Winery
<i>Candida tropicalis</i>	QA44	E & J Gallo Winery
<i>Hanseniaspora uvarum</i>	GS014	Gallo of Sonoma
<i>Pinchia anomala</i>	34080	ATCC ⁴
<i>Saccharomyces cerevisiae</i>	GS084	Gallo of Sonoma
<i>Zygosaccharomyces bisporus</i>	Y476	Gallo of Sonoma
<i>Zygosaccharomyces cidri</i>	36238	ATCC

¹ Centraalbureau voor Schimmelcultures, The Netherlands

² Gallo of Sonoma Winery, Healdsburg, CA, USA

³ E & J Gallo Winery, Modesto, CA, USA

⁴ American Type Culture Collection, Rockville, MD, USA

Table 2		
Target Organism	Primer Name	Primer Sequence
5 <i>Metschnikowia pulcherrima</i>	MXL 258	5'-AAGCAGGACCAAACCGGAGG-3' (SEQ ID NO: 16)
<i>Metschnikowia pulcherrima</i>	MXL 268	5'-TATTAGGCCGAAGCAGGACC-3' (SEQ ID NO: 17)
10 <i>Candida stellata</i>	CsF	5'-TTTGCCAAAACCACTGTGAACA-3' (SEQ ID NO: 18)
<i>Candida stellata</i>	CsR	5'-TTTAAAGATTGGGCGCCTTTC-3' (SEQ ID NO: 19)
<i>Zygosaccharomyces bailii</i>	QAZ165	5'-TGGGAGGATGGGTTCGTC-3' (SEQ ID NO: 20)
15 <i>Zygosaccharomyces bailii</i>	QAZ 622	5'-GCTATCACTCACCCAATCTC-3' (SEQ ID NO: 21)
<i>Kluyveromyces thermotolerans</i>	QSZ 450	5'-CCTCAGTCAGCAACAGCC-3' (SEQ ID NO: 22)
<i>Torulaspora delbrueckii</i>	TdF1	5'-CTATATGAATGAAGTTAGAGGACGTCTAAAGAT-3' (SEQ ID NO: 23)
20 <i>Torulaspora delbrueckii</i>	TdR1	5'-GGAAGCACGCACAAGACGTATC-3' (SEQ ID NO: 24)

We claim:

1. An isolated double stranded nucleic acid selected from the group consisting of SEQ ID NOS: 6 to 15 and their complementary sequences.
5
2. An isolated nucleic acid which specifically hybridizes with the nucleic acid of claim 1.
3. An oligonucleotide sequence for identification of a
10 fermentation-related microorganism, wherein said sequence is selected from the group consisting of SEQ ID NOS: 16 to 24.
4. An oligonucleotide primer which is a fragment of the sequences according to claim 3, and which specifically hybridizes to the ITS1 or ITS2 of
15 *Metschnikowia pulcherrima*; *Zygosaccharomyces bailii*; *Kluyveromyces thermotolerans*; *Candida stellata* and *Torulaspora delbrueckii*.
5. An oligonucleotide primer for identification of a fermentation-related microorganism, wherein said primer is selected from the group
20 consisting of primers which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24, primers which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15, and primers of 10 nucleotide bases or
25 longer which contain at least 5 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with from 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15.

- 18 -

6. A pair of oligonucleotide primers for use in the amplification-based detection of an internal transcribed spacer sequence of a fermentation-related microorganism, wherein the primers are selected from the group consisting of primers which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24, primers which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15, and primers of 10 nucleotide bases or longer which contain at least 5 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with from 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15.

7. The pair of oligonucleotide primers according to claim 6 which comprises SEQ ID NO: 16 and SEQ ID NO: 17.

8. The pair of oligonucleotide primers according to claim 6 which comprises SEQ ID NO: 18 and SEQ ID NO: 19.

9. The pair of oligonucleotide primers according to claim 6 which comprises SEQ ID NO: 20 and SEQ ID NO: 21.

10. The pair of oligonucleotide primers according to claim 6 which comprises SEQ ID NO: 23 and SEQ ID NO: 24.

11. A method for detection of a fermentation-related microorganism comprising:

- (a) obtaining DNA from a fungal culture or colony isolated from a fermentation, or from an organism present in a fermentation beverage;
- (b) amplifying a part of the internal transcribed spacer sequence of said fermentation-related microorganism using said DNA as a template in a polymerase chain reaction with a pair of oligonucleotide primers according to claim 6; and

- 19 -

(c) detecting said amplified part of the internal transcribed spacer sequence to determine whether said fermentation-related microorganism is present.

12. The method according to claim 11, wherein said
5 fermentation-related microorganism is selected from the group consisting of *Metschnikowia pulcherrima*; *Zygosaccharomyces bailii*; *Kluyveromyces thermotolerans*; *Candida stellata* and *Torulaspora delbrueckii*.

13. The method according to claim 12, wherein said fermentation
10 culture or fermentation beverage is a wine fermentation culture or wine fermentation beverage.

14. The method according to claim 12, wherein the pair of
15 oligonucleotide primers comprises SEQ ID NO: 16 and SEQ ID NO: 17.

15. The method according to claim 12, wherein the pair of
oligonucleotide primers comprises SEQ ID NO: 18 and SEQ ID NO: 19.

16. The method according to claim 12, wherein the pair of
20 oligonucleotide primers comprises SEQ ID NO: 20 and SEQ ID NO: 21.

17. The method according to claim 12, wherein the pair of
oligonucleotide primers comprises SEQ ID NO: 23 and SEQ ID NO: 24.

25 18. A kit comprising a carrier to receive therein one or more
containers, at least one of said containers comprising an oligonucleotide primer according
to claim 5.

30 19. A kit comprising a carrier to receive therein one or more
containers, at least one of said containers comprising a pair of oligonucleotide primers
according to claim 6.

AMENDED CLAIMS

[received by the International Bureau on 27 April 2001 (27.04.01);
original claims 1 and 2 amended; remaining claims unchanged (1 page)]

1. An isolated double stranded nucleic acid which consists of a member from the group consisting of SEQ ID NOS: 6 to 9 and 12 to 13 and its complementary sequence.
2. An isolated nucleic acid which specifically hybridizes with a nucleic acid selected from the group consisting of SEQ ID NOS: 6 to 9 and their complementary sequences.
3. An oligonucleotide sequence for identification of a fermentation-related microorganism, wherein said sequence is selected from the group consisting of SEQ ID NOS: 16 to 24.
4. An oligonucleotide primer which is a fragment of the sequences according to claim 3, and which specifically hybridizes to the ITS1 or ITS2 of *Metschnikowia pulcherrima*; *Zygosaccharomyces bailii*; *Kluyveromyces thermotolerans*; *Candida stellata* and *Torulaspora delbrueckii*.
5. An oligonucleotide primer for identification of a fermentation-related microorganism, wherein said primer is selected from the group consisting of primers which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24, primers which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15, and primers of 10 nucleotide bases or longer which contain at least 5 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with from 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15.

6. A pair of oligonucleotide primers for use in the amplification-based detection of an internal transcribed spacer sequence of a fermentation-related microorganism, wherein the primers are selected from the group consisting of primers which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24, primers which contain at least 10 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15, and primers of 10 nucleotide bases or longer which contain at least 5 contiguous nucleotide bases from one of SEQ ID NOS: 16 to 24 contiguous with from 1 to 15 flanking nucleotide bases in the 5' and/or 3' direction of SEQ ID NOS: 6 to 15.

7. The pair of oligonucleotide primers according to claim 6 which comprises SEQ ID NO: 16 and SEQ ID NO: 17.

8. The pair of oligonucleotide primers according to claim 6 which comprises SEQ ID NO: 18 and SEQ ID NO: 19.

9. The pair of oligonucleotide primers according to claim 6 which comprises SEQ ID NO: 20 and SEQ ID NO: 21.

10. The pair of oligonucleotide primers according to claim 6 which comprises SEQ ID NO: 23 and SEQ ID NO: 24.

11. A method for detection of a fermentation-related microorganism comprising:

- (a) obtaining DNA from a fungal culture or colony isolated from a fermentation, or from an organism present in a fermentation beverage;
- (b) amplifying a part of the internal transcribed spacer sequence of said fermentation-related microorganism using said DNA as a template in a polymerase chain reaction with a pair of oligonucleotide primers according to claim 6; and
- (c) detecting said amplified part of the internal transcribed spacer sequence to determine whether said fermentation-related microorganism is present.

12. The method according to claim 11, wherein said fermentation-related microorganism is selected from the group consisting of *Metschnikowia pulcherrima*; *Zygosaccharomyces bailii*; *Kluyveromyces thermotolerans*; *Candida stellata* and *Torulaspora delbrueckii*.

13. The method according to claim 12, wherein said fermentation culture or fermentation beverage is a wine fermentation culture or wine fermentation beverage.

14. The method according to claim 12, wherein the pair of oligonucleotide primers comprises SEQ ID NO: 16 and SEQ ID NO: 17.

15. The method according to claim 12, wherein the pair of oligonucleotide primers comprises SEQ ID NO: 18 and SEQ ID NO: 19.

16. The method according to claim 12, wherein the pair of oligonucleotide primers comprises SEQ ID NO: 20 and SEQ ID NO: 21.

17. The method according to claim 12, wherein the pair of oligonucleotide primers comprises SEQ ID NO: 23 and SEQ ID NO: 24.

18. A kit comprising a carrier to receive therein one or more containers, at least one of said containers comprising an oligonucleotide primer according to claim 5.

19. A kit comprising a carrier to receive therein one or more containers, at least one of said containers comprising a pair of oligonucleotide primers according to claim 6.

- 1 -

SEQUENCE LISTING

<110> E. & J. Gallo Winery

<120> DETECTION OF FERMENTATION RELATED MICROORGANISMS

<130> 29520/35994

<140>

<141>

<160> 26

<170> PatentIn Ver. 2.0

<210> 1

<211> 344

<212> DNA

<213> Metschnikowia pulcherrima

<220>

<221> misc_feature

<222> (24)

<223> At position 24, Y = C or T

<220>

<221> misc_feature

<222> (56)

<223> At position 56, D = A, G, or T

<220>

<221> misc_feature

<222> (58)

<223> At position 58, W = A or T

<220>

<221> misc_feature

<222> (60)

<223> At position 60, R = A or G

<220>

<221> misc_feature

<222> (61)

<223> At position 61, Y = C or T

<220>

<221> misc_feature

<222> (70)

<223> At position 70, M = A or C

<220>

<221> misc_feature

<222> (74)

<223> At position 74, Y = C or T

<220>

<221> misc_feature

<222> (76)

<223> At position 76, R = A or G

<220>

<221> misc_feature

- 2 -

<222> (100)
 <223> At position 100, N = A, G, C, or T

<220>
 <221> misc_feature
 <222> (108)
 <223> At position 108, Y = C or T

<220>
 <221> misc_feature
 <222> (109)
 <223> At position 109, W = A or T

<220>
 <221> misc_feature
 <222> (110)
 <223> At position 110, Y = C or T

<220>
 <221> misc_feature
 <222> (111)
 <223> At position 111, W = A or T

<220>
 <221> misc_feature
 <222> (112)
 <223> At position 112, W = A or T

<220>
 <221> misc_feature
 <222> (114)
 <223> At position 114, H = A, C, or T

<220>
 <221> misc_feature
 <222> (115)
 <223> At position 115, M = A or C

<220>
 <221> misc_feature
 <222> (341)
 <223> At position 341, K = G or T

<220>
 <221> misc_feature
 <222> (342)
 <223> At position 342, Y = C or T

<400> 1
 ggaagtaaaaa gtcgtaacaa ggtytccgta ggtgaacctg cggaaggatc attaadawtr 60
 ytattacacm cttytrggca caaactctaa atcttaaccn tcaataaywy wwthmaaaaa 120
 ctttcaacaa cggatctctt gggtctcgca tcatgaaga acgcagcgaa ttgcgatacg 180
 taatatgact tgcagacgtg aatcattgaa tctttgaacg cacattgcgc cccgggggat 240
 tccccagggc atgcgtgggt gagcgatatt tactctcaaa cctccgggtt ggtcctgctt 300
 cggcctaata tcaacggcgc tagaataagt ttagcccca kyct 344

<210> 2
 <211> 372
 <212> DNA
 <213> Candida stellata

<400> 2

- 3 -

```

ctgaaggctt tttgccaaaa ccactgtgaa cagcttagac ttcgggtcttt gcaattgctt 60
gggtgtcgaa aggcgccccaa tctttaaaaac ttttatattt gttctgaaac aatgaaaatt 120
taaaactttc aacaacggat ctcttggttc tcgtatcgat gaagaacgca gcaaagcgcg 180
ataggtaatg cgaattgcag acgtgagtca ttgaattttt gaacgcata tgcgctatta 240
gtttgtctaa tagcatgctt gttggagtga taatcttcct ctcaaccatt tttggtatga 300
ggcttgcctc cttttaggag ttaaaatcat ggaagtgcac acgttaatta actctgtgca 360
gttatacact tt

```

372

<210> 3

<211> 697

<212> DNA

<213> *Zygosaccharomyces bailii*

<400> 3

```

tagaacaatt tctcgattga ctatctggac cagttgtggt ctttctgttt tttttaaggc 60
ctgcgcttaa ttgcgcggtc tagagcggag ggagttaagc atagttgctt tggctttcaa 120
tttacacaca gtggagtttc tacttttttt attcttcttt gggaggatgg gttcgtcccg 180
cttcagaggg taaacacaaa caattttttt tattttattt tattttatta ttataataat 240
aataatacag tcaaaacgaa tactaaaaca aaaaatattc aaaactttca acaacggatc 300
tcttggttct cgcacgatg aagaacgcag cgaaatgcga tacgtaattg gaattgcaga 360
attccgtgaa tcatcgaatc tttgaacgca cattgcgccc cttggtattc cagggggcat 420
gcctggttga gcgtcatttc cttctcaaac attcgtgttt ggtagttagt gatactctgt 480
tttattattt ggggttaactt gaaattgcaa gccttttggg gacgcgtgtg ggtgagtttt 540
aggcggaaac gtcttgctct cctctttcct aaccaaattg cgtattaggt tttaccgact 600
ccgacagacg ggactaggag attgggtgag tgatagcaat atcgagctct gcctaatttt 660
ttttttgcgc gccttgggca aacaatactc tcaaagt

```

697

<210> 4

<211> 588

<212> DNA

<213> *Kluyveromyces thermotolerans*

<220>

<221> misc_feature

<222> (38)

<223> At position 38, N = A, G, C or T

<220>

<221> misc_feature

<222> (530)

<223> At position 530, N = A, G, C, or T

<400> 4

```

aagaattttg ttagagcagc cgggaaagtt cagaagcntg cgcttgattg cgcggccgat 60
gatgctttct gttaacgact gtctctctac acacacactg tggagtaatt tattttacaa 120
cgcttcttct ttgggcttta cggcccaagg gttacaaaca caaacaacta ttgtatttta 180
aacattgtca attatttttc atttttagaaa aaaaatattt aaaactttca acaacggatc 240
tcttggttct cgcacgatg aagaacgcag cgaaatgcga taagtattgt gaattgcaga 300
tattcgtgaa tcatcgaatc tttgaacgca cattgcgccc tctggtattc cagggggcat 360
gcctggttga gcgtcatttc cttctcaaac cctcgggttt ggtagttagt ggtactcttt 420
ctgggttaac ttgaaaatgc tggccatctg gctggtgctg actgaggttt tagtccagtc 480
cgctgatact ctgcgtatta ggttttacca actcgtagtg gcgttagtan gcgtttttaa 540
ggcttttact gaaagtacag acagtctggc aaacagtatt cataaagt

```

588

<210> 5

<211> 710

<212> DNA

<213> *Torulaspora delbrueckii*

<400> 5

```

gagaaatcta tatgaatgaa gttagaggac gtctaaagat actgtaagag aggatcaggt 60
tcaagaccag cgcttaattg cgcggttgcg gcttggttcg ctttttgcg aacatgtctt 120

```

- 4 -

```

ttctcgttgt taactctact tcaacttcta caacactgtg gagttttcta cacaactttt 180
cttctttggg aagatacgtc ttgtgcgtgc ttcccagagg tgacaaacac aaacaacttt 240
ttattattat aaaccagtca aaaccaattt cgttatgaaa ttaaaaatat ttaaaacttt 300
caacaacgga tctcttgggt ctcgcatcga tgaagaacgc agcgaaatgc gatacgtaat 360
gtgaattgca gaattccgtg aatcatcgaa tctttgaacg cacattgcgc cccttggtat 420
tccagggggc atgectgttt gagcgtcatt tccttctcaa acaatcatgt ttggtagtga 480
gtgatactct gtcaagggtt aacttgaaat tgctagcctg ttatttgggt gtgattttgc 540
tggcttggat gactttgtcc agtctagcta ataccgaatt gtcgtattag gttttaccaa 600
cttcggcgaga ctgtgtgttg gctcggggcg tttaaagact ttgtcgtaaa cgatttatcg 660
tttgtttgag cttttcgcac acgcaatccg gcgaacaata ctctcaaagt 710

```

<210> 6

<211> 63

<212> DNA

<213> *Metschnikowia pulcherrima*

<220>

<221> misc_feature

<222> (2)

<223> At position 2, D= A, G, or T

<220>

<221> misc_feature

<222> (4)

<223> At position 4, W = A or T

<220>

<221> misc_feature

<222> (6)

<223> At position 6, R = A or G

<220>

<221> misc_feature

<222> (7)

<223> At position 7, Y = C or T

<220>

<221> misc_feature

<222> (16)

<223> At position 16, M = A or C

<220>

<221> misc_feature

<222> (20)

<223> At position 20, Y = C or T

<220>

<221> misc_feature

<222> (22)

<223> At position 22, R = A or G

<220>

<221> misc_feature

<222> (46)

<223> At position 46, N = A, G, C, or T

<220>

<221> misc_feature

<222> (54)

<223> At position 54, Y = C or T

<220>

- 5 -

<221> misc_feature
 <222> (55)
 <223> At position 55, W = A or T

<220>
 <221> misc_feature
 <222> (56)
 <223> At position 56, Y = C or T

<220>
 <221> misc_feature
 <222> (57)
 <223> At position 57, W = A or T

<220>
 <221> misc_feature
 <222> (58)
 <223> At position 58, W = A or T

<220>
 <221> misc_feature
 <222> (60)
 <223> At position 60, H = A, C, or T

<220>
 <221> misc_feature
 <222> (61)
 <223> At position 61, M = A or C

<400> 6
 adawtrytat tacacmctty trggcacaaa ctctaaatct taaccntcaa taaywywwth 60
 maa 63

<210> 7
 <211> 73
 <212> DNA
 <213> Metschnikowia pulcherrima

<220>
 <221> misc_feature
 <222> (70)
 <223> At position 70, K = G or T

<220>
 <221> misc_feature
 <222> (71)
 <223> At position 71, Y = C or T

<400> 7
 actctcaaac ctccggtttg gtcttgcttc ggcctaatat caacggcgct agaataagtt 60
 ttagccccak yct 73

<210> 8
 <211> 122
 <212> DNA
 <213> Candida stellata

<400> 8
 ctgaaggctt tttgccaaaa ccactgtgaa cagcttagac ttcgggtcttt gcaattgctt 60
 ggggtgtcgaa aggcgcccac tctttaaaac ttttatattt gttctgaaac aatgaaaatt 120
 ta 122

- 6 -

<210> 9
 <211> 89
 <212> DNA
 <213> *Candida stellata*

<400> 9
 aaccattttt ggtatgaggt cttgctcctt ttaggagtta aaatcatgga agtgcacacg 60
 ttaattaact ctgtgcagtt atacacttt 89

<210> 10
 <211> 281
 <212> DNA
 <213> *Zygosaccharomyces bailii*

<400> 10
 tagaacaatt tctcgattga ctatctggac cagttgtggt ctttctgttt tttttaaggc 60
 ctgcgcttaa ttgcgcggtc tagagcggag ggagttaagc atagttgctt tggctttcaa 120
 tttacacaca gtggagtttc tacttttttt attcttcttt gggaggatgg gttcgtcccg 180
 cttccagagg taaacacaaa caattttttt tattttattt tattttatta ttataataat 240
 aataatacag tcaaaacgaa tactaaaaca aaaaatattc a 281

<210> 11
 <211> 258
 <212> DNA
 <213> *Zygosaccharomyces bailii*

<400> 11
 ccttctcaaa cattcgtggt tggtagtgag tgatactctg ttttattatt tgggttaact 60
 tgaaattgca agccttttgg ggacgcgtgt gggtagtatt taggcggaaa cgtcttgctc 120
 tcctctttcc taaccaaagt tctgattagg ttttacggac tccgacagac gggactagga 180
 gattgggtga gtgatagcaa tatcgagctc tgcctaattt tttttttgcg cgccttgggc 240
 aaacaatact ctcaaagt 258

<210> 12
 <211> 221
 <212> DNA
 <213> *Kluyveromyces thermotolerans*

<220>
 <221> misc_feature
 <222> (38)
 <223> At position 38, N = A, G, C, or T

<400> 12
 aagaattttg ttagagcagc cgggaaagtt cagaagcntg cgcttgattg cgcggccgat 60
 gatgctttct gttaacgact gtctctctac acacacactg tggagtaatt tattttacaa 120
 cgcttcttct ttgggcttta cggcccaagg gttacaaaca caaacaacta ttgtatttta 180
 aacattgtca attatttttc atttttagaaa aaaaatattt a 221

<210> 13
 <211> 209
 <212> DNA
 <213> *Kluyveromyces thermotolerans*

<220>
 <221> misc_feature
 <222> (151)
 <223> At position 151, N = A, G, C or T

<400> 13
 ccttctcaaa ccctcggggt tggtagtgag tggtagtctt tctgggttaa cttgaaaatg 60
 ctggccatct ggctgttgct gactgaggtt ttagtccagt ccgctgatac tctgctgatt 120

- 7 -

agggttttacc aactcgtagt ggcgttagta ngcgttttta aggctttttac tgaaagtaca 180
 gacagtctgg caaacagtat tcataaagt 209

<210> 14

<211> 293

<212> DNA

<213> *Torulaspora delbrueckii*

<400> 14

gagaaatcta tatgaatgaa gtttagaggac gtctaaagat actgtaagag aggatcaggt 60
 tcaagaccag cgcttaattg cgcggttgcg gcttggttcg ccttttgccg aacatgtctt 120
 ttctcgttgt taactctact tcaacttcta caacactgtg gagttttcta cacaactttt 180
 cttctttggg aagatacgtc ttgtgcgtgc ttcccagagg tgacaaacac aaacaacttt 240
 ttattattat aaaccagtca aaaccaattt cgttatgaaa ttaaaaaatat tta 293

<210> 15

<211> 259

<212> DNA

<213> *Torulaspora delbrueckii*

<400> 15

ccttctcaaa caatcatggt tggtagttag tgatactctg tcaaggggta acttgaaatt 60
 gctagcctgt tatttggttg tgattttgct ggcttggtatg actttgtcca gtctagctaa 120
 taccgaattg tcgtattagg ttttaccacac ttccgcagac tgtgtggttg ctccggcgct 180
 ttaaagactt tgcgtgaaac gatttatcgt ttgtttgagc ttttcgcata cgcaatccgg 240
 cgaacaatac tctcaaagt 259

<210> 16

<211> 20

<212> DNA

<213> Oligonucleotide MXL258

<400> 16

aagcaggacc aaaccggagg 20

<210> 17

<211> 20

<212> DNA

<213> Oligonucleotide MXL 268

<400> 17

tattaggccg aagcaggacc 20

<210> 18

<211> 22

<212> DNA

<213> Oligonucleotide CsF

<400> 18

tttgccaaaa ccactgtgaa ca 22

<210> 19

<211> 21

<212> DNA

<213> Oligonucleotide CsR

<400> 19

tttaaagatt gggcgcttt c 21

<210> 20

<211> 18

<212> DNA

- 8 -

<213> Oligonucleotide QAZ165

<400> 20
tgggaggatg ggttcgtc 18

<210> 21
<211> 21
<212> DNA
<213> Oligonucleotide QAZ622

<400> 21
gctatcactc acccaatctc c 21

<210> 22
<211> 19
<212> DNA
<213> Oligonucleotide GSZ450

<400> 22
acctcagtca gcaacagcc 19

<210> 23
<211> 33
<212> DNA
<213> Oligonucleotide TdF1

<400> 23
ctatatgaat gaagtttagag gacgtctaaa gat 33

<210> 24
<211> 22
<212> DNA
<213> Oligonucleotide TdR1

<400> 24
ggaagcacgc acaagacgta tc 22

<210> 25
<211> 22
<212> DNA
<213> Oligonucleotide ITS5

<400> 25
ggaagtaaaa gtcgtaacaa gg 22

<210> 26
<211> 20
<212> DNA
<213> Oligonucleotide ITS4

<220>

<400> 26
tcctccgctt attgatatgc 20

INTERNATIONAL SEARCH REPORT

International application No.

PCT/US01/01458

A. CLASSIFICATION OF SUBJECT MATTER

IPC(7) : C07H 21/04; C12Q 1/70
 US CL : 536/23.1, 24.32, 24.3; 435/6

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)
 U.S. : 536/23.1, 24.32, 24.3; 435/6

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	JAMES et al. "Use of an rRNA Internal Transcribed Spacer Region to distinguish phylogenetically close related species of the general Zygosaccharomyces and Torulaspora" International J. of Systematic Bacteriology. January 1996, Vol 46, No. 1, pages 189-194.	1-2
Y		3-6, 9-13, 16-19
X	GODDARD et al. "Recurrent invasion and extinction of a selfish gene" PNAS, November 1999, Vol 96, No. 24, pages 13880-13885.	2
X	JAMES et al. Genbank Accession Number AJ224310 or AJ224312, January 1999.	1-2
Y		3-6, 10, 18-19
X	VAN NUES et al. "Separate structural elements within internal transcribed spacer 1 of Saccharomyces cerevisiae precursor rRNA" Nucleic Acids Research, 1994, Vol 22, No. 6, pages 912-919.	1-2
Y		3-6, 10, 18-19
X	VAN der SNADE et al. "Functional Analysis of ITS2 of Saccharomyces cerevisiae ribosomal DNA" J. Mol. Biol. 1992, Vol 22, pages 899-910.	1-2
Y		3-4
Y	Stratagene Catalog 1988, pg 39.	18-19

☐ Further documents are listed in the continuation of Box C.☐ See patent family annex.

* Special categories of cited documents:

"A" document defining the general state of the art which is not considered to be of particular relevance

"E" earlier application or patent published on or after the international filing date

"L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)

"O" document referring to an oral disclosure, use, exhibition or other means

"P" document published prior to the international filing date but later than the priority date claimed

"T"

later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention

"X"

document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone

"Y"

document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art

"&"

document member of the same patent family

Date of the actual completion of the international search

20 February 2001 (20.02.2001)

Date of mailing of the international search report

11 APR 2001

Name and mailing address of the ISA/US

Commissioner of Patents and Trademarks
 Box PCT
 Washington, D.C. 20231

Facsimile No. (703)305-3230

Authorized officer

Jeanine Enewold Goldberg

Telephone No. (703) 308-0196