



## INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

<b>(51) International Patent Classification <sup>6</sup> :</b>  <b>C12Q 1/68</b>	<b>A1</b>	<b>(11) International Publication Number:</b> <b>WO 98/55649</b>  <b>(43) International Publication Date:</b> 10 December 1998 (10.12.98)
<b>(21) International Application Number:</b> PCT/IB98/00865  <b>(22) International Filing Date:</b> 4 June 1998 (04.06.98)  <b>(30) Priority Data:</b> 08/871,678                      6 June 1997 (06.06.97)                      US  <b>(63) Related by Continuation (CON) or Continuation-in-Part (CIP) to Earlier Application</b> US    08/871,678 (CON) Filed on    6 June 1997 (06.06.97)  <b>(71) Applicant (for all designated States except US):</b> CHIRON DIAGNOSTICS CORPORATION [US/US]; 333 Coney Street, East Walpole, MA 02032 (US).  <b>(72) Inventors; and</b> <b>(75) Inventors/Applicants (for US only):</b> SANDHU, Gurpreet, S. [IN/US]; Apartment 212, 10320 Devonshire Circle, Bloomington, MN 55431 (US). KLINE, Bruce, C. [US/US]; 2315 Kline Lane, S.W., Rochester, MN 55902 (US).  <b>(74) Agents:</b> MORGENSTERN, Arthur, S.; Chiron Diagnostics Corporation, 63 North Street, Medfield, MA 02052 (US) et al.		<b>(81) Designated States:</b> AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GE, GH, GM, GW, HU, ID, IL, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZW, ARIPO patent (GH, GM, KE, LS, MW, SD, SZ, UG, ZW), Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG).  <b>Published</b> <i>With international search report.</i>
<b>(54) Title:</b> NUCLEIC ACID PROBES FOR THE DETECTION AND IDENTIFICATION OF FUNGI  <b>(57) Abstract</b>  Nucleic acid probes and primers are described for detecting fungi that cause disease in humans and animals, as well as spoilage of food and beverages. These probes can detect rRNA, rDNA or polymerase chain reaction products from a majority of fungi in clinical, environmental or food samples. Nucleic acid hybridization assay probes specific for <i>Paracoccidioides brasiliensis</i> and <i>Pneumocystis carinii</i> are also described.		

**FOR THE PURPOSES OF INFORMATION ONLY**

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
AT	Austria	FR	France	LU	Luxembourg	SN	Senegal
AU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav Republic of Macedonia	TM	Turkmenistan
BF	Burkina Faso	GR	Greece			TR	Turkey
BG	Bulgaria	HU	Hungary	ML	Mali	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MN	Mongolia	UA	Ukraine
BR	Brazil	IL	Israel	MR	Mauritania	UG	Uganda
BY	Belarus	IS	Iceland	MW	Malawi	US	United States of America
CA	Canada	IT	Italy	MX	Mexico	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NE	Niger	VN	Viet Nam
CG	Congo	KE	Kenya	NL	Netherlands	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NO	Norway	ZW	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's Republic of Korea	NZ	New Zealand		
CM	Cameroon			PL	Poland		
CN	China	KR	Republic of Korea	PT	Portugal		
CU	Cuba	KZ	Kazakstan	RO	Romania		
CZ	Czech Republic	LC	Saint Lucia	RU	Russian Federation		
DE	Germany	LI	Liechtenstein	SD	Sudan		
DK	Denmark	LK	Sri Lanka	SE	Sweden		
EE	Estonia	LR	Liberia	SG	Singapore		

5 NUCLEIC ACID PROBES FOR THE DETECTION AND IDENTIFICATION  
OF FUNGI

FIELD OF INVENTION

The inventions described and claimed herein relate to the design and  
10 composition of two nucleic acid probes capable of detecting many different fungal  
organisms in clinical, food, environmental and other samples. The inventions  
described and claimed herein also relate to the design and composition of probes  
capable of specifically detecting and identifying *Acremonium* sp., *Aspergillus*  
*clavatus*, *Aspergillus flavus*, *Aspergillus fumigatus*, *Aspergillus glaucus*, *Aspergillus*  
15 *nidulans*, *Aspergillus niger*, *Aspergillus ochraceus*, *Aspergillus terreus*, *Aspergillus*  
*unguis*, *Aspergillus ustus*, *Beauveria* sp., *Bipolaris* sp., *Blastoschizomyces* sp.,  
*Blastomyces dermatitidis*, *Candida albicans*, *Candida glabrata*, *Candida*  
*guilliermondii*, *Candida kefyr*, *Candida krusei*, *Candida lusitanae*, *Candida*  
*parapsilosis*, *Candida tropicalis*, *Chrysosporium* sp., *Cladosporium* sp., *Coccidioides*  
20 *immitis*, *Cryptococcus neoformans* var *gattii* serotype B, *Cryptococcus neoformans*  
serotype A, *Cryptococcus laurentii*, *Cryptococcus terreus*, *Curvularia* sp., *Fusarium*  
sp., *Filobasidium capsuligenum*, *Filobasidiella* (*Cryptococcus*) *neoformans* var  
*bacillispora* serotype C, *Filobasidiella* (*Cryptococcus*) *neoformans* var *neoformans*  
serotype D, *Filobasidium uniguttulatum*, *Geotrichum* sp., *Histoplasma capsulatum*,  
25 *Malbranchea* sp., *Mucor* sp., *Paecilomyces* sp., *Paracoccidioides brasiliensis*,  
*Penicillium species*, *Pneumocystis carinii*, *Pseudallescheria boydii*, *Rhizopus* sp.,  
*Sporothrix schenkii*, *Scopulariopsis brevicaulis* sp., *Scopulariopsis brumpti*,  
*Saccharomyces cerevisiae*, and *Trichosporon beigelii* in clinical, food, environmental  
and other samples.

30 Fungi are eukaryotic microorganisms that are universally distributed. While in  
nature fungi play a major role in the decomposition of plant materials, they are also  
responsible for spoilage of food, beverage and pharmaceutical preparations. Out of an

estimated 100,000 species of fungi described by mycologists, approximately 150 species are pathogenic to man and animals. The increasing incidence of AIDS and the development of newer treatments for hematologic malignancies and organ transplants has lead to an increase in the number of immunocompromised patients. These patients  
5 have a high risk of developing fungal infections, which if not rapidly diagnosed and treated are capable of causing death in a matter of days. The number of antifungal drugs is limited and their toxic side effects on the patient are much higher than that of comparable antibacterial therapy. A rapid diagnosis of fungal infection and start of treatment is critical in these patients. Books by Kwon-Chung and Bennett, along with  
10 Sarosi and Davies, provide an overview into the medical importance of fungi.

Fungal organisms are identified by morphology and nutritional characteristics. Fungi may take anywhere from two days to several weeks to grow in culture and often the same organism can take radically different forms depending on the growth conditions. This makes timely identification difficult even for the classically trained  
15 expert and impedes the treatment of patients where rapid identification of genus and species is of medical advantage.

The incidence and distribution of major pathogenic fungi varies by geographic location. *Aspergillus fumigatus*, *Blastomyces dermatitidis*, *Candida albicans*, *Coccidioides immitis*, *Cryptococcus neoformans*, *Histoplasma capsulatum*,  
20 *Paracoccidioides brasiliensis*, *Pseudallescheria boydii* and *Sporothrix schenckii* represent some of the leading causes of mycotic infections.

*Aspergillus fumigatus* is among the top three causes of systemic fungal infection treated in hospitals. It usually affects patients with organ transplants, acute leukemias and burns and can be rapidly fatal if not diagnosed quickly. With over 150  
25 species of *Aspergillus* present in the soil, air and water, accurate detection of *Aspergillus fumigatus* becomes extremely important. *Aspergillus clavatus*, *Aspergillus flavus*, *Aspergillus fumigatus*, *Aspergillus glaucus*, *Aspergillus nidulans*, *Aspergillus niger*, *Aspergillus ochraceus*, *Aspergillus terreus*, *Aspergillus unguis* and *Aspergillus ustus* represent a majority of *Aspergillus* species seen in clinical specimens and their  
30 presence can cause diagnostic difficulties. *Aspergillus flavus*, *Aspergillus fumigatus* and *Aspergillus niger* have been linked with disease in humans, with *Aspergillus fumigatus* being the predominant pathogen in North America. A few immunologic

tests exist for *Aspergillus fumigatus* but these have limited sensitivity and specificity. There are also reports of development of polymerase chain reaction based tests for *Aspergillus fumigatus* based on the amplification of the *Asp fl* antigen gene and a ribosomal intergenic spacer (Spreadbury et. al.). The Spreadbury technique is based on the PCR amplification of a 401 bp fragment spanning the large subunit rRNA/intergenic spacer region. This relies on a pair of primers to specifically amplify DNA from *Aspergillus fumigatus* only, and is of no utility in identifying other fungi.

*Blastomyces dermatitidis* is present in the soil, usually in bird droppings and animal feces. Infections often occur at construction sites and the ensuing lung infiltration and pneumonitis are usually fatal in immunocompromised patients. Diagnosis by culture may take weeks, and the organism is occasionally mistaken for other fungi. Existing immunological diagnostic tests are unreliable, and there is a need for rapid and reliable DNA based diagnostic tests. Similarly, *Histoplasma capsulatum* exists in the soil and is known to have infected at least 20% of the population of North America. Most infections start in the lung and resolve spontaneously, but may occasionally spread to other organs. AIDS patients represent a growing number of cases of Histoplasmosis. Diagnosis is difficult as immunological tests are often negative during the first 4-6 weeks of infection. *Coccidioides immitis* is found in abundance in the soil in Southwestern United States. Dust storms, farming, building construction, earthquakes and even hiking have been linked with outbreaks of disease. Lung infection followed by cavitation and disseminated miliary coccidioidomycosis are seen. Meningitis is usually lethal, and as with other fungi, mortality is highest in debilitated hosts. Four serotypes of *Cryptococcus neoformans* cause disease in humans. These are *Cryptococcus neoformans* serotype A, *Cryptococcus neoformans* var *gatti* serotype B, *Filobasidiella (Cryptococcus) neoformans* var *bacillispora* serotype C and *Filobasidiella (Cryptococcus) neoformans* var. *neoformans* serotype D. The incidence of this disease is growing rapidly, with up to 10% of HIV infected people developing cryptococcosis. DNA probes capable of detecting all 4 serotypes are required for the early diagnosis and treatment for life threatening infections like cryptococcal meningitis. A report by Stockman et. al. discusses commercial tests for *Histoplasma*, *Blastomyces*, *Coccidioides*, and *Cryptococcus* based on the 18S rRNA (Gen-Probe, Inc., San Diego, CA). The authors report sensitivities ranging from 87.8

to 100% and a specificity of 100%. One drawback of these probes is that these are used on rRNA extracted from fungal cultures. As some fungi may require up to 3 weeks to grow in culture, this technique cannot be used to expedite diagnosis until a culture becomes available.

5        *Candida albicans* is one of the most common causes of fungal infection in humans. It is present in the respiratory, gastrointestinal and female genital tract of healthy individuals, and acts as an opportunistic pathogen in debilitated individuals on steroid or chemotherapy. Diabetes mellitus and indwelling catheters are other predisposing causes. Immunocompromised hosts show rapid hematogenous spread of  
10 fungi. Morbidity and mortality in untreated cases is high. *Candida glabrata*, *Candida guilliermondii*, *Candida kefyr*, *Candida krusei*, *Candida lusitanae*, *Candida parapsilosis* and *Candida tropicalis* are also known to cause disease in humans. DNA probes capable of identifying these individual species would eliminate the need for multiple blood cultures and lengthy biochemical speciation.

15        Paracoccidioidomycosis, which is a deep-seated systemic infection of humans, is a major health problem in Central and South America. The disease is caused by *Paracoccidioides brasiliensis*, a thermally dimorphic fungus. Classically, diagnosis has been made by detecting yeast cells with multiple buds, the distinctive "pilots wheel" morphology. Unfortunately, such forms are not always observed in clinical materials.  
20 Histologically, the yeast form of this species is sometimes difficult to distinguish from that of *Histoplasma capsulatum*, *Blastomyces dermatitidis*, *Coccidioides immitis*, and *P. loboi*. While *P. brasiliensis* produces a glycoprotein antigen than is also diagnostic, the antigen is also produced by *P. loboi*. *P. brasiliensis* can be cultured in vitro; however *P. loboi* cannot be cultured. *P. brasiliensis* is a slow growing organism that can take one to  
25 three weeks to culture for identification. Clearly, the existence of a species specific probe against *P. brasiliensis* would be of significant value for rapid detection and identification.

*Pneumocystis carinii* is a common and major opportunistic pathogenic fungus that infects immunocompromised patients. More than 80% of HIV infected patients  
30 will develop *P. carinii* pneumonia if not prophylactically treated to prevent infection. Therefore, this fungus is widely recognized as one of the most important causes of morbidity and mortality in AIDS patients. To date, polymerase chain reaction primers

directed at six different gene targets have been developed. These include 5S and 18 rDNA genes, an intragenic transcribed spacer sequence between ribosomal genes, thymidylate synthase, dihydrofolate reductase, and large ribosomal rDNA gene of mitochondria. However, each set of primers requires a unique and specific set of PCR conditions to function optimally (Lu et al., J. Clin. Microbial. 33:2785-2788, 1995).

Recent advances in molecular techniques have led to the approach of microbe detection and identification based upon the DNA sequence of ribosomal genes. Commonly used detection techniques include either direct amplification of the ribosomal DNA (rDNA) genes by the polymerase chain reaction, or reverse transcription of the ribosomal RNA (rRNA) into complementary DNA (cDNA) followed by polymerase chain reaction amplification of the cDNA. Ribosomes are composites of unique rRNA and protein species that function in the translation of messenger RNA into protein. Evolutionary studies are consistent with the interpretation that all extant life has evolved from a single organism. Thus, all cellular organisms contain rRNA and these rRNAs are related by evolution. The evolutionary process is such that each species of organism appears to have unique regions of sequence in its ribosomal genes. The presence of these unique species specific regions allows one to design DNA probes that under conditions of hybridization will specifically bind to, and identify the polymerase chain reaction amplified DNA from only one species of fungus. For the purposes of this application, the word "primer" is used to mean a nucleotide sequence which can be extended by template-directed polymerization, and "probe" is used to mean a nucleotide sequence capable of detecting its complementary sequence by hybridization. Also, for the purpose of this application, the phrase "nucleotide sequence" is intended to include either DNA or RNA forms or modification thereof. Furthermore, those versed in the art will recognize that primer sequences can be used as probes and vice versa. The use of nucleic acid hybridization to detect specific nucleic acid sequences of interest is also described by Kohne (U.S. Patent 4,851,330, 7/1989).

In prokaryotes and eukaryotes, ribosomal RNA and the corresponding rDNA genes are identified by the size of the RNA. The sizes are related in terms of sedimentation velocity or S values. Thus, for prokaryotes the values are 5S, 16S, and 23S; and for eukaryotes the values are 5S, 5.8S, 18S and 28S. Because all ribosomes

perform the same function which is essential for cell viability, ribosomal sequences are largely conserved, yet certain regions of each ribosomal species are subject to more variation without consequence to function. It is these hypervariable regions that allow one to identify different species amongst members of the same genus. As noted  
5 in the references, there are several reports where 5S, 18S and the intergenic spacer between 5.8S and 28S rDNA have been used for the detection and identification of fungi (Holmes et. al., Hopfer et. al., Lott et. al., Maiwald et. al., Makimura et. al., Mitchell et. al., Nakamura et. al.). Holmes et. al. describe a PCR test based on the co-amplification of the 5S rDNA and an adjacent nontranscribed spacer region. This  
10 identifies only *Candida albicans* and detects other *Candida* species without identifying individual organisms. Hopfer et. al. and Maiwald et. al. both use universal primers to amplify 18S rDNA from several fungi including *Candida* sp., *Aspergillus fumigatus*, *Cryptococcus neoformans* and *Trichosporon* sp. These amplicons are digested with restriction enzymes and the cut fragments are sized by gel  
15 electrophoresis. This restriction fragment length polymorphism pattern enables them to identify most but not all organisms. This technique can be used on amplified DNA from a pure fungal culture. As clinical samples such as sputum usually contain multiple fungal organisms, this technique has little utility in diagnosis as multiple overlapping fragments obtained from a mix of fungi would be nearly impossible to  
20 interpret. Lott et. al. use the 5.8S RNA and the internal transcribed spacer (ITS2) to identify and speciate *Candida albicans* and related *Candida* species. Makimura amplifies a 687 bp fragment from the 18S rDNA of 25 medically important fungi and uses these in the diagnosis of *Candida albicans* in clinical samples. Mitchell uses nested PCR to amplify 5.8S and internal transcribed spacer (ITS) to identify  
25 *Cryptococcus neoformans*. No subsequent testing is done to verify the identity of the amplified DNA. Nakamura et. al. use 18S primers to detect *Aspergillus fumigatus* infections of the lung. Most protocols given in these references can only be used to detect an extremely limited number of fungi from a clinical specimen. Hopfer et. al. and Maiwald et. al. can identify multiple organisms from pure cultures, but their  
30 utility for clinical specimens containing multiple fungal species is limited at best.

United States patents have been issued to Weisburg et. al. for probes developed for the detection of 18S small subunit ribosomal RNA sequences in fungi.



These probes will detect fungi from many species, but cannot be used easily to identify any single species. United States patents have also been issued to Milliman for probes developed for the specific detection of the bacteria *Staphylococcus aureus* based on the 16S ribosomal sequences. Hogan et. al. (European Pat. App. 0,272,009) describe one fungal probe for 18S rRNA and three fungal probes for 28S rRNA sequences. Two of these 28S probes detect several different fungi while the third probe detects *Candida krusei* from a limited panel of 10 fungi. None of the 28S probes described by Hogan et. al. is related to any of the probes described in our invention. All probes claimed in our invention can be mapped within the first 900 base pairs of a 28S gene. The probes described by Hogan et. al. are located further 3' on the 28S sequence, between base pairs 1000 and 2000 (these numbers are comparable to the primary sequence of *Saccharomyces cerevisiae* 28S rRNA gene. Genbank accession number: J01355). Leclerc et. al. have published reports analyzing the phylogenetic relationship between fungi based on partial DNA sequences of several fungal 28S genes sequenced by them. Some of the organisms claimed to have been sequenced by Leclerc are the same as some organisms sequenced by us. These are *Sporothrix schenckii*, *Pseudallescheria boydii*, *Blastomyces dermatitidis*, *Histoplasma capsulatum* and *Chrysosporium* sp. Leclerc et. al. have not published any sequence data in their report, and to the best of our knowledge, they have not made these sequences publically available. The reverse-complement sequence of their sequencing primer 401 (TCCCTTTCAA CAATTTCACG) overlaps our SEQ ID NO: 1 (GTGAAATTGT TGAAAGGGAA) by 19 nucleotides and their sequencing primer 636 (GGTCCGTGTT TCAAGACGG) overlaps our SEQ ID NO: 2 (GACTCCTTGG TCCGTGTT) by 10 nucleotides. We are aware of no reports in the literature of variable regions from 28S rRNA genes of fungi being used as targets for the development of species specific diagnostic probes.

As discussed above, most present techniques for the molecular detection of fungi rely on the use of highly specific primers for the PCR amplification of only one fungal species. Those that employ "Universal" primers for a PCR amplification of DNA from multiple organisms, use post-PCR amplicon identification techniques that are useful only on pure cultures of fungi. These are not be able to identify fungi from a clinical specimen containing multiple fungal organisms. Our first aim was to

develop "Universal" primers for the 28S gene. These primers would be capable of amplifying in a PCR, 28S rDNA from most fungi. Our subsequent aim was to develop species specific probes for fungi of interest, that would be used to analyze our "Universal" 28S amplicon. These species specific probes would be able to detect the presence of fungi of interest even in situations containing mixed fungal species.

One aspect of this invention is to provide nucleic acid primers capable of detecting 28S sequences from DNA or RNA of most fungi. These would be used as "Universal" primers in a polymerase chain reaction to amplify 28S sequences from any fungus present in clinical, food, environmental or other samples. These "Universal" primers would also be used to sequence the amplified DNA. The sequence obtained would be used to identify the fungus by comparing with a database of known fungal sequences.

A second aspect of this invention is to provide nucleic acid probes capable of detecting and identifying, by nucleic acid hybridization, the pathogens *Aspergillus fumigatus*, *Blastomyces dermatitidis*, *Candida albicans*, *Coccidioides immitis*, *Cryptococcus neoformans*, *Histoplasma capsulatum*, *Aspergillus flavus*, *Aspergillus glaucus*, *Aspergillus niger*, *Aspergillus terreus*, *Candida glabrata*, *Candida guilliermondii*, *Candida kefyr*, *Candida krusei*, *Candida lusitanae*, *Candida parapsilosis*, *Candida tropicalis*, *Paracoccidioides brasiliensis*, *Pneumocystis carinii*, *Pseudallescheria boydii*, *Sporothrix schenckii* and other species by use of any of several different formats. Additionally, nucleotide sequence information is provided to identify these pathogens and other fungi by DNA sequence comparison (Table 3) or by the construction of additional probes.

## SUMMARY OF THE INVENTION

Nucleic acid probes and primers are described for detecting fungi that cause disease in humans and animals, as well as spoilage of food and beverages. These probes can detect rRNA, rDNA or polymerase chain reaction products from a majority of fungi in clinical, environmental or food samples. Nucleic acid hybridization assay probes specific for *Aspergillus fumigatus*, *Blastomyces dermatitidis*, *Candida albicans*, *Coccidioides immitis*, *Cryptococcus neoformans*, *Histoplasma capsulatum*, *Aspergillus flavus*, *Aspergillus glaucus*, *Aspergillus niger*, *Aspergillus terreus*, *Candida glabrata*, *Candida guilliermondii*, *Candida kefyr*, *Candida krusei*, *Candida lusitanae*, *Candida parapsilosis*, *Candida tropicalis*, *Paracoccidioides brasiliensis*, *Pneumocystis carinii*, *Pseudallescheria boydii*, *Sporothrix schenckii* and other species (Table 1 and Table 3) are also described.

Figure 1 represents the relative position of the sequences described on the 28S subunit of fungi.

## DETAILS OF THE INVENTION

Our first objective was to develop nucleic acid primers for use in a polymerase chain reaction to amplify 28S genes from all fungi likely to be present in a clinical sample. This amplified DNA would then be amenable to probing with several different species specific probes. Each one of these species specific probes would, under conditions of hybridization, anneal to 28S ribosomal DNA from only one species of fungus, thereby detecting and identifying the species of fungus present in the clinical sample. The 28S gene was selected as a target because it had regions that were conserved among fungi and these would provide potential annealing sites for "universal" fungal probes. The ribosomal 28S genes were also expected to have hypervariable regions that would be unique enough to provide sites for species specific probes. The large rRNA gene is called the 23S rRNA gene in prokaryotes and 28S in eukaryotes. This designation is based on the length and therefore the sedimentation coefficient of these rRNA molecules. Fungal large subunit rRNAs vary in size among different organisms and are often referred to as being 25S, 26S or 28S.

Since fungi are eukaryotes, and to maintain uniformity in this application, we shall refer to fungal large subunit rRNA as 28S rRNA.

Published sequences from *Cryptococcus neoformans*, two *Candida albicans*, *Saccharomyces cerevisiae* and two *Schizosaccharomyces pombe* 28S genes are approximately 3.5 kilobases in length (Genbank accession numbers: L14068, L28817, X70659, J01355, Z19136 & Z19578). These four sequences were aligned, and a region of sequence variability was found clustered between coordinates 200 and 700 from the 5' end of these genes. As an initial starting point, two nucleic acid primers P1 (ATCAATAAGC GGAGGAAAAG) (SEQ ID NO:79) and P2 (CTCTGGCTTC ACCCTATTC) (SEQ ID NO:80) (see figure 1), capable of hybridizing to all 4 of the above mentioned organisms and not to human 28S sequences (GenBank accession number: M11167), were designed and used under low stringency hybridization conditions in a polymerase chain reaction to amplify approximately 800 base pairs of DNA spanning this hypervariable region from the following 34 fungi that were obtained from the Mayo Clinic fungal collection: *Acremonium* sp., *Aspergillus clavatus*, *Aspergillus fumigatus*, *Aspergillus glaucus*, *Aspergillus nidulans*, *Aspergillus niger*, *Aspergillus ochraceus*, *Aspergillus terreus*, *Aspergillus unguis*, *Aspergillus ustus*, *Beauveria* sp., *Bipolaris* sp., *Blastomyces dermatitidis*, *Candida albicans*, *Candida glabrata*, *Candida guilliermondii*, *Candida kefyr*, *Candida krusei*, *Candida lusitanae*, *Candida parapsilosis*, *Candida tropicalis*, *Chrysosporium* sp., *Cladosporium* sp., *Coccidioides immitis*, *Cryptococcus neoformans* serotype A, *Curvularia* sp., *Geotrichum* sp., *Histoplasma capsulatum*, *Mucor* sp., *Penicillium* sp., *Pseudallescheria boydii*, *Saccharomyces cerevisiae*, *Sporothrix schenckii* and *Trichosporon beigelii*.

DNA was extracted from the fungi listed above by the following method. A loopful of fungal culture was scraped off a culture plate using a sterile inoculation loop. The fungus was added one milliliter of sterile water in a 1.5 ml Sarsted (Newton, North Carolina) screw cap microcentrifuge tube. This tube was placed in a boiling water bath for 20 minutes in order to lyse the fungus and release DNA from the cells. Two microliters of this whole cell lysate was used in a PCR to amplify 28S rDNA. All PCR amplifications were carried out as hot-start reactions in a 50 ul reaction volume using Perkin-Elmer (Norwalk, CT) 0.5 ml thin-wall polypropylene tubes and a Perkin-

Elmer thermal cycler. Reagents added to the tube initially were 2.5 ul of 10X PCR buffer (100 mM tris pH 8.3, 500 mM KCl, 15 mM MgCl<sub>2</sub>), 5.0 ul of 50% glycerol/1 mM cresol red, 8.0 ul of dNTP mix (1.25 mM each of dATP, dGTP, dTTP and dCTP), 12 picomoles of each nucleic acid primer and sterile water to make up a  
5 volume of 25 ul. A wax bead (Ampliwax Gem-100, Perkin-Elmer) was added and the tubes heated to 77°C for 1 minute and cooled to room temperature to form a wax barrier. 2.5 ul of 10X PCR buffer, 5.0 ul of 50% glycerol/1 mM cresol red, 0.2 ul Taq polymerase (AmpliTaq 5U/ul, Perkin-Elmer) and 15.3 ul of sterile water was added to the tube along with 2.0 ul of DNA from the fungal whole cell lysate described above.  
10 50 cycles of thermal cycling was carried out at 94°C - 30 sec, 40°C - 1 min, 72°C - 2 min. The amplified DNA was electrophoresed and purified from a low melt agarose gel by tris buffered phenol pH 8.0, phenol/chloroform/isoamyl alcohol (25:24:1 by vol.) and 3 ether extractions, followed by isopropanol precipitation and 70% ethanol wash.

15 We completely sequenced both strands of DNA amplified from the organisms listed above. All sequencing was carried out on an Applied Biosystems 373A sequencer. Every nucleotide in the sequences generated was verified and confirmed by examining the complementary nucleotide from the second strand sequence. We had now created a novel database consisting of nucleic acid sequences spanning a  
20 variable region of the 28S rDNA from a diverse collection of medically important fungi.

While the complete sequences for *Candida albicans*, *Cryptococcus neoformans* and *Saccharomyces cerevisiae* 28S genes had previously been published and deposited in GenBank, it was not obvious, nor had it been defined, whether any  
25 regions of sequence identity among these three organisms would also be conserved among all fungi of interest. DNA sequences from all the fungi in our novel 28S database had to be analyzed in order to develop "Universal" 28S probes. All sequences were subjected to extensive manipulation to identify optimal relative alignments in order to identify regions of similarity for use as "Universal" probes.  
30 The selected probe sequences had to meet several important criteria besides the condition of being present in 28S genes from most fungal species. Each probe sequence required an appropriate thermal profile, secondary structure and utility in a

DNA amplification reaction. These probes were optimized to work for PCR amplification in pure cultures of fungus, as well as in the presence of DNA from multiple sources as in the case of clinical specimens. The probes were also designed to facilitate direct sequencing of the amplified DNA. Our analysis led to the discovery of the oligonucleotide probes listed in (SEQ ID NO:1) and (SEQ ID NO:2). (For their location, see Figure 1.) The successful identification of these two probes ((SEQ ID NO:1) and (SEQ ID NO:2)) completed our first objective to develop nucleic acid probes that would hybridize to, and detect 28S rRNA and rDNA from a majority of fungi (Figure 1 and Table 1). As shown later in this application, the novel sequence information generated by the use of our "Universal" probes allowed us to develop species-specific probes ((SEQ ID NO:3) to (SEQ ID NO:23), (SEQ ID NO:75) and (SEQ ID NO:76)) capable of identifying 21 different disease-causing fungi.

Table 1:

Presence of hybridization sites for probes SEQ ID NO: 1 and SEQ ID NO: 2 in 28S nucleic acid sequences.

	SEQ ID NO: 1	SEQ ID NO: 2
<i>Acremonium</i> sp.	+	+
<i>Aspergillus clavatus</i>	+	+
<i>Aspergillus flavus</i>	+	+
<i>Aspergillus fumigatus</i>	+	+
<i>Aspergillus glaucus</i>	+	+
<i>Aspergillus nidulans</i>	+	+
<i>Aspergillus niger</i>	+	+
<i>Aspergillus ochraceus</i>	+	+
<i>Aspergillus terreus</i>	+	+
<i>Aspergillus unguis</i>	+	+
<i>Aspergillus ustus</i>	+	+
<i>Beauveria</i> sp.	+	+

<i>Bipolaris</i> sp.	+	+
<i>Blastomyces dermatitidis</i>	+	+
<i>Blastoschizomyces</i> sp.	+	+
<i>Candida albicans</i>	+	+
<i>Candida glabrata</i>	+	+
<i>Candida guilliermondii</i>	+	+
<i>Candida kefyr</i>	+	+
<i>Candida krusei</i>	+	+
<i>Candida lusitaniae</i>	+	+
<i>Candida parapsilosis</i>	+	+
<i>Candida tropicalis</i>	+	+
<i>Chrysosporium</i> sp.	+	+
<i>Cladosporium</i> sp.	+	+
<i>Coccidioides immitis</i>	+	+
<i>Cryptococcus laurentii</i>	+	+
<i>Cryptococcus neoformans</i> serotype A	+	+
<i>Cryptococcus neoformans</i> var. <i>gattii</i> serotype B	+	+
<i>Cryptococcus terreus</i>	+	+
<i>Curvularia</i> sp.	+	+
<i>Filobasidiella</i> ( <i>Cryptococcus</i> ) <i>neoformans</i> var <i>bacillispora</i> serotype C	+	+
<i>Filobasidiella</i> ( <i>Cryptococcus</i> ) <i>neoformans</i> var <i>neoformans</i> serotype D	+	+
<i>Filobasidium capsuligenum</i>	+	+
<i>Filobasidium uniguttulatum</i>	+	+
<i>Fusarium</i> sp.	+	+
<i>Geotrichum</i> sp.	+	+

<i>Histoplasma capsulatum</i>	+	+
<i>Malbranchea</i> sp.	+	+
<i>Mucor</i> sp.	+	+
<i>Paecilomyces</i> sp.	+	+
<i>Penicillium</i> sp.	+	+
<i>Pseudallescheria boydii</i>	+	+
<i>Rhizopus</i> sp.	+	+
<i>Saccharomyces cerevisiae</i>	+	+
<i>Scopulariopsis brevicaulis</i>	+	+
<i>Scopulariopsis brumptii</i>	+	+
<i>Sporothrix schenckii</i>	+	+
<i>Trichosporon beigelii</i>	+	+
<i>Paracoccidioides brasiliensis</i>	+	+
<i>Pneumocystis carinii</i>	+	+
Human	-	+

Probes SEQ ID NO: 1 and SEQ ID NO: 2 were used to successfully amplify (Table 2) and sequence DNA (Table 3) spanning this variable region from the following 51 organisms: *Acremonium* sp., *Aspergillus clavatus*, *Aspergillus flavus*,  
5 *Aspergillus fumigatus*, *Aspergillus glaucus*, *Aspergillus nidulans*, *Aspergillus niger*,  
*Aspergillus ochraceus*, *Aspergillus terreus*, *Aspergillus unguis*, *Aspergillus ustus*,  
*Beauveria* sp., *Bipolaris* sp., *Blastomyces dermatitidis*, *Blastoschizomyces* sp.,  
*Candida albicans*, *Candida glabrata*, *Candida guilliermondii*, *Candida kefyr*,  
*Candida krusei*, *Candida lusitanae*, *Candida parapsilosis*, *Candida tropicalis*,  
10 *Chrysosporium* sp., *Cladosporium* sp., *Coccidioides immitis*, *Cryptococcus*  
*neoformans* serotype A, *Cryptococcus neoformans* var. *gattii* serotype B,  
*Cryptococcus terreus*, *Cryptococcus laurentii*, *Curvularia* sp., *Filobasidiella*  
*(Cryptococcus) neoformans* var. *bacillispora* serotype C, *Filobasidiella*  
*(Cryptococcus) neoformans* var. *neoformans* serotype D, *Filobasidium capsuligenum*,  
15 *Filobasidium uniguttulatum*, *Fusarium* sp., *Geotrichum* sp., *Histoplasma capsulatum*,  
*Malbranchea* sp., *Mucor* sp., *Paecilomyces* sp., *Penicillium* sp., *Pseudallescheria*



*boydii*, *Rhizopus* sp., *Saccharomyces cerevisiae*, *Scopulariopsis brevicaulis*, *Scopulariopsis brumptii*, *Sporothrix schenckii*, *Paracoccidioides brasiliensis*, *Pneumocystis carinii*, and *Trichosporon beigelii*. This list contains all 4 serotypes (A, B, C and D) of *Cryptococcus neoformans*. This sequence information generated by the use of probes SEQ ID NO: 1 and SEQ ID NO: 2 expanded the size of our database consisting of fungal 28S sequences. All amplified DNA was sequenced across both strands from a minimum of two different isolates of each organism to ensure accuracy of the data generated.

Table 2:

Polymerase chain reaction amplification of 28S rDNA with probes SEQ ID NO: 1 and SEQ ID NO: 2.

	PCR with SEQ ID NO: 1 & NO: 2
<i>Acremonium</i> sp.	+
<i>Aspergillus clavatus</i>	+
<i>Aspergillus flavus</i>	+
<i>Aspergillus fumigatus</i>	+
<i>Aspergillus glaucus</i>	+
<i>Aspergillus nidulans</i>	+
<i>Aspergillus niger</i>	+
<i>Aspergillus ochraceus</i>	+
<i>Aspergillus terreus</i>	+
<i>Aspergillus unguis</i>	+
<i>Aspergillus ustus</i>	+
<i>Beauveria</i> sp.	+
<i>Bipolaris</i> sp.	+
<i>Blastomyces dermatitidis</i>	+
<i>Blastoschizomyces</i> sp.	+
<i>Candida albicans</i>	+

<i>Candida glabrata</i>	+
<i>Candida guilliermondii</i>	+
<i>Candida kefyr</i>	+
<i>Candida krusei</i>	+
<i>Candida lusitaniae</i>	+
<i>Candida parapsilosis</i>	+
<i>Candida tropicalis</i>	+
<i>Chrysosporium</i> sp.	+
<i>Cladosporium</i> sp.	+
<i>Coccidioides immitis</i>	+
<i>Cryptococcus laurentii</i>	+
<i>Cryptococcus neoformans</i> serotype A	+
<i>Cryptococcus neoformans</i> var. <i>gattii</i> serotype B	+
<i>Cryptococcus terreus</i>	+
<i>Curvularia</i> sp.	+
<i>Filobasidiella</i> ( <i>Cryptococcus</i> ) <i>neoformans</i> var. <i>bacillispora</i> serotype C	+
<i>Filobasidiella</i> ( <i>Cryptococcus</i> ) <i>neoformans</i> var. <i>neoformans</i> serotype D	+
<i>Filobasidium capsuligenum</i>	+
<i>Filobasidium uniguttulatum</i>	+
<i>Fusarium</i> sp.	+
<i>Geotrichum</i> sp.	+
<i>Histoplasma capsulatum</i>	+
<i>Malbranchea</i> sp.	+
<i>Mucor</i> sp.	+
<i>Paecilomyces</i> sp.	+

<i>Penicillium</i> sp.	+
<i>Pseudallescheria boydii</i>	+
<i>Rhizopus</i> sp.	+
<i>Saccharomyces cerevisiae</i>	+
<i>Scopulariopsis brevicaulis</i>	+
<i>Scopulariopsis brumptii</i>	+
<i>Sporothrix schenckii</i>	+
<i>Trichosporon beigelii</i>	+
<i>Paracoccidioides brasiliensis</i>	+
<i>Pneumocystis carinii</i>	+
Human	-

This list of fungi sequenced by us represents organisms responsible for most cases of subcutaneous and deep mycotic infections in humans and also includes saprophytes (non-pathogenic fungi) commonly encountered in clinical isolates. Since the two probes (SEQ ID NO: 1 and SEQ ID NO: 2) hybridize to 28S rDNA from all the fungi listed above, they are capable of diagnosing the presence of a majority of fungi that are likely to be present in a clinical specimen. They are believed to be primers for universally detecting fungi.

Probes listed in SEQ ID NO: 1 and SEQ ID NO: 2 were also checked for their potential ability to hybridize to, and amplify (in a polymerase chain reaction) 23S sequences from bacteria by searching for hybridization sites among the 539 bacterial 23S genes listed in GenBank. Bacterial 23S rDNAs do not have suitable hybridization sites for SEQ ID NO: 1 and SEQ ID NO: 2 and these two probes should not be able to amplify bacterial DNA under stringent conditions.

Our second objective was to develop species specific probes, which under hybridization conditions, would detect *Aspergillus fumigatus*, *Blastomyces dermatitidis*, *Candida albicans*, *Coccidioides immitis*, *Cryptococcus neoformans*, *Histoplasma capsulatum*, *Aspergillus flavus*, *Aspergillus glaucus*, *Aspergillus niger*, *Aspergillus terreus*, *Candida glabrata*, *Candida guilliermondii*, *Candida kefyr*, *Candida krusei*, *Candida lusitanae*, *Candida parapsilosis*, *Candida tropicalis*,

*Paracoccidioides brasiliensis*, *Pneumocystis carinii*, *Pseudallescheria boydii*, and *Sporothrix schenckii*. We used our database of fungal 28S nucleic acid sequences to create a multiple sequence alignment of all the organisms that we had sequenced.

Every individual sequence was subjected to intensive comparison with all other sequences in our database in order to discover unique regions of sequence that would be present only in the fungus of interest, and would be absent in all other fungi. When unique stretches of sequence were identified, these were further analyzed for thermal profile and secondary structure. Each probe constructed by us will, under conditions of hybridization, specifically hybridize to and detect, nucleic acid sequence from the unique region of only one specific target fungus. Those versed in the art will recognize that specification of a single-stranded DNA sequence implies the utility of the complementary DNA sequence, as well as the two equivalent RNA sequences. Furthermore, sequences incorporating modification of any of the moieties comprising the nucleic acid (i.e., the base, the sugar or the backbone) are functional equivalents of the sequence. It should also be recognized that these additional sequences can potentially serve as probes or primers. Finally, those versed in the art recognize that comparisons of extensive DNA sequences provides enough variability and uniqueness to speciate organisms (Table 3).

The nucleic acid sequences for these species specific synthetic probes are listed in SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76.

There are two probes specific for *Cryptococcus neoformans*, two probes specific for *Sporothrix schenckii*, and one probe each for *Aspergillus fumigatus*, *Blastomyces dermatitidis*, *Candida albicans*, *Coccidioides immitis*, *Histoplasma capsulatum*, *Aspergillus flavus*, *Aspergillus glaucus*, *Aspergillus niger*, *Aspergillus terreus*, *Candida glabrata*, *Candida guilliermondii*, *Candida kefyr*, *Candida krusei*, *Candida lusitanae*, *Candida parapsilosis*, *Candida tropicalis*, *Paracoccidioides brasiliensis*, *Pneumocystis carinii*, and *Pseudallescheria boydii* 28S rRNA and rDNA. (See Tables 4 - 8 and further discussion below.)

All species specific probes developed by us are novel and to the best of our knowledge have not been reported in the literature. While all 28S genes sequenced by us had several regions that were different among the various species analyzed, the regions that would function best as species specific probes under conditions of

hybridization were not obvious. Extensive analysis of each 28S sequence yielded several potential probe sites. These were studied in detail to enable the selection of optimal unique sites for each probe, based on the need to obtain optimal hybridization characteristics under the test conditions. The highly specific hybridization characteristics of all probe sequences developed by us were then validated by experimental results. The prior existence in GenBank of sequences for *Candida albicans* and serotypes A and B (GenBank accession numbers L14067 and L14068) of *Cryptococcus neoformans* 28S genes was in itself not sufficient to enable even an individual versed in this field to develop specific probes for either of these two organisms. We had to obtain novel 28S sequence from *Candida albicans*, *Candida glabrata*, *Candida guilliermondii*, *Candida kefyr*, *Candida krusei*, *Candida lusitanae*, *Candida parapsilosis*, *Candida tropicalis*, *Cryptococcus neoformans* serotype A, *Cryptococcus neoformans* var. *gattii* serotype B, *Cryptococcus terreus*, *Cryptococcus laurentii*, *Filobasidiella (Cryptococcus) neoformans* var *bacillispora* serotype C, *Filobasidiella (Cryptococcus) neoformans* var *neoformans* serotype D, *Filobasidium capsuligenum* and *Filobasidium uniguttulatum* before we were able to identify potential regions for the development of species specific probes for these two fungal organisms that would not cross react with the others listed above.

Our modification of the Chomczynski technique (see Example 2, below) allows us to obtain DNA from any clinical specimen, irrespective of source (see Table 10 for a variety of clinical specimens tested), within a 3 hour period. The PCR amplification and subsequent probing can be accomplished with ease within a 24 hour period. The final identification is therefore possible in a day as opposed to several days or weeks required by traditional methods. This speed and sensitivity of diagnosis can make a difference between life and death in debilitated patients battling fungal diseases of undetermined cause. Rapid diagnosis will allow physicians to immediately direct their therapy towards curing the identified causative fungus, rather than wait for days or weeks while the patient succumbs to an unknown fungus.

Our probes have the ability to pick out the correct target organism even in a mixed fungal infection because of their high level of specificity. The methods of Hopfer et. al. and Maiwald et. al., do not allow identification of individual species in a mixed fungal infection because restriction fragment length polymorphism results are

nearly impossible to interpret when multiple organisms contribute to the restriction fragments. Their method can therefore only be used on a pure culture, and this also does not save any diagnostic time, because the fungus first has to be grown in culture.

The probes developed by us allow rapid species identification of a large  
5 number of pathogenic fungi by using multiple probes against only one PCR amplified fragment of DNA. Coupled with our modified DNA extraction technique and our ability to accurately diagnose in the case of mixed organisms, this strategy can provide the greatest amount of diagnostic information in the shortest amount of time. This diagnostic strategy is also amenable to automation, which can result in even  
10 greater savings in time, money and effort.

The sequences and the complement of the sequences claimed in this disclosure, along with any modifications to these sequences, may potentially be utilized in assays for the identification of fungi based on several existing methodologies, as well as future improvements and alterations of this technology.  
15 These techniques include, but are not limited to, assays based on hybridization, ligation, polymerization, depolymerization, sequencing, chemical degradation, enzymatic digestion, electrophoresis, chromatography and amplification. Furthermore, all such variations ultimately are based in some selection or amplification process, some ligand or some nucleic acid moiety that recognizes or  
20 utilizes the sequences (SEQ ID NO: 1) to (SEQ ID NO:23), (SEQ ID NO:75) and (SEQ ID NO:76) claimed in this application. Such variations include but are not limited to use of a variety of linear or exponential target amplification schemes, such as, any of the myriad forms of PCR, the ligase chain reaction, Q-beta repliase, etc.; direct detection of species-specific nucleic acid purified or extracted from pure fungal  
25 culture using a probe selected from the group (SEQ ID NO: 3) to (SEQ ID NO: 23), (SEQ ID NO:75) and (SEQ ID NO:76); use of the complementary DNA forms of (SEQ ID NO:1) to (SEQ ID NO:23), (SEQ ID NO:75) and (SEQ ID NO:76); use of the RNA forms of these sequences and their complements; and use of derivatives of these DNA or RNA sequences by the addition of one or more reporter moieties from a  
30 variety of labels including nucleic acid sequences, proteins, signal generating ligands such as acridinium esters, and/or paramagnetic particles. These techniques may be

utilized with DNA, RNA or modified derivatives used as either the target or the detection molecule.

In addition to the 25 sequences SEQ ID NO: 1 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76, we also describe an additional 53 sequences SEQ ID NO: 24 to SEQ ID NO: 74, SEQ ID NO:77 and SEQ ID NO:78. These 53 sequences are inclusive of SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76 and are shown as a multiple sequence alignment (Table 3) with coordinate 1 corresponding to base # 431 of a reference *S. cerevisiae* 28S rRNA gene. (The numbers are comparable to the primary sequence of *S. cerevisiae* 28S rRNA gene. Genbank accession number: J01355). These sequences were obtained by amplifying and sequencing 28S rDNA from various fungi with primers SEQ ID NO: 1 and SEQ ID NO: 2. (SEQ ID NO: 1 corresponds to coordinates 403-422 and the SEQ ID NO: 2 corresponds to coordinates 645-662 of the reference *S. cerevisiae* gene).

An analysis of these aligned sequences enabled us to develop the species specific probes SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76 and sites for these probes are shown underlined. These 53 aligned sequences contain sufficient variability, to enable a person versed in this art, to develop additional species specific hybridization probes in the 10-50 nucleotide length. Similarly, longer species specific hybridization probes encompassing the entire 200+ nucleotide length can also be envisioned. Species identification may also be accomplished by direct DNA sequence determination of any DNA amplified with primers SEQ ID NO: 1 and SEQ ID NO: 2. If the derived sequence matches approximately 98% or more of any sequence in SEQ ID NO: 24 to SEQ ID NO: 74, SEQ ID NO:77 and SEQ ID NO:78, then the identity of the organism can be ascertained. Additionally, we recognize that parts of SEQ ID NO: 24 to SEQ ID NO: 74, SEQ ID NO:77 and SEQ ID NO:78 may be specific for groups of fungi arranged phylogenetically at the level of genus or higher. SEQ ID NO: 24 to SEQ ID NO: 74, SEQ ID NO:77 and SEQ ID NO:78, their complements, along with any modification to these sequences may also potentially be utilized in assays for the identification of fungi based on existing methodologies and future technologies as noted above for SEQ ID NO: 1 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76.

Table 3:

Multiple sequence alignment for (SEQ ID NO: 24) through  
(SEQ ID NO: 74), (SEQ ID NO: 77) and (SEQ ID NO: 78)

```

1
(Rhizo2) AGCCAGACTG GCTTGTCTGT AATCAATCTA GGCTTTCG GC CTGGATGCAC TTGCAGGCTA ...TGCTTGC
(Rhizo3) AGCCAGACTG GCTTGTCTGT AATCAGICTA AGCTTTCG GC TTGGATGCAC TTGCAGGCTA ...TGCTTGC
(Rhizo1) AGCCAGACTG GCTTGTCTGT AATCAATCTA GGTTTCGTGC CTGGATGCAC TTGCAGACTA TTGCTTGC
(Mucor_) AGCCAGACTG GTTTCAGTGT AATCAAOCTA GAATTCGTTC .TGGGTGCAC TTGCAGTCTA ...TAOCTGC
(C_Terr) AGTCAGTCAT GTCTATTGGA CTCAGOOGGT TCT.....G COGGTGTACT TOCTTTAGAT GGGGTCAAC.
(F_Caps) AGTCAGTCAT GTCTATTGGA CTCAGOOGGT TCT.....G COGGTGTATT TOCTTTAGAT GGGGTCAAC.
(F_Unig) AGTCAGTOGT GTCTAATGGA CTCAGOOG.. TTC.....T GOGGTGTATT TOCAITGGGT GGGGTCAAC.
(C_Neob) AGTCAGTOGT GTCTATTGGG TTCAGOCAGC TCT.....G CTGGTGTATT OOCITTAGA. OGGGTCAAC.
(F_Neoc) AGTCAGTOGT GTCTATTGGG TTCAGOCAGC TCT.....G CTGGTGTATT OOCITTAGA. OGGGTCAAC.
(F_Neod) AGTCAGTOGT GTCTATTGGG TTCAGOCAGT TCT.....G CTGGTGTATT OOCITTAGA. OGGGTCAAC.
(C_Neof) AGTCAGTOGT GTCTATTGGG TTCAGOCAGT TCT.....G CTGGTGTATT OOCITTAGA. OGGGTCAAC.
(T_Beig) AGTCAGTOGT GTTCTTTGGA TTCAGOCAGT TCT.....G CTGGTGTACT TOCTTGGAA. OGGGTCAAC.
(C_Laur) AGTCAGTOGT GTCTGGGAGG CTCAGOOGGT TCT.....G COGGTGTATT OCTCTAGA. OGGGTCAAC.
(Beauve) GAOCAGACTT GGGCTTGGTT GATCATOOGG GGTTTC.TOC. OGGTGCAC CTTC.GGOC CAGGOCAGC.
(Fusari) GAOCAGACTT GGGCTTGGTT AATCATOOGG GGTTTC.TCY. OCAGTGCAC TTTC.AGTC CAGGOCAGC.
(Acremo) GAOCAGACTT GGGCTOOGGT AATCATOOGG OGTTTC.TOG. OGGTGCAC TTGOC.GTTC CAGGOCAGC.
(Paecil) GAOCAGACTT GGGOOOGGTG GATCATOCAG OGTTTC.TOG. CTGGTGCAC OGGOOOGGT CAGGOCAGC.
(P_Boyd) GAOCAGACTT GTGOOOGTGG AATCAGOOGC OGCTOGTGG. GOGGOGCACT TOGGOGGGCT CAGGOCAGC.
(S_Brum) GAOCAGACTC GOGOOOGTGG GATCAGOOGT OGCTOGTGG. GOGGOGCACT OGGOGGGGT CAGGOCAGC.
(S_Brev) GAOCAGACTT GOGOOOGTGG GATCAAOOGT OGCTTTC.GG. GOGGOGCACT OGGOGGGGT CAGGOCAGC.
(Sporot) GAOCAGACTT GOGOCYOGGG GAOTAOOGGG OGTTTC.TOG. OGGTGCAC CTGOGKGG CAGGOCAGC.
(B_Derm) GAOCAGACTC GGOOGTGGGG GTTCAGOGGG CATTOGT.TG OOGTGCAC OOOOAOOGG OGGGOCAGC.
(H_Caps) GAYCAGACTC GGOOGYGGGG GTTCAGOGGG CATTOGT.TG OOGTGCAC OOOOAOOGG OGGGOCAGC.
(A_Nidu) GAOCAGACTC GGOOOC.GGG GTTCAROCAG CACTOG..TG CTGGTGTACT TCOOOGGGGG OGGGOCAGC.
(A_Ungu) GAOCAGACTC GGOCTC.GGG GTTCAGOCAG CACTOG..TG CTGGTGTACT TCOOOGGGGG OGGGOCAGC.
(A_Ustu) GAOCAGACTC GGOOOC.GGG GTTCAGOCAG CACTOG..TG CTGGTGTACT TCOOOGGGGG OGGGOCAGC.
(A_Clav) GAOCAGACTC GGOOOC.GGG GTTCAGOOGG CATTOG..TG OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(A_Fumi) GAOCAGACTC GGOOOC.GGG GTTCAGOOGG CATTOG..TG OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(A_Flav) GAOCAGACTC GGOOOC.GGG GTTCAGOOGG CATTOG..TG OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(A_Ochr) GAOCAGACTC GGOOOC.GGG GTTCAGOOGG CATTOG..TG OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(A_Nige) GAOCAGACTC GGOOOC.GGG GTTCAGOOGG CATTOG..TG OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(A_Terr) AAOAGACTC GCTOOC.GGG GTTCAGOOGG CATTOG..TG OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(A_Glau) GAOCAGACTC GCTOOC.GGG GTTCAGOOGG CATTOG..TG OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(Penici) GAOCAGACTC GGOOCAC.GGG GTTCAGOOGG CATTOG..TG OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(C_Immi) AAOAGACTC GGTGTGGGG GCTCAGOGGG CATGAGT.GC OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(Bipola) AAOAGACTT GCTTGCAGTT GCTCATOOGG GCTT.T.GC OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(Curvtl) AAOAGACTT GCTTGCAGTT GCTCATOOGG GCTT.T.GC OGGTGTACT TCOOOGTGGG OGGGOCAGC.
(Chryso) AAOAGACTT GOGOGOGGG GATCATOOGG TCTTTC.T.GA OGGTGTACT OGGOGGTGG CAGGOCAGC.
(Cladoc) AAOAGACTT GCTOOGGGT. GTTCAGOOGG TCTTTC.T.GA OGGTGTACT OGGOGGTGG CAGGOCAGC.
(Malbra) AGACAGACTC GAGOGOGGG GCTCAGOGGG TATTGTTATG OOGTGTACT OOOOGOGGG OGGGOCAGC.
(C_Para) GATCAGACTT GGTATTGTTG ATG..TTACT CTCTOOGGG. ..GTGGGCTC TACAGTTTAC OGGGOCAGC.
(C_Trop) GATCAGACTT GGTATTGTTG ATG..TTACT CTCTOOGGG. ..GTGGGCTC TACAGTTTAC OGGGOCAGC.
(C_Albi) GATCAGACTT GGTATTGTTG ATG..CTGCT. CTCTOOGGG. ..GGGGOOGC TGGGTTTAC OGGGOCAGC.
(C_Guil) GATCAGACTC GATATTGTTG GAGGCTTGGC TCTGGGGG. ..GGGTTGOC OGCAGCTTAC OGGGOCAGC.
(C_Glab) GATCAGACTT GGTATTGTTG GGOOCTTGGC TCTGGGGG. TGGGACTCT OGCAGCTTAC TGGGOCAGC.
(S_Cere) GATCAGACTT GGTATTGTTG GGOOCTTGGC TCTGGGGG. TGGGACTCT OGCAGCTTAC TGGGOCAGC.
(C_Kefy) GATCAGACTT GGTATTGTTG GGOOCTTGGC TCTGGGGG. TGGGACTCT OGCAGCTTAC TGGGOCAGC.
(Geotri) AATCAGACTT GGTGCTGT.. .TGTTCACCT RTGTTTGGC ATAGTGTACT CAGCAGTACT AGGOCAGG.
(C_Lusi) AAGCAGACAC GGT..... .TTTAC OGGGOCAGC.
(C_Krus) OGOOOGACTT GGGGTTGGG CAGGCTTGGC TCTGGGGG. ..GGGCTCT GGGCTTTCOC TGGGOCAGC.
(Blasch) .....
(P_Braz) ACCAGAGTCG GCCGCGGGGG CTCAGCGGGC ACTCGT.TGC CCGTGCACTC CCCCCTG... .GTCGGGGCCA
(Pc_Hum) ATCAGACATG CCTTTATAGG AGATGCCATT GTT..TCGGC ATTGGCAGTA TTATCCGAAT TGGCAGGCCA

```



71

140

(Rhizo2) AACGACAAATT TGACTTGAGG GAAAAAATA GGGGAAATGT GGCC..... CACTTGTGGG TGTATATAGTC  
 (Rhizo3) AACGACAAATT TGGCTTGAGG GAAAAAATA AGGGAAATGT GGCC..... CATCGTGGG TGTATATAGTC  
 (Rhizo1) AACGACAAATT TTTTITGAGT GTAAAAATA TTGGAAATGT GGCCAAATAT TATTATITGG TGTATATAGTC  
 (Mucor) AACMACAGTT TGATTTGGAG GAAAAAATA GTAGGAATGT AGCC..... TCTOGA GGTGTATAGTC  
 (C\_Terr) .ATCAGTTTT .GATOGCTGG AAAAGGGCAG GAGGAATGTA GCACTC.TOG GGTGAACCTA TAGCCTTCTG  
 (F\_Caps) .ATCAGTTTT .GACOGTTGG ATAAAGGCAG GAGGAATGTA GCACTC.TOG GGTGAACCTA TAGCCTTCTG  
 (F\_Unig) .ATCAGTTTT .GATOGCTGG ATAAAGGCAG GAGGAATGTA GCACTC.TOG GGTGAACCTA TAGCCTTCTG  
 (C\_Neob) .ATCAGTTCT .GATOGGTGG ATAAAGGCAG GAGGAATGTA GCACTC.TOG GGTGTGTGTA TAGCCTTCTG  
 (F\_Neoc) .ATCAGTTCT .GATOGGTGG ATAAAGGCAG GAGGAATGTA GCACTC.TOG GGTGTGTGTA TAGCCTTCTG  
 (F\_Neod) .ATCAGTTCT .GATOGGTGG ATAAAGGCAG GAGGAATGTA GCACTC.TOG GGTGTGTGTA TAGCCTTCTG  
 (C\_Neof) .ATCAGTTCT .GATOGGTGG ATAAAGGCAG GAGGAATGTA GCACTC.TOG GGTGTGTGTA TAGCCTTCTG  
 (T\_Beig) .ATCAGTTTT .GTOGGTGG ATAAAGGTAG TAGGAATGTA .ACTTCTOC GGAAGTGTA TAGCCTTCTG  
 (C\_Laur) .ATCAGTTTT .GTOGGAOGG ATAAAGGTAG TAGGAATGTA GCACTC.TOC GGTGTGTGTA TAGCCTTCTG  
 (Beauve) .ATCAGTTGG CCGT.GGGGG ACAAGGGCTT CCGGAAGTGT GCTCTCTOC. ....GGGA .....  
 (Fusari) .ATCAGTTTT CSOC.GGGGG ATAAAGGCTT CCGGAATGTA GCTCTCTOC. ....GGGA .....  
 (Acemo) .ATCAGTTGG CGOC.GGGGG ATAAAGGCTT CCGGAATGTA GCTCTCTOC. ....GGGA .....  
 (Paecil) .ATCAGTTGG CGOC.GGGGG ATAAAGGCTT CCGGAAGTGT GCTCTCTOC. ....GGGA .....  
 (P\_Boyd) .ATCAGTTGG CTGCAGGGGG ACAAGGGCTT CCGGAATGTA GCTCTCTOC. ....GGGA .....  
 (S\_Brum) .ATCAGTTGG CTOGGGGGG AGAAGGGCTT CCGGAATGTA GCTCTCTOC. ....GGGA .....  
 (S\_Brev) .ATCAGTTGG .TOGGGGGG AGAAGGGCTT CCGGAATGTA GCTCTCTOC. ....GGGA .....  
 (Sporot) .ATOGTTTCT C.OCAGGGGG ACAAGGGCTT CCGGAATGTA GCTCTCTOC. ....GGGA .....  
 (B\_Derm) .GTGGTTTC .GACGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (H\_Caps) .GTGGTTTC .GACGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Nidu) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Ungu) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Ustu) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Clav) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Fumi) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Flav) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Ochr) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Nige) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Terr) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (A\_Glau) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (Penici) .GTGGTTTC .GGGGGCGG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (C\_Immi) .ATCAGTTCT .GGGGGCTG TTAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (Bipola) .ATCAGTTTC .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (Curvul) .ATCAGTTTC .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (Chryso) .ATOGTTTT .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (Clados) .ATOGTTTT .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (Halbra) .ATCAGTTTC .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (C\_Para) .ATCAGTTTC .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (C\_Trop) .ATCAGTTTC .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (C\_Albi) .ATOGTTTT .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (C\_Guil) .ATOGTTTT .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (C\_Glab) .ATOGTTTT .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (S\_Cere) .ATCAGTTTC .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (C\_Kefy) .ATCAGTTTC .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (Geotri) .TGGGGTGT .TGGGGTGT GAAAAAGTAG TAGGAATGTA ACTCTCTOC. ....GGG.C .....  
 (C\_Lusi) .GTC.GAAAA .GGGGGCTG ATAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (C\_Krus) .ATOGTTTT .TGGGGTGT GAAAAAGTAG TAGGAATGTA ACTCTCTOC. ....GGG.C .....  
 (Blasch) .....  
 (P\_Braz) GCGTCGGTTT CGACGGCCCG TCAAAGGCTT CCGGAATGTA TOGCTCTOC. ....GGG.C .....  
 (Pc\_Hum) GCATCGGTTT CAGTTACTGG ATAAAGGCTT GAGGAATGTA GCTCTCTOC. ....GGG.C .....

141 210

(Rhizo2) CCTTAGAAAA TACCTTGGGT TGGATTGAGG AACGCAGCGA ATG..... CTTATTG  
 (Rhizo3) CCTTAGAAAA TACCTTGGGC TGGATTGAGG TACGCAGCGA ATG..... CTATTG  
 (Rhizo) CTTTAGAAAA TACCTTGAAT TGGATTGAGG AACGCAGCGA ATGCTTCTCT TTnGAGGCAA AGTCTTTTAT  
 (Mucor) OCTACTATCA TACTCTGGAT TGGACTGAGG AAGGCAGCGA ATGCTTCTCT GCRAGATTGC TGGGTGCTTT  
 (C\_Terr) TGGTATACAG TGGTTGGGAC TGAGGAAAGC AGCATGOCIT TATGGGCGGG GTTGGGOCAC GTACATGCTT  
 (F\_Caps) TCACATACAA TGGTTGGGAC TGAGGAAAGC AGCATGOCIT TATGGGCGGG ATTGGTGCAC GTACATGCTT  
 (F\_Unig) TCACATACAG TGGTTGGGAC TGAGGAAAGC AGCATGOCIT TATGGGCGGG ATTGGTGCAC GTACATGCTT  
 (C\_Neob) TOGCATACAC TGGTTGGGAC TGAGGAAAGC AGCTGOCIT TATGGGCGGG GTTGGGOCAC GTTGGAGCTT  
 (F\_Neoc) TOGCATACAC TGGTTGGGAC TGAGGAAAGC AGCTGOCIT TATGGGCGGG GTTGGGOCAC GTTGGAGCTT  
 (F\_Neod) TOGCATACAC TGGTTGGGAC TGAGGAAAGC AGCTGOCIT TATGGGCGGG GTTGGGOCAC GTTGGAGCTT  
 (C\_Neof) TOGCATACAC TGGTTGGGAC TGAGGAAAGC AGCTGOCIT TATGGGCGGG GTTGGGOCAC GTTGGAGCTT  
 (T\_Beig) TCACATACAC TGGGTGAGAC TGAGGACTGC AGCTGOCIT TATGGGCGGG GTTGGGOCAC GTTGGAGCTT  
 (C\_Laur) TOGCATACAC TGGGTGAGAC TGAGGACTGC AGCTGOCIT TATGGGCGGG GTTGGGOCAC GTTGGAGCTT  
 (Beauve) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CATTOGCA  
 (Fusari) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CWTCTGCA  
 (Acremo) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... C.TCTGCA  
 (Paecil) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... C.TGCGCA  
 (P\_Boyd) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CATCTGCA  
 (S\_Brum) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... OGTTATGCA  
 (S\_Brev) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... OGTTATGCA  
 (Sporot) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (B\_Derm) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (H\_Caps) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Nidu) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Ungu) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Ustu) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Clav) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Fumi) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Flav) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Ochr) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Nige) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Terr) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (A\_Glau) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (Penici) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (C\_Immi) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (Bipola) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (Curvul) .....G TGTATAGOC OGTTGGTAA TACC.CTGG GGGACTGAG GTTGGG... CTTGGGCA  
 (Chryso) .....T GTTATAGC.C TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (Clados) .....T TATA.G OCTCTGTA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (Malbra) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (C\_Para) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (C\_Trop) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (C\_Albi) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (C\_Guil) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (C\_Glab) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (S\_Cere) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (C\_Kefy) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (Geotri) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (C\_Lusi) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (C\_Krus) .....G TGTATAGOC TGGGTGCAA TGGGOCAGC TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (Blasch) .....TGAAT TGTGTAAG GGAAGGOCAT GGTAGGAATA AGAGGCTGOG GTTGGGAATA ATTGTTTTTC  
 (P\_Braz) .....CG TGTATAGCC GGGGGTGCAA TGGGOCAGC GGGGOCAG GTTGGG... CTTGGGCA  
 (Pc\_Hum) TAGCTGCAGT GACCGGGACC GGAAGGGAAA TTGGGTCTTT GAAGACCTTA TGATGTTGCC AGAAATGGTC

	211	250		
{Rhizo2}	GCGAGTTTTC CAGGAAGGT. ....TTTCT	GAGGTACTAC	SEQ	ID NO: 68
{Rhizo3}	GCGAGTTGGC TGGGAATAT. ....TTTCT	GAGGTGCTTT	SEQ	ID NO: 69
{Rhizo1}	TGGGATTAC GGATCAGAC. ....TGTGG	CATTGTCACA	SEQ	ID NO: 67
{Mucor}	CGCTAATAAA TGTTAGAATT TCTGCTTOGG	GTGGTGCTAA	SEQ	ID NO: 63
{C_Terr}	AGG..ATGTT GACATAATGG CTTTAAAOGA	COOGTCTTGA	SEQ	ID NO: 53
{F_Caps}	AGG..ATGTT GACATAATGG CTTTAAAOGA	COOGTCTTGA	SEQ	ID NO: 56
{F_Unig}	AGG..ATGTT GACATAATGG CTTTAAAOGA	COOGTCTTGA	SEQ	ID NO: 59
{C_Neob}	AGG..ATGTT GACAAAATGG CTTTAAAOGA	COOGTCTTGA	SEQ	ID NO: 50
{F_Neoc}	AGG..ATGTT GACAAAATGG CTTTAAAOGA	COOGTCTTGA	SEQ	ID NO: 57
{F_Neod}	AGG..ATGTT GACAAAATGG CTTTAAAOGA	COOGTCTTGA	SEQ	ID NO: 58
{C_Neof}	AGG..ATGTT GACAAAATGG CTTTAAAOGA	COOGTCTTGA	SEQ	ID NO: 51
{T_Beig}	AGG..ATGTT GACATAATGG CTTTAAAOGA	COOGTCTTGA	SEQ	ID NO: 74
{C_Laur}	AGG..ATGTT GAOGTAATGG CTTTAAAOGA	COOGTCTTGA	SEQ	ID NO: 48
{Beauve}	AGG..ATGCT GGOGTAATGG TCATCAGTGA	COOGTCT...	SEQ	ID NO: 35
{Fusari}	AGG..ATGCT GGOGTAATGG TCATCAGTGA	COOGTCTTGA	SEQ	ID NO: 55
{Acreno}	AGG..ATGCT GGOGTAATGG TCATCAGTGA	COOGTCTTGA	SEQ	ID NO: 24
{Paecil}	AGG..ATGCT GGOGTAATGG TCATCAGOGA	COOGTCTTGA	SEQ	ID NO: 64
{P_Boyd}	AGG..ATGCT GGOGTAATGG TOGTCAGOGA	COOGTCTTGA	SEQ	ID NO: 66
{S_Brum}	AGG..ATGCT GGOGTAATGG TOGTCAGOGA	COOGTCTTGA	SEQ	ID NO: 72
{S_Brev}	AGG..ATGCT GGOGTAATGG TOGTCAGOGA	COOGTCTTGA	SEQ	ID NO: 71
{Sporot}	AGG..ATGCT GGOGTAATGG TCAOCAGOGA	ACOGTCTTGA	SEQ	ID NO: 70
{B_Derm}	CGG..AOGCT GGCCTAATGG TOGTAAOGGA	COOGTCTTGA	SEQ	ID NO: 38
{H_Caps}	CGG..AOGCT GGCCTAATGG TOGTCAOGGA	COOGTCTTGA	SEQ	ID NO: 61
{A_Nidu}	CGG..AOGCT GGOGTAATGG TOGCAAOGA	COOGTCTTGA	SEQ	ID NO: 29
{A_Ungu}	CGG..AOGCT GGCATAATGG TTGCAAOGA	COOGTCTTGA	SEQ	ID NO: 33
{A_Ustu}	CGG..AOGCT GGOGTAATGG TOGCAAOGA	COOGTCTTGA	SEQ	ID NO: 34
{A_Clav}	CGG..AOGCT GGOGTAATGG TOGTAAATGA	COOGTCTTGA	SEQ	ID NO: 25
{A_Fumi}	CGG..AOGCT GGOGTAATGG TOGTAAATGA	COOGTCTTGA	SEQ	ID NO: 27
{A_Flav}	CGG..AOGCT GGCATAATGG TOGYAAOGA	COOGTCTTGA	SEQ	ID NO: 26
{A_Ochr}	CGG..AOGCT GGCATAATGG TOGTAAOGA	COOGTCTTGA	SEQ	ID NO: 31
{A_Nige}	CGG..AOGCT GGCATAATGG TOGTAAOGA	COOGTCTTGA	SEQ	ID NO: 30
{A_Terr}	CGG..AOGCT GGCATAATGG TTGTAAOGA	COOGTCTTGA	SEQ	ID NO: 32
{A_Glau}	CGG..AOGCT GGCATAATGG TOGTAAOGA	COOGTCTTGA	SEQ	ID NO: 28
{Penici}	CGG..AOGCT GGCATAATGG TOGTAA... ..		SEQ	ID NO: 65
{C_Immi}	CGG..ATGCT GGCATAATGG TTGTAAOGG	COOGTCTTGA	SEQ	ID NO: 45
{Bipola}	AGG..ATGCT GGOGTAATGG CTGTAAOGG	COOGTCTTGA	SEQ	ID NO: 36
{Curvul}	AGG..ATGCT GGOGTAATGG CTGTAAOGG	COOGTCTTGA	SEQ	ID NO: 41
{Chryso}	AGG..ATGCT GGOGTAATGG TTGTAAOGG	COOGTCTTGA	SEQ	ID NO: 39
{Clados}	AGG..ATGCT GGOGTAATGG TOGTAAOOG	COOGTCTTGA	SEQ	ID NO: 40
{Malbra}	CGG..ATGCT GGOGTAATGG CTGTAAOGG	COOGTCTTGA	SEQ	ID NO: 62
{C_Para}	AGG..ATGTT GGCATAATGA TCTTAAGTGG	COOGTCTTGA	SEQ	ID NO: 52
{C_Trop}	AGG..ATGTT GGCATAATGA TCTTAAGTGG	COOGTCTTGA	SEQ	ID NO: 54
{C_Albi}	AGG..ATGTT GGCATAATGA TCTTAAGTGG	COOGTCTTGA	SEQ	ID NO: 42
{C_Guil}	AGG..ATGCT GGCATAATGA TCCCAAAOOG	COOGTCTTGA	SEQ	ID NO: 44
{C_Glab}	AGG..ATGCT GGCATAATGG TTATATGOOG	COOGTCTTGA	SEQ	ID NO: 43
{S_Cere}	AGG..ATGCT GGCATAATGG TTATATGOOG	COOGTCTTGA	SEQ	ID NO: 73
{C_Kefy}	AGG..ATGCT GGOGTAATGG TTAAATGOOG	COOGTCTTGA	SEQ	ID NO: 46
{Geotri}	AGG..AOGCT GGCATAATGA TTCTATAOOG	COOGTCTTGA	SEQ	ID NO: 60
{C_Iusi}	AGG..AOGCT GGOGTAATGG TTGCAAGOOG	COOGTCTTGA	SEQ	ID NO: 49
{C_Krus}	CGG..ATGCT GGCAGAAOOG OGCAACOOG	COOGTCTTGA	SEQ	ID NO: 47
{Blasch}	GGGCAOOGT CTCTGAGOC TGCTTTGCA	COOGTCTTGA	SEQ	ID NO: 37
{P_Braz}	CGG..ACGCT GGCTTAATGG TCGTAAGCGA	CCCGTCTTG-	SEQ	ID NO: 77
{Pc_Hum}	CTAAGCGACC CGTCTTGAAA CACGGACCAA	GGAGTCTAAT	SEQ	ID NO: 78

## Legend to Table 3:

The multiple sequence alignment shows the sequence of 28S ribosomal RNA genes amplified with primers SEQ ID NO: 1 and SEQ ID NO: 2. 23 species specific probes (SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76) are shown underlined. Minor sequence variation among two isolate of the same organism are represented by the appropriate code (see key below). Major differences among *Rhizopus* species are depicted by including 3 separate *Rhizopus* sequences in the alignment. (The organisms in this figure are listed according to their sequence relatedness.)

## Key to symbols:

(.) gap in sequence to facilitate alignment

(R) A or G

(W) A or T

(Y) T or C

(M) A or C

(K) T or G

(S) G or C

(B) T,G or C

Acremo	<i>Acremonium</i> species
A_clav	<i>Aspergillus clavatus</i>
25 A_flav	<i>Aspergillus flavus</i>
A_fumi	<i>Aspergillus fumigatus</i>
A_glau	<i>Aspergillus glaucus</i>
A_nidu	<i>Aspergillus nidulans</i>
A_nige	<i>Aspergillus niger</i>
30 A_ochr	<i>Aspergillus ochraceus</i>
A_terr	<i>Aspergillus terreus</i>
A_ungu	<i>Aspergillus unguis</i>

	A_ustu	<i>Aspergillus ustus</i>
	Beauve	<i>Beauveria</i> species
	Bipola	<i>Bipolaris</i> species
	Blasch	<i>Blastoschizomyces</i> species
5	B_derm	<i>Blastomyces dermatitidis</i>
	Chryso	<i>Chrysosporium</i> species
	Clados	<i>Cladosporium</i> species
	Curvul	<i>Curvularia</i> species
	C_albi	<i>Candida albicans</i>
10	C_glab	<i>Candida glabrata</i>
	C_guil	<i>Candida guilliermondii</i>
	C_immi	<i>Coccidioides immitis</i>
	C_kefy	<i>Candida kefyi</i>
	C_krus	<i>Candida krusei</i>
15	C_laur	<i>Cryptococcus laurentii</i>
	C_lusi	<i>Candida lusitaniae</i>
	C_neob	<i>Cryptococcus neoformans</i> var <i>gattii</i> serotype B
	C_neof	<i>Cryptococcus neoformans</i> serotype A
	C_para	<i>Candida parapsilosis</i>
20	C_terr	<i>Cryptococcus terreus</i>
	C_trop	<i>Candida tropicalis</i>
	Fusari	<i>Fusarium</i> species
	F_caps	<i>Filobasidium capsuligenum</i>
	F_neoc	<i>Filobasidiella</i> ( <i>Cryptococcus</i> ) <i>neoformans</i> var <i>bacillispora</i>
25		serotype C
	F_neod	<i>Filobasidiella</i> ( <i>Cryptococcus</i> ) <i>neoformans</i> var <i>neoformans</i> serotype D
	F_unig	<i>Filobasidium uniguttulatum</i>
	Geotri	<i>Geotrichum</i> species
	H_caps	<i>Histoplasma capsulatum</i>
30	Malbra	<i>Malbranchea</i> species
	Mucor_	<i>Mucor</i> species
	Paecil	<i>Paecilomyces</i> species

Penici	<i>Penicillium</i> species
P_boyd	<i>Pseudallescheria boydii</i>
P_braz	<i>Paracoccidioides brasiliensis</i>
Pc-Hum	<i>Pneumocystis carinii</i> , human isolate
5 Rhizo1	<i>Rhizopus</i> species isolate #1
Rhizo2	<i>Rhizopus</i> species isolate #2
Rhizo3	<i>Rhizopus</i> species isolate #3
Sporot	<i>Sporothrix schenckii</i>
S_brev	<i>Scopulariopsis brevicaulis</i>
10 S_brum	<i>Scopulariopsis brumptii</i>
S_cere	<i>Saccharomyces cerevisiae</i>
T_beig	<i>Trichosporon beigelii</i>
	<i>Pneumocystis carinii</i>

15 Further variations of the invention that utilize any of the named sequences will be apparent to those with ordinary skill in the art. The following examples illustrate various aspects of the invention but are not intended to limit its usefulness.

EXAMPLE 1. Testing probes SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and  
20 SEQ ID NO:76 for hybridization specificity.

Probes listed in SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76 were tested for specificity against their target organisms. Probe SEQ ID NO: 5 for *Candida albicans* was the first one tested against a panel of fungi taken  
25 from the Mayo Clinic collection. 28S rDNA from *Acremonium* sp., *Aspergillus clavatus*, *Aspergillus flavus*, *Aspergillus fumigatus*, *Aspergillus glaucus*, *Aspergillus nidulans*, *Aspergillus niger*, *Aspergillus ochraceus*, *Aspergillus terreus*, *Aspergillus unguis*, *Aspergillus ustus*, *Aspergillus* sp., *Beauveria* sp., *Bipolaris* sp., *Blastomyces dermatitidis*, *Candida albicans*, *Candida glabrata*, *Candida guilliermondii*, *Candida*  
30 *kefyr*, *Candida krusei*, *Candida lusitanae*, *Candida parapsilosis*, *Candida tropicalis*, *Chrysosporium* sp., *Cladosporium* sp., *Coccidioides immitis*, *Cryptococcus neoformans* serotype A, *Curvularia* sp., *Fusarium* sp., *Geotrichum* sp., *Histoplasma*

*capsulatum*, *Mucor* sp., *Penicillium* sp., *Pseudallescheria boydii*, *Rhizopus* sp., *Saccharomyces cerevisiae*, *Scopulariopsis brevicaulis*, *Sporothrix schenckii* and *Trichosporon beigelii* was amplified in a polymerase chain reaction using oligonucleotide probes SEQ ID NO: 1 and SEQ ID NO: 2. All PCR amplifications were carried out as hot-start reactions in a 50 ul reaction volume using Perkin-Elmer (Norwalk, CT) 0.5 ml thin-wall polypropylene tubes and a Perkin-Elmer thermal cycler. Reagents added to the tube initially were 2.5 ul of 10X PCR buffer (100 mM tris pH 8.3, 500 mM KCl, 15 mM MgCl<sub>2</sub>), 5.0 ul of 50% glycerol/1 mM cresol red, 8.0 ul of dNTP mix (1.25 mM each of dATP, dGTP, dTTP and dCTP), 11 picomoles of each nucleic acid primer and sterile water to make up a volume of 25 ul. A wax bead (Ampliwax Gem=100, Perkin-Elmer) was added and the tubes heated to 77°C for 30 seconds and cooled to room temperature to form a wax barrier. 2.5 ul of 10X PCR buffer, 5.0 ul of 50% glycerol/1 mM cresol red, 0.2 ul Taq polymerase (AmpliTaq 5U/ul, Perkin-Elmer) and 15.3 ul of sterile water was added to the tube along with 2.0 ul of DNA from the fungal whole cell boiled lysate described above. 50 cycles of thermal cycling was carried out at 94°C - 30 sec, 50°C - 1 min, 72°C - 2 min. Five microliters of polymerase chain reaction mix from each sample was run on a 5% polyacrylamide gel to visually confirm the successful amplification of 28S rDNA from each fungus listed above. 40 ul of the remaining amplified 28S rDNA was denatured in 1 N NaOH, and half of this denatured rDNA was slot blotted on to a positively charged polysulphone based membrane equilibrated in 0.5 N NaOH. The membrane was air dried for 15 minutes and baked in a vacuum oven at 80°C for 30 minutes. Amplified rDNA from each species was now bound and immobilized at a separate spot on the membrane. The free binding sites on the membrane were blocked by incubating the membrane for 3 hours at 40°C in hybridization buffer (100 ml of hybridization buffer was made using 1g non-fat milk powder, 6g NaH<sub>2</sub>PO<sub>4</sub>, 7g SDS, 200 ul 0.5M EDTA and adjusted to pH 7.2 with NaOH). The specific probe for *Candida albicans* (SEQ ID NO: 5) was end-labeled with radioactive phosphorus using <sup>32</sup>P ATP and T4 polynucleotide kinase. 50 picomoles of this probe was added to 70 milliliters of hybridization buffer and the membrane was probed at 40°C overnight. The membrane was washed in hybridization buffer at 40°C for 15 minutes followed by a wash in 2X SSC at 40°C for 15 minutes. The membrane was then exposed on x-

ray film for at least 1 hour. The oligonucleotide probe SEQ ID NO: 5 only hybridized to amplified 28S rDNA from *Candida albicans* (see Table 4) Under these hybridization conditions, probe SEQ ID NO: 5 is extremely specific for *Candida albicans*. The sequence of oligonucleotide probe SEQ ID NO: 5 differs from the sequences of other species of *Candida* by as few as 1 or 2 bases, but these mismatches are sufficient to prevent stable hybrids from forming with the other *Candida* species.

Probes SEQ ID NO: 3 to SEQ ID NO: 23 were tested for specificity, as described above for the *Candida albicans* probe SEQ ID NO: 5, against the same panel of fungi listed in the preceding paragraph. The positively charged polysulphone based membrane probed with *Candida albicans* probe SEQ ID NO: 5 was washed in 0.5 N NaOH at 40°C for 10 minutes to remove all bound *Candida albicans* probe. The membrane was sequentially probed with all probes listed in SEQ ID NO: 3 to SEQ ID NO: 23. For each subsequently tested probe, the membrane was blocked for at least 30 minutes, probe hybridization was carried out at 40-42°C for at least 3 hours, and post-hybridization washes were done in 2X SSC for 20 minutes. The membrane was stripped between probings by washing in 0.5 to 1.0 N NaOH at 40-42°C. Results are listed in Tables 4 to 7.

Probes SEQ ID NO:75 and SEQ ID NO:76 were tested for specificity against a modified panel of fungi, using the procedure just described. The specific organisms tested and the results obtained are listed in Table 8.

As shown in Tables 4 to 8, each probe listed in SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76 specifically hybridizes to only one target fungal 28S nucleic acid sequence. This specificity is essential for identifying a given species of fungus in clinical specimens containing mixed fungal organisms with a high level of reliability. The organisms listed in these Tables represent a majority of organisms that are commonly isolated from clinical samples. While we have developed 23 species specific probes (SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO: 76) that identify a total of 21 individual organisms, the additional organisms listed in the test panel were used to ensure that our probes did not have any cross-reactivity with other fungi likely to be present in a clinical specimen. The ability to accurately and reliably diagnose, and identify to a species level, this large a number of pathogens is unmatched by any other report. The fact that



we can achieve this by probing DNA amplified by a single pair of "Universal" probes (SEQ ID NO: 1 and SEQ ID NO: 2) is highly advantageous as it saves time, money and effort by providing the ability to test a single amplified target with 23 different probes (SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76).

5           A GenBank search was carried out with all probes listed in SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76 in order to determine whether similar gene sequences were present in the database. 28S sequences for *Candida albicans* *Cryptococcus neoformans* are already present in GenBank, and as expected, the probes for *Candida albicans* and *Cryptococcus neoformans* correctly identified the  
10 28S sequences from these two organisms. Ten other probes also matched DNA sequences from a variety of genes not related to the 28S gene (Table 9). This was expected because short stretches of sequence identity can often be found for any query sequence in unrelated genes from the same or a different organism. This observation is known to those versed in this art. In all cases, sequences that matched a probe  
15 sequence were not located within the 28S rRNA genes. Our probes are used to analyze 28S DNA that has been previously amplified in a polymerase chain reaction with our probes SEQ ID NO: 1 and SEQ ID NO: 2. Under stringent conditions, these two probes only amplify DNA from fungal 28S rRNA genes. Therefore no amplified DNA from the non-28S genes listed in Table 9 will be available for the hybridization  
20 of probes SEQ ID NO: 3 to SEQ ID NO: 23. The presence of related sequences in non-28S, unamplified genes will not be detected and will, thus, not have any effect on the sensitivity or the specificity of our detection and identification strategy.

Table 4:

Detection of species specific 28S sequence with probes SEQ ID NO: 3 to SEQ ID NO: 8

FUNGUS	SEQ ID: 3	SEQ ID: 4	SEQ ID: 5	SEQ ID: 6	SEQ ID: 7	SEQ ID: 8
<i>Acremonium</i> sp.	-	-	-	-	-	-
<i>Aspergillus clavatus</i>	-	-	-	-	-	-
<i>Aspergillus flavus</i>	-	-	-	-	-	-
<i>Aspergillus fumigatus</i>	+	-	-	-	-	-
<i>Aspergillus glaucus</i>	-	-	-	-	-	-
<i>Aspergillus nidulans</i>	-	-	-	-	-	-
<i>Aspergillus niger</i>	-	-	-	-	-	-
<i>Aspergillus ochraceus</i>	-	-	-	-	-	-
<i>Aspergillus terreus</i>	-	-	-	-	-	-
<i>Aspergillus unguis</i>	-	-	-	-	-	-
<i>Aspergillus ustus</i>	-	-	-	-	-	-
<i>Aspergillus</i> sp.	-	-	-	-	-	-
<i>Beauvaria</i> sp.	-	-	-	-	-	-
<i>Bipolaris</i> sp.	-	-	-	-	-	-
<i>Blastomyces dermatitidis</i>	-	+	-	-	-	-
<i>Candida albicans</i>	-	-	+	-	-	-
<i>Candida glabrata</i>	-	-	-	-	-	-
<i>Candida guilliermondii</i>	-	-	-	-	-	-
<i>Candida kefyr</i>	-	-	-	-	-	-
<i>Candida krusei</i>	-	-	-	-	-	-
<i>Candida lusitanae</i>	-	-	-	-	-	-
<i>Candida parapsilosis</i>	-	-	-	-	-	-
<i>Candida tropicalis</i>	-	-	-	-	-	-

<i>Chrysosporium</i> sp.	-	-	-	-	-	-
<i>Cladosporium</i> sp.	-	-	-	-	-	-
<i>Coccidioides immitis</i>	-	-	-	+	-	-
<i>Cryptococcus neoformans</i>	-	-	-	-	+	+
<i>Curvularia</i> sp.	-	-	-	-	-	-
<i>Fusarium</i> sp.	-	-	-	-	-	-
<i>Geotrichum</i> sp.	-	-	-	-	-	-
<i>Histoplasma capsulatum</i>	-	-	-	-	-	-
<i>Mucor</i> sp.	-	-	-	-	-	-
<i>Penicillium</i> sp.	-	-	-	-	-	-
<i>Pseudallescheria boydii</i>	-	-	-	-	-	-
<i>P. brasiliensis</i>	-	-	-	-	-	-
<i>Pneumocystis carinii</i>	-	-	-	-	-	-
<i>Rhizopus</i> sp.	-	-	-	-	-	-
<i>Saccharomyces cerevisiae</i>	-	-	-	-	-	-
<i>Scopulariopsis brevicaulis</i>	-	-	-	-	-	-
<i>Sporothrix schenckii</i>	-	-	-	-	-	-
<i>Trichosporon beigelii</i>	-	-	-	-	-	-

+	Positive
-	Negative after 20 minute wash in 2X SSC

Table 5:

Detection of species specific 28S sequence with probes SEQ ID NO: 9 to SEQ ID

NO: 14

FUNGUS	SEQ ID: 9	SEQ ID: 10	SEQ ID: 11	SEQ ID: 12	SEQ ID: 13	SEQ ID: 14
<i>Acremonium</i> sp.	-	-	-	-	-	-
<i>Aspergillus clavatus</i>	-	-	-	-	-	-
<i>Aspergillus flavus</i>	-	-	-	-	-	-
<i>Aspergillus fumigatus</i>	-	-	-	-	-	-
<i>Aspergillus glaucus</i>	-	+	-	-	-	-
<i>Aspergillus nidulans</i>	-	-	-	-	-	-
<i>Aspergillus niger</i>	-	-	+	-	-	-
<i>Aspergillus ochraceus</i>	-	-	-	-	-	-
<i>Aspergillus terreus</i>	-	-	-	+	-	-
<i>Aspergillus unguis</i>	-	-	-	-	-	-
<i>Aspergillus ustus</i>	-	-	-	-	-	-
<i>Aspergillus</i> sp.	-	-	-	-	-	-
<i>Beauvaria</i> sp.	-	-	-	-	-	-
<i>Bipolaris</i> sp.	-	-	-	-	-	-
<i>Blastomyces dermatitidis</i>	-	-	-	-	-	-
<i>Candida albicans</i>	-	-	-	-	-	-
<i>Candida glabrata</i>	-	-	-	-	+	-
<i>Candida guilliermondii</i>	-	-	-	-	-	+
<i>Candida kefyr</i>	-	-	-	-	-	-
<i>Candida krusei</i>	-	-	-	-	-	-
<i>Candida lusitanae</i>	-	-	-	-	-	-
<i>Candida parapsilosis</i>	-	-	-	-	-	-
<i>Candida tropicalis</i>	-	-	-	-	-	-

<i>Chrysosporium</i> sp.	-	-	-	-	-	-
<i>Cladosporium</i> sp.	-	-	-	-	-	-
<i>Coccidioides immitis</i>	-	-	-	-	-	-
<i>Cryptococcus neoformans</i>	-	-	-	-	-	-
<i>Curvularia</i> sp.	-	-	-	-	-	-
<i>Fusarium</i> sp.	-	-	-	-	-	-
<i>Geotrichum</i> sp.	-	-	-	-	-	-
<i>Histoplasma capsulatum</i>	+	-	-	-	-	-
<i>Mucor</i> sp.	-	-	-	-	-	-
<i>Penicillium</i> sp.	-	-	-	-	-	-
<i>P. brasiliensis</i>	-	-	-	-	-	-
<i>Pneumocystis carinii</i>	-	-	-	-	-	-
<i>Pseudallescheria boydii</i>	-	-	-	-	-	-
<i>Rhizopus</i> sp.	-	-	-	-	-	-
<i>Saccharomyces cerevisiae</i>	-	-	-	-	-	-
<i>Scopulariopsis brevicaulis</i>	-	-	-	-	-	-
<i>Sporothrix schenckii</i>	-	-	-	-	-	-
<i>Trichosporon beigeli</i>	-	-	-	-	-	-

+	Positive
-	Negative after 20 minute wash in 2X SSC

Table 6:

Detection of species specific 28S sequence with probes SEQ ID NO: 15 to SEQ ID

NO: 20

FUNGUS	SEQ ID: 15	SEQ ID: 16	SEQ ID: 17	SEQ ID: 18	SEQ ID: 19	SEQ ID: 20
<i>Acremonium</i> sp.	-	-	-	-	-	-
<i>Aspergillus clavatus</i>	-	-	-	-	-	-
<i>Aspergillus flavus</i>	-	-	-	-	-	-
<i>Aspergillus fumigatus</i>	-	-	-	-	-	-
<i>Aspergillus glaucus</i>	-	-	-	-	-	-
<i>Aspergillus nidulans</i>	-	-	-	-	-	-
<i>Aspergillus niger</i>	-	-	-	-	-	-
<i>Aspergillus ochraceus</i>	-	-	-	-	-	-
<i>Aspergillus terreus</i>	-	-	-	-	-	-
<i>Aspergillus unguis</i>	-	-	-	-	-	-
<i>Aspergillus ustus</i>	-	-	-	-	-	-
<i>Aspergillus</i> sp.	-	-	-	-	-	-
<i>Beauveria</i> sp.	-	-	-	-	-	-
<i>Bipolaris</i> sp.	-	-	-	-	-	-
<i>Blastomyces dermatitidis</i>	-	-	-	-	-	-
<i>Candida albicans</i>	-	-	-	-	-	-
<i>Candida glabrata</i>	-	-	-	-	-	-
<i>Candida guilliermondii</i>	-	-	-	-	-	-
<i>Candida kefyr</i>	+	-	-	-	-	-
<i>Candida krusei</i>	-	+	-	-	-	-
<i>Candida lusitanae</i>	-	-	+	-	-	-
<i>Candida parapsilosis</i>	-	-	-	+	-	-
<i>Candida tropicalis</i>	-	-	-	-	+	-

<i>Chrysosporium</i> sp.	-	-	-	-	-	-
<i>Cladosporium</i> sp.	-	-	-	-	-	-
<i>Coccidioides immitis</i>	-	-	-	-	-	-
<i>Cryptococcus neoformans</i>	-	-	-	-	-	-
<i>Curvularia</i> sp.	-	-	-	-	-	-
<i>Fusarium</i> sp.	-	-	-	-	-	-
<i>Geotrichum</i> sp.	-	-	-	-	-	-
<i>Histoplasma capsulatum</i>	-	-	-	-	-	-
<i>Mucor</i> sp.	-	-	-	-	-	-
<i>Penicillium</i> sp.	-	-	-	-	-	-
<i>P. brasiliensis</i>	-	-	-	-	-	-
<i>Pneumocystis carinii</i>	-	-	-	-	-	-
<i>Pseudallescheria boydii</i>	-	-	-	-	-	+
<i>Rhizopus</i> sp.	-	-	-	-	-	-
<i>Saccharomyces cerevisiae</i>	-	-	-	-	-	-
<i>Scopulariopsis brevicaulis</i>	-	-	-	-	-	-
<i>Sporothrix schenckii</i>	-	-	-	-	-	-
<i>Trichosporon beigelii</i>	-	-	-	-	-	-

+	Positive
-	Negative after 20 minute wash in 2X SSC

Table 7:

Detection of species specific 28S sequence with probes SEQ ID NO: 21 to SEQ ID NO: 23

FUNGUS	SEQ ID: 21	SEQ ID: 22	SEQ ID: 23
<i>Acremonium</i> sp.	-	-	-
<i>Aspergillus clavatus</i>	-	-	-
<i>Aspergillus flavus</i>	+	-	-
<i>Aspergillus fumigatus</i>	-	-	-
<i>Aspergillus glaucus</i>	-	-	-
<i>Aspergillus nidulans</i>	-	-	-
<i>Aspergillus niger</i>	-	-	-
<i>Aspergillus ochraceus</i>	-	-	-
<i>Aspergillus terreus</i>	-	-	-
<i>Aspergillus unguis</i>	-	-	-
<i>Aspergillus ustus</i>	-	-	-
<i>Aspergillus</i> sp.	-	-	-
<i>Beauveria</i> sp.	-	-	-
<i>Bipolaris</i> sp.	-	-	-
<i>Blastomyces dermatitidis</i>	-	-	-
<i>Candida albicans</i>	-	-	-
<i>Candida glabrata</i>	-	-	-
<i>Candida guilliermondii</i>	-	-	-
<i>Candida kefyr</i>	-	-	-
<i>Candida krusei</i>	-	-	-
<i>Candida lusitanae</i>	-	-	-
<i>Candida parapsilosis</i>	-	-	-
<i>Candida tropicalis</i>	-	-	-
<i>Chrysosporium</i> sp.	-	-	-
<i>Cladosporium</i> sp.	-	-	-



<i>Coccidioides immitis</i>	-	-	-
<i>Cryptococcus neoformans</i>	-	-	-
<i>Curvularia</i> sp.	-	-	-
<i>Fusarium</i> sp.	-	-	-
<i>Geotrichum</i> sp.	-	-	-
<i>Histoplasma capsulatum</i>	-	-	-
<i>Mucor</i> sp.	-	-	-
<i>Penicillium</i> sp.	-	-	-
<i>Paracoccidioides brasiliensis</i>	-	-	-
<i>Pneumocystis carinii</i>	-	-	-
<i>Pseudallescheria boydii</i>	-	-	-
<i>Rhizopus</i> sp.	-	-	-
<i>Saccharomyces cerevisiae</i>	-	-	-
<i>Scopulariopsis brevicaulis</i>	-	-	-
<i>Sporothrix schenckii</i>	-	+	+
<i>Trichosporon beigelii</i>	-	-	-

+	Positive
-	Negative after 20 minute wash in 2X SSC

Table 8:

Detection of species specific 28S sequence with probes SEQ ID NO: 75 to SEQ ID NO: 76

FUNGUS	SEQ ID NO: 75	SEQ ID NO: 76
<i>Acremonium sp.</i>	-	-
<i>Aspergillus clavatus</i>	-	-
<i>Aspergillus flavus</i>	-	-
<i>Aspergillus fumigatus</i>	-	-
<i>Aspergillus glaucus</i>	-	-
<i>Aspergillus nidulans</i>	-	-
<i>Aspergillus niger</i>	-	-
<i>Aspergillus ochraceus</i>	-	-
<i>Aspergillus terreus</i>	-	-
<i>Aspergillus unguis</i>	-	-
<i>Aspergillus ustus</i>	-	-
<i>Beauveria sp.</i>	-	-
<i>Bipolaris sp.</i>	-	-
<i>Blastomyces dermatitidis</i>	-	-
<i>Blasotschizomyces capitatus</i>	-	-
<i>Candida albicans</i>	-	-
<i>Candida glabrata</i>	-	-
<i>Candida guilliermondii</i>	-	-
<i>Candida krusei</i>	-	-
<i>Candida lusitaniae</i>	-	-
<i>Candida parapsilosis</i>	-	-
<i>Candida tropicalis</i>	-	-
<i>Chrysosporium sp.</i>	-	-
<i>Cladosporium sp.</i>	-	-
<i>Coccidioides immitis</i>	-	-

<i>Cryptococcus laurentii</i>	-	-
<i>Cryptococcus neoformans</i>	-	-
var <i>neoformans</i> (sero A)	-	-
var <i>neoformans</i> (sero D)	-	-
var <i>gatti</i> (sero B)	-	-
<i>Cryptococcus terreus</i>	-	-
<i>Curvularia</i> sp.	-	-
<i>Filobasidium capsuligenum</i>	-	-
<i>Filobasidium uniguttulatum</i>	-	-
<i>Fusarium</i> sp.	-	-
<i>Geotrichum</i> sp.	-	-
<i>Histoplasma capsulatum</i>	-	-
<i>Malbranchea filamentosum</i>	-	-
<i>Mucor</i> sp.	-	-
<i>Paecilomyces</i> sp.	-	-
<i>Paracoccidioides brasiliensis</i>		
isolate no. 135	+	-
isolate no. 262	+	-
isolate no. 265	+	-
isolate no. 927	+	-
isolate no. 9919	+	-
isolate no. 9894	+	-
<i>Prototheca</i> sp.	-	-
<i>Penicillium</i> sp.	-	-
<i>Pneumocystis carinii</i>		
isolate no. Mayo ref. cult.	-	+
isolate no. 56	-	+
isolate no. 62	-	+
isolate no. 69	-	+
isolate no. 85	-	+
<i>Pseudallesheria boydii</i>	-	-

<i>Rhizopus sp.</i>	-	-
<i>Saccharomyces cerevisiae</i>	-	-
<i>Scopulariopsis brevicaulis</i>	-	-
<i>Scopulariopsis brumptii</i>	-	-
<i>Sporothrix schenckii</i>	-	-
<i>Trichosporon beigeli</i>	-	-

Table 9:

GenBank search results listing genes from other organisms having 100% identity to probes SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76

	PROBE SEQ ID NO:	ORGANISM MATCHED	GENE MATCHED* (see note below)	ACCESSION NUMBER
<i>Aspergillus fumigatus</i>	3	-	-	-
<i>Blastomyces dermatitidis</i>	4	<i>Streptomyces verticillus</i>	bleomycin acetyl transferase	L26955
	4	<i>Giardia muris</i>	upstream of rRNA genes	X65063, S53320
	4	<i>Aspergillus nidulans</i>	uric acid-xanthine permease	X71807
	4	<i>Homo sapiens</i>	T-cell surface glycoprotein	X16996
	4	<i>Homo sapiens</i>	MIC2	M16279, M22557, J03841, M22556
<i>Candida albicans</i>	5	<i>Candida albicans</i>	28S rRNA	L28817
<i>Coccidioides immitis</i>	6	-	-	-
<i>Cryptococcus neoformans</i>	7	<i>Cryptococcus neoformans</i>	28S rRNA	L14067, L14068,
<i>Cryptococcus neoformans</i>	8	<i>Cryptococcus neoformans</i>	28S rRNA	L14067, L14068, L20964
	8	<i>Escherichia coli</i>	0111 cld	Z17241
<i>Histoplasma capsulatum</i>	9	-	-	-

<i>Aspergillus glaucus</i>	10	<i>Pseudomonas denitrificans</i>	cob genes	M62866
<i>Aspergillus niger</i>	11	-	-	-
<i>Aspergillus terreus</i>	12	Human cytomegalovirus	genome	X17403
	12	<i>Homo sapiens</i>	GABA receptor	L08485
<i>Candida glabrata</i>	13	<i>Homo sapiens</i>	Class I MHC	X03664, X03665
<i>Candida guilliermondii</i>	14	-	-	-
<i>Candida kefyr</i>	15	-	-	-
<i>Candida krusei</i>	16	<i>Pseudomonas syringae</i>	penicillin binding protein	L28837
<i>Candida lusitanae</i>	17	Chicken	AK1	D00251
	17	Mouse	IL10	M84340
<i>Candida parapsilosis</i>	18	<i>Polytomella agilis</i>	beta-2 tubulin	M33373
	18	Tobacco chloroplast	genome	Z00044, S54304
	18	<i>Aedes aegypti</i>	amylase	L03640
	18	<i>Homo sapiens</i>	chromosome 13q14	L14473
<i>Candida tropicalis</i>	19	-	-	-
<i>Pseudallescheria boydii</i>	20	<i>Drosophila melanogaster</i>	AcTr66B	X71789
		Cow	actin 2	D12816
<i>Aspergillus flavus</i>	21	-	-	-
<i>Sporothrix schenckii</i>	22	-	-	-
<i>Sporothrix schenckii</i>	23	Sulfate reducing bacteria	FMN binding protein	D21804
	23	Equine herpesvirus 1	genome	M86664
<i>Paracoccidioides brasiliensis</i>	75	European rabbit	BI-1 calcium channel gene	X57476

<i>Pneumocystis carinii</i>	76	rice	heat shock protein	D1775
	76	<i>C. elegans</i>	repetitive DNA	X61259
	76	<i>Synechocystis</i> sp.	genome	D64000
	76	<i>C. elegans</i>	various cosmids	Z81063, U59749, Z69662, Z49936, U88182, U80931, Z47358, Z81142, Z79602, Z82084, Z70269

\* Note: As discussed earlier in this document, the presence of sequences similar to probes SEQ ID NO:3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76 in genes not related to 28S does not have any effect on the specificity or sensitivity of our diagnostic strategy. Our species specific probes are used to analyze 28S DNA that has been previously amplified in a polymerase chain reaction with our probes SEQ ID NO: 1 and SEQ ID NO:2. These two probes will not amplify DNA from any gene other than 28S in column #4 (GENE MATCHED), and therefore no amplified DNA from these non-28S genes will be available for the hybridization of probes SEQ ID NO: 3 to SEQ ID NO: 23, SEQ ID NO:75 and SEQ ID NO:76.

EXAMPLE 2. Use of method in example 1 to test clinical specimens for specific fungal organisms.

Clinical samples taken from the respiratory and gastrointestinal tract of healthy individuals almost always contain some fungal flora. Most of these fungi are non-

pathogenic, but may give false positives on traditional immunochemical diagnostic tests for pathogenic fungi.

We obtained 44 clinical specimens from diverse sources ranging from sputum and incision drainage tubes, to intervertebral disc and lung biopsies. Traditional smear and culture results showed that all 44 specimens contained at least 1 type of fungus. In order to test the efficacy of our probes, we extracted DNA from all 44 clinical samples and used probes SEQ ID NO: 1 & 2 in a polymerase chain reaction to amplify fungal 28S sequences present in these samples.

DNA was extracted from all clinical samples by our modification of the technique of Chomczynski and Sacchi which originally described the use of acid guanidinium thiocyanate-phenol-chloroform to preferentially extract RNA from cells and tissues. We replaced room temperature cell lysis by boiling lysis, and acid guanidinium thiocyanate-phenol-chloroform extraction by alkaline phenol-guanidine thiocyanate to preferentially extract DNA from cells. 1.5 ml Sarsted (Newton, North Carolina) polypropylene screw cap tubes with o-ring seals were used for the extractions. 200 ul of specimen was added to 500 ul of GPT reagent (6 M guanidine thiocyanate dissolved in 50 mM tris pH 8.3 mixed with an equal volume of phenol buffered in tris pH 8.0). This was mixed by vortexing and immediately placed in a boiling water bath for 15 minutes. The tubes were spun in a microcentrifuge for 5 seconds and 250 ul of chloroform/iso-amyl alcohol (24:1 by volume) was added and mixed by vortexing. The liquid phases were separated by centrifugation for 10 minutes and 450 ul of aqueous (upper) phase was transferred to a fresh tube. The aqueous phase was mixed with 500 ul of 100% isopropanol and placed at -20°C for at least 1 hour. At the end of this period the tubes were centrifuged for 15 minutes and the supernatant removed without disturbing the nucleic acid pellet. The pellet was washed with 500 ul of ice-cold 70% ethanol to remove traces of GPT reagent by gently inverting 2 times and then centrifuged for 5 minutes. The ethanol was removed and the pellet dried in a speed vac for 10 minutes. The pellet was resuspended in 25 ul of sterile deionized water and 5 ul was used in a 50 ul PCR amplification. The PCR was carried out as a hot-start reaction using the thermal cycling conditions for probes SEQ ID NO: 1 and SEQ ID NO: 2 described in example 1. Gel electrophoresis



showed that probes SEQ ID NO: 1 and SEQ ID NO: 2 successfully amplified DNA from all 44 specimens.

The amplified DNA from each specimen was transferred to a positively charged polysulphone based membrane. We radioactively labeled our species specific probes SEQ ID NO: 3, SEQ ID NO: 5, SEQ ID NO: 6 and SEQ ID NO: 7, and sequentially probed the membrane to test for the presence of 28S rDNA from *Aspergillus fumigatus*, *Candida albicans*, *Coccidioides immitis* and *Cryptococcus neoformans* respectively. Membrane blocking, probe hybridization and washes were done exactly as described in Example 1. The results are shown in Table 10.

No false positives were observed, indicating a specificity of 100% for these 4 probes in the clinical specimens tested. 10 out of 12 culture positive samples for *Aspergillus fumigatus*, and 11 out of 13 samples of *Candida albicans* were identified, indicating a detection sensitivity of about 85% for these two probes. Additionally, two out of two *Coccidioides immitis* and two out of two *Cryptococcus neoformans* were correctly identified (detection sensitivity of 100%). As seen by these results, the probes described in this invention can be used on a diverse variety of clinical specimens with excellent efficacy.

Table 10.

Detection of *Aspergillus fumigatus*, *Candida albicans*, *Coccidioides immitis* and *Cryptococcus neoformans* in clinical specimens using species specific probes.

Specimen type	Smear and culture results	PCR with SEQ ID: 1, SEQ ID: 2	SEQ ID: 3	SEQ ID: 5	SEQ ID: 6	SEQ ID: 7
U035 sputum	A. flavus	+	-	-	-	-
U069 pleura	A. fumigatus	+	+	-	-	-
U070 bronchial wash	A. flavus	+	-	-	-	-

M019 bronchial wash	A. fumigatus	+	+	-	-	-
M020 sputum	mixed fungal flora	+	-	+	-	-
X35254 sputum	C. albicans	+	-	+	-	-
M20910 sputum	A. fumigatus	+	+	-	-	-
M055 sputum	C. albicans	+	-	+	-	-
M056 abdominal	mixed fungal flora	+	-	-	-	-
M057 drainage tube	C. albicans	+	-	(-)	-	-
M059 ind. sputum	C. albicans	+	-	+	-	-
M060 ind. sputum	mixed fungal flora	+	-	-	-	-
M083 bronchial wash	C. albicans	+	-	+	-	-
M084 sputum	A. fumigatus	+	(-)	-	-	-
M085 throat	C. albicans	+	-	(-)	-	-
A001 sputum	A. fumigatus	+	(-)	-	-	-
A002 leg	Blastomyces	+	-	-	-	-
A003 leg	Blastomyces	+	-	-	-	-
A005 disc	A. fumigatus	+	+	-	-	-
A037 disc	A. fumigatus	+	+	-	-	-
A039 trachea	C. albicans	+	-	+	-	-
A040 trachea	C. albicans	+	-	+	-	-
A102 empyema	A. fumigatus	+	+	-	-	-
Y004 sputum	C. albicans	+	-	+	-	-
Y016 induced sputum	Coccidioides	+	-	-	+	-
Y028 sputum	Coccidioides	+	-	-	+	-
J003 chest	Aspergillus sp.	+	-	-	-	-

J045 bronchial wash	C. albicans	+	-	+	-	-
J046 ethmoid	yeast	+	-	-	-	-
J047 chest	A. fumigatus	+	+	-	-	-
J048 sputum	C. albicans	+	-	+	-	-
J073 lung	Aspergillus sp.	+	-	-	-	-
J074 lung	A. fumigatus	+	+	-	-	-
U017 lip	A. fumigatus	+	+	-	-	-
U033 sputum	mixed fungal flora	+	-	-	-	-
U071 sputum	C. albicans	+	-	+	-	-
U072 BA lavage	Sporothrix	+	-	-	-	-
U073 knee	Histoplasma	+	-	-	-	-
U074 mandible	Cryptococcus	+	-	-	-	+
U075 CSF	Cryptococcus	+	-	-	-	+
U076 knee	Histoplasma	+	-	-	-	-
U077 soft tissue	Histoplasma	+	-	-	-	-
U051 buccal	A. fumigatus	+	+	-	-	-
Y055 sputum	mixed fungal flora	+	-	-	-	-
+ Positive    - Negative    (-) Missed						

### EXAMPLE 3. DNA sequence based identification of unknown fungal organisms.

Another utility of our probes is in the rapid DNA sequence based identification of a pure culture of fungus. Probes SEQ ID NO: 1 and SEQ ID NO: 2 are used in a polymerase chain reaction to amplify 28S rDNA from an unknown fungus. Probes SEQ ID NO: 1 or SEQ ID NO: 2 are then used as sequencing primers to obtain DNA sequence from this amplified 28S DNA belonging to the unknown fungus. This DNA sequence is compared to the fungal 28S DNA sequences in our database, and a sequence match at, or overlapping any one of the probe sequences in SEQ ID NO: 3 to

SEQ ID NO: 78 will confirm the identity of the fungus. This technique cannot be used directly on clinical samples, as these usually contain DNA from more than one fungus, and the DNA sequence generated will consist of overlapping sequences of several organisms. This technique has utility in rapidly and reliably identifying colonies of a single fungus on culture plates, clinical specimens, food, pharmaceutical, environmental or other samples containing only one species of fungus.

#### EXAMPLE 4. Capture and identification of target DNA or RNA

All primers and probes described in this invention disclosure may be labeled with any detectable reporter or signal moiety including, but not limited to radioisotopes, enzymes, antigens, antibodies, chemiluminescent reagents and fluorescent chemicals. Additionally, these probes may be modified without changing the substance of their purpose by terminal addition of nucleotides designed to incorporate restriction sites or other useful sequences. These probes may also be modified by the addition of a capture moiety (including, but not limited to paramagnetic particles, biotin, fluorescein, dioxigenin, antigens, antibodies) or attached to the walls of microtiter trays to assist in the solid phase capture and purification of these probes and any DNA or RNA hybridized to these probes. Fluorescein may be used as a signal moiety as well as a capture moiety, the latter by interacting with an anti-fluorescein antibody.

A typical utility of these modifications would be as follows. Primers SEQ ID NO: 1 and SEQ ID NO: 2 would be utilized to amplify 28S rDNA from a sample, if present, as described previously. Primers would be modified so as to contain a biotin moiety at their 5' ends. A streptavidin solid phase, such as a paramagnetic particle, would be used to separate PCR products, if present, from the reaction mixture. The amplified target may be subsequently hybridized to a third probe ((SEQ ID NO: 3) to (SEQ ID NO: 78) or their complements) attached to a detectable moiety to determine which species of fungus is present in the given sample. Multiple probes, each labeled with a different detectable moiety may be used at one time to analyze the amplified target.

Alternatively, Primers SEQ ID NO: 1 and SEQ ID NO: 2 would be utilized to amplify 28S rDNA from a sample, if present, as above. In a separate reaction, individually, either SEQ ID NO: 1 or SEQ ID NO: 2 would be modified by attachment to a solid phase capture moiety, such as a paramagnetic particle, and SEQ ID NO: 3 to SEQ ID NO: 78 (or their complements) would be modified by addition of a detectable moiety. Alternately, in the amplicon, any sequences delimited by SEQ ID NO: 1 and SEQ ID NO: 2, including but not limited to SEQ ID NO: 3 to SEQ ID NO: 78, may be used in the design of a capture probe. One of the probes attached to a solid phase (SEQ ID NO: 1 and SEQ ID NO: 2) or any other appropriately designed sequences and one of the probes modified by attachment to a detectable moiety (SEQ ID NO: 3 to SEQ ID NO: 78 or their complements) would be hybridized together, in solution, to products of the PCR, if they had been generated. The hybrids, if present, would be captured from the solution, and analyzed by a method appropriate to the detection moiety. Detection of the hybridized probe would indicate which species of fungus was present in the given sample. Multiple probes, each labeled with a different detectable moiety may be used at one time to analyze the amplified target.

#### EXAMPLE 5. Species-specific amplification of fungal DNA

Another utility of the probes described in this invention is their usage as primers in the direct detection of a specific fungal species by virtue of a nucleic acid amplification reaction. In this embodiment, one primer is a universal one, such as (SEQ ID NO:1) or (SEQ ID NO:2), and the other is a species-specific primer selected from the group consisting of (SEQ ID NO:3) to (SEQ ID NO: 23), (SEQ ID NO:75) and (SEQ ID NO:76) or the complements thereof. One variation of this approach is the substitution of (SEQ ID NO:1) or (SEQ ID NO:2) with any functional sequence located in proximity to the species-specific primer. Another variation of this approach is the selection of any appropriate species specific primer pair from SEQ ID NO: 24 to SEQ ID NO: 74, SEQ ID NO:77 and SEQ ID NO:78.

## SEQUENCE LISTING

## (1) GENERAL INFORMATION:

5

## (i) APPLICANT:

(A) Sandhu, Gurpreet S.

(B) Kline, Bruce C.

10

## (ii) TITLE OF INVENTION:

Nucleic Acid Probes for the Detection and Identification of Fungi

## (iii) NUMBER OF SEQUENCES: 80

15

## (iv) CORRESPONDENCE ADDRESS:

(A) ADDRESSEE: Chiron Diagnostics Corporation

(B) STREET: 63 North Street

(C) CITY: Medfield

(D) STATE: Massachusetts

20

(E) COUNTRY: USA

(F) ZIP: 02052

## (v) COMPUTER READABLE FORM:

(A) MEDIUM TYPE: Diskette 3.5 inch, 1.44 Mb storage

25

(B) COMPUTER: IBM PS/2

(C) OPERATING SYSTEM: MS-DOS 6.2

(D) SOFTWARE: Word 6.0

## (vi) CURRENT APPLICATION DATA:

30

(A) APPLICATION NUMBER:

(B) FILING DATE:

(C) CLASSIFICATION:

## (vii) PRIOR APPLICATION DATA:

(A) APPLICATION NUMBER: 08/373,127

(B) FILING DATE: 13-JAN-1995

(C) CLASSIFICATION

## (vii) PRIOR APPLICATION DATA:

(A) APPLICATION NUMBER: 08/435,684

(B) FILING DATE: 05-MAY-1995

(C) CLASSIFICATION

## (viii) ATTORNEY INFORMATION:

(A) NAME: Morgenstern, Arthur S.

(B) REGISTRATION NUMBER: 28,244

(C) DOCKET NUMBER: CCD-180CIPII

## (ix) TELECOMMUNICATION INFORMATION:

(A) TELEPHONE: 508 359-3836

(B) TELEFAX: 508 359-3885

## (2) INFORMATION FOR SEQ ID NO 1:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 20

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

## (ii) MOLECULE TYPE: Nucleic acid probe for fungal organisms

## (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 1:

5 GTGAAATTGT TGAAAGGGAA 20

(3) INFORMATION FOR SEQ ID NO 2:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 18

10 (B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for fungal organisms

15

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

20

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 2:

GACTCCTTGG TCCGTGTT 18

(4) INFORMATION FOR SEQ ID NO 3:

25

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

30

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for *Aspergillus fumigatus*



(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

5

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 3:

CTCGGAATGT ATCA

14

10 (5) INFORMATION FOR SEQ ID NO 4:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 13

(B) TYPE: nucleic acid

15

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for Blastomyces dermatitidis

20

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 4:

25

ACTCCCCAC GGG

13

(6) INFORMATION FOR SEQ ID NO 5:

30

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for *Candida albicans*

5

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

10

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 5:

CCTCTGACGA TGCT

14

(7) INFORMATION FOR SEQ ID NO 6:

15

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

20

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for *Coccidioides immitis*

(iii) HYPOTHETICAL: No

25

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 6:

30

TCTGGCGGTT GGTT

14

(8) INFORMATION FOR SEQ ID NO 7:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 14
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

5

(ii) MOLECULE TYPE: Nucleic acid probe for *Cryptococcus neoformans*

(iii) HYPOTHETICAL: No

10

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 7:

CTCCTGTCGC ATAC

14

15

(9) INFORMATION FOR SEQ ID NO 8:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 14
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

20

(ii) MOLECULE TYPE: Nucleic acid probe for *Cryptococcus neoformans*

25

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 8:

30

AGTTCTGATC GGTG

14

## (10) INFORMATION FOR SEQ ID NO 9:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for Histoplasma capsulatum

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 9:

CAATCCCCCG CGGC

14

## (11) INFORMATION FOR SEQ ID NO 10:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for Aspergillus glaucus

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 10:

GTGTCATGCG GCCA

14

5 (12) INFORMATION FOR SEQ ID NO 11:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

10 (C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for *Aspergillus niger*

15 (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 11:

20

CCCTGGAATG TAGT

14

(13) INFORMATION FOR SEQ ID NO 12:

(i) SEQUENCE CHARACTERISTICS:

25 (A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

30 (ii) MOLECULE TYPE: Nucleic acid probe for *Aspergillus terreus*

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 12:

5

GCTTCGGCCC GGTG 14

(14) INFORMATION FOR SEQ ID NO 13:

10

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

15

(ii) MOLECULE TYPE: Nucleic acid probe for *Candida glabrata*

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

20

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 13:

CTTGGGACTC TCGC 14

25

(15) INFORMATION FOR SEQ ID NO 14:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

30

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for *Candida guilliermondii*

(iii) HYPOTHETICAL: No

5 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 14:

ATATTTTGTG AGCC 14

10

(16) INFORMATION FOR SEQ ID NO 15:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

15

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for *Candida kefyr*

20

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

25

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 15:

TTCGGCTTTC GCTG 14

(17) INFORMATION FOR SEQ ID NO 16:

30

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

5 (ii) MOLECULE TYPE: Nucleic acid probe for *Candida krusei*

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

10

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 16:

GGGATTGCGC ACCG 14

(18) INFORMATION FOR SEQ ID NO 17:

15

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 14
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

20

(ii) MOLECULE TYPE: Nucleic acid probe for *Candida lusitanae*

(iii) HYPOTHETICAL: No

25

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 17:

30 GCCTCCATCC CTTT 14

(19) INFORMATION FOR SEQ ID NO 18:



## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

5 (C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for Candida parapsilosis

10 (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 18:

15 ATAAGTGCAA AGAA

14

(20) INFORMATION FOR SEQ ID NO 19:

## (i) SEQUENCE CHARACTERISTICS:

20 (A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

25 (ii) MOLECULE TYPE: Nucleic acid probe for Candida tropicalis

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

30

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 19:

AGAATTGCGT TGGA

14

## (21) INFORMATION FOR SEQ ID NO 20:

## 5 (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

10 (ii) MOLECULE TYPE: Nucleic acid probe for *Pseudallescheria boydii*

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

15

## (v) SEQUENCE DESCRIPTION: SEQ ID NO: 20:

GCGATGGGAA TGTG

14

## 20 (22) INFORMATION FOR SEQ ID NO 21:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

## 25 (C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for *Aspergillus flavus*

## 30 (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 21:

AGACTCGCCT CCAG 14

5 (23) INFORMATION FOR SEQ ID NO 22:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

10 (C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for *Sporothrix schenckii*

15 (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 22:

20

CGGACCACCC GGCG 14

(24) INFORMATION FOR SEQ ID NO 23:

25 (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

30

(ii) MOLECULE TYPE: Nucleic acid probe for *Sporothrix schenckii*

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 23:

5 CGGCGGCATG CCCC 14

(25) INFORMATION FOR SEQ ID NO 24:

(i) SEQUENCE CHARACTERISTICS:

10 (A) LENGTH: 208  
(B) TYPE: nucleic acid  
(C) STRANDEDNESS: single  
(D) TOPOLOGY: linear

15 (ii) MOLECULE TYPE: Acremonium species specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

20 (v) SEQUENCE DESCRIPTION: SEQ ID NO: 24:

GACCAGACTT GGGCTCGGTG AATCATCCGG CGTTCTCGCC  
25 GGTGCACTTT 50  
GCCGTCCCAG GCCAGCATCA GTTCGCGCCG GGGGATAAAG  
GTTTCGGGAA 100  
TGTAGCTCCT TCGGGAGTGT TATAGCCCGT TGCGTAATAC  
CCTGGCGTGG 150  
30 ACTGAGGTCC GCGCTCTGCA AGGATGCTGG CGTAATGGTC  
ATCAGTGACC 200  
CGTCTTGA 208

## (26) INFORMATION FOR SEQ ID NO 25:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 212

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Aspergillus clavatus* specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 25:

GACCAGACTC GCTCGCGGGG TTCAGCCGGC ATTCGTGCCG  
GTGTACTTCC 50

CCGTGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG  
GCCTCCGGAA 100

TGTATCACCT CTCGGGGTGT CTTATAGCCG GGGGTGCAAT  
GCGGCCTGCC 150

TGGACCGAGG AACGCGCTTC GGCTCGGACG CTGGCGTAAT  
GGTCGTAAAT 200

GACCCGTCTT GA

212

## (27) INFORMATION FOR SEQ ID NO 26:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 212

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Aspergillus flavus* specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 26:

GACCAGACTC GCCTCCAGGG TTCAGCCGGC ATTCGTGCCG

GTGTACTTCC 50

CTGGGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG

GCTCCCGGAA 100

TGTA GTGCCC TYCGGGGCAC CTTATAGCCG GGAGTGCAAT

GCGGCCAGCC 150

TGGACCGAGG AACGCGCTTC GGCACGGACG CTGGCATAAT

GGTCGYAAAC 200

GACCCGTCTT GA

212

(28) INFORMATION FOR SEQ ID NO 27:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 212

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Aspergillus fumigatus* specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

5

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 27:

GACCAGACTC GCCCGCGGGG TTCAGCCGGC ATTCGTGCCG  
10 GTGTACTTCC 50  
CCGTGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG  
GCCCTCGGAA 100  
TGTATCACCT CTCGGGGTGT CTTATAGCCG AGGGTGCAAT  
GCGGCCTGCC 150  
15 TGGACCGAGG AACGCGCTTC GGCTCGGACG CTGGCGTAAT  
GGTCGTAAAT 200  
GACCCGTCTT GA 212

20 (29) INFORMATION FOR SEQ ID NO 28:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 212

(B) TYPE: nucleic acid

25

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Aspergillus glaucus specific region of 28S gene.

(iii) HYPOTHETICAL: No

30

(iv) ANTISENSE: No

## (v) SEQUENCE DESCRIPTION: SEQ ID NO: 28:

GACCAGACTC GCTTCCGGGG TTCAGCCGGC TTTCGGGCCG  
5 GTGTACTTCC 50  
CCGGGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG  
GCCCCTGGAA 100  
TGTAACGCCT CTCGGGGCGC CTTATAGCCA GGGGTGTCAT  
GCGGCCAGCC 150  
10 TGGACCGAGG AACGCGCTTC GGCACGGACG CTGGCATAAT  
GGTCGTAAAC 200  
GACCCGTCTT GA 212

## 15 (30) INFORMATION FOR SEQ ID NO 29:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 213  
(B) TYPE: nucleic acid  
20 (C) STRANDEDNESS: single  
(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Aspergillus nidulans specific region of 28S gene.

25 (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 29:

30



GACCAGACTC GGCCCCGGGG TTCARCCAGC ACTCGTGCTG  
GTGTACTTCC 50

CCGGGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG  
GCCCCAGGAA 100

5 TGTATCGCCC TCCGGGGTTG TCTTATAGCC TGGGGTGCAA  
TGC GGCCAGC 150

CCG GACCGAG GAACGCGCTT CGGCACGGAC GCTGGCGTAA  
TGGTCGCAA 200

CGACCCGTCT TGA 213

10

(31) INFORMATION FOR SEQ ID NO 30:

(i) SEQUENCE CHARACTERISTICS:

15

(A) LENGTH: 212

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

20

(ii) MOLECULE TYPE: Aspergillus niger specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

25

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 30:

GACCAGACTC GCCCGCGGGG TTCAGCCGGC ATTCGTGCCG  
GTGTACTTCC 50

30

CCGTGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG  
GCCCCTGGAA 100

TGTAGTRCCC TCCGGGGYAC CTTATAGCCA GGGGTGCAAT  
GCGGCCAGCC 150

TGGACCGAGG AACGCGCTTC GGCACGGACG CTGGCATAAT  
GGTCGTAAAC 200

5 GACCCGTCTT GA 212

(32) INFORMATION FOR SEQ ID NO 31:

10 (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 212

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

15

(ii) MOLECULE TYPE: Aspergillus ochraceus specific region of 28S gene.

(iii) HYPOTHETICAL: No

20 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 31:

GACCAGACTC GCCCGCGGGG TTCAGCCGGC ATTCGTGCCG

25 GTGTACTTCC 50

CCGCGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG  
GCCCCCGGAA 100

TGTAGCACCC TTCGGGGTGC CTTATAGCCG GGGGTGCAAT  
GCGGCCAGCC 150

30 TGGACCGAGG AACGCGCTTC GGCACGGACG CTGGCATAAT  
GGTCGTAAAC 200

GACCCGTCTT GA 212

## (33) INFORMATION FOR SEQ ID NO 32:

## 5 (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 212

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

10

(ii) MOLECULE TYPE: Aspergillus terreus specific region of 28S gene.

(iii) HYPOTHETICAL: No

15

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 32:

AACCAGACTC GCTCGCGGGG TTCAGCCGGG CTTCGGCCCCG  
20 GTGTACTTCC 50  
CCGCGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG  
GCCTCCGGAA 100  
TGTAGCGCCC TTCGGGGCGC CTTATAGCCG GGGGTGCAAT  
GCGGCCAGCC 150  
25 TGGACCGAGG AACGCGCTTC GGCACGGACG CTGGCATAAT  
GGTTGTAAAC 200  
GACCCGTCTT GA 212

## 30 (34) INFORMATION FOR SEQ ID NO 33:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 213
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

5

(ii) MOLECULE TYPE: *Aspergillus unguis* specific region of 28S gene.

(iii) HYPOTHETICAL: No

10

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 33:

15

GACCAGACTC GGCCTCGGGG TTCAGCCAGC ACTCGTGCTG  
GTGTACTTCC 50

CCGGGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG  
GCCCCAGGAA 100

20

TGTATCACCC TCCGGGGTTG TCTTATAGCC TGGGGTGCAA  
TGCGGCCAGC 50

CTGGACCGAG GAACGCGCTT CGGCACGGAC GCTGGCATAA  
TGTTTGCAAA 200

CGACCCGTCT TGA 213

25

(35) INFORMATION FOR SEQ ID NO 34:

(i) SEQUENCE CHARACTERISTICS:

30

- (A) LENGTH: 212
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Aspergillus ustus* specific region of 28S gene.

(iii) HYPOTHETICAL: No

5 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 34:

10 GACCAGACTC GGCCCCGGGG TTCAGCCAGC ACTCGTGCTG  
GTGTACTTCC 50  
CCGGGGGCGG GCCAGCGTCG GTTTGGGCGG CCGGTCAAAG  
GCCCCAGGAA 100  
TGTGTCGCCC TCCGGGGCGT CTTATAGCCT GGGGTGCAAT  
15 GCGGCCAGCC 150  
CGGACCGAGG AACGCGCTTC GGCACGGACG CTGGCGTAAT  
GGTCGCAAAC 200  
GACCCGTCTT GA 212

20

(36) INFORMATION FOR SEQ ID NO 35:

(i) SEQUENCE CHARACTERISTICS:

25 (A) LENGTH: 208  
(B) TYPE: nucleic acid  
(C) STRANDEDNESS: single  
(D) TOPOLOGY: linear

30 (ii) MOLECULE TYPE: *Beauveria* species specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 35:

5

GACCAGACTT GGGCTTGGTT GATCATCCGG GGTTCTCCCC  
GGTGCACTCT 50

TCCGGCCCAG GCCAGCATCA GTTCGCCCTG GGGGACAAAG  
10 GCTTCGGGAA 100

CGTGGCTCTC TCCGGGGAGT GTTATAGCCC GTTGCGTAAT  
ACCCTGTGGC 150

GGACTGAGGT TCGCGCATTC GCAAGGATGC TGGCGTAATG  
GTCATCAGTG 200

15 ACCCGTCT 208

(37) INFORMATION FOR SEQ ID NO 36:

20 (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 213

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

25

(ii) MOLECULE TYPE: Bipolaris species specific region of 28S gene.

(iii) HYPOTHETICAL: No

30 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 36:

AGCCAGACTT GCTTGCAGTT GCTCATCCGG GCTTTTGCCC  
GGTGC ACTCT 50  
5 TCTGCAGGCA GGCCAGCATC AGTTTGGGCG GTGGGATAAA  
GGTCTCTGTC 100  
ACGTACCTTC CTTCGGGTTG GCCATATAGG GGAGACGTCA  
TACCACCAGC 150  
CTGGACTGAG GTCCGCGCAT CTGCTAGGAT GCTGGCGTAA  
10 TGGCTGTAAG 200  
CGGCCCCGTCT TGA 213

(38) INFORMATION FOR SEQ ID NO 37:

15

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 105

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

20

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Blastoschizomyces species specific region of 28S  
gene.

25

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 37:

30

TGAAATTGTT GAAAGGGAAG GCGATGGTAG GAATAAGAGG  
CTGCGGTTTG 50

AAATAATTGT TTTTCGGGCC ACGGTCTCCT GAGCCTGCTT  
TCGCACCCGT 100

5 CTTGA 105

(39) INFORMATION FOR SEQ ID NO 38:

10 (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 214

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

15

(ii) MOLECULE TYPE: Blastomyces dermatitidis specific region of 28S  
gene.

(iii) HYPOTHETICAL: No

20

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 38:

25

GACCAGAGTC GGCCGTGGGG GTTCAGCGGG CATTCGTTGC  
CCGTGCACTC 50

CCCCACGGGC GGGCCAGCGT CGGTTTCGAC GGCCGGTCAA  
AGGCCCCCGG 100

30

AATGTGTCGC CTCTCGGGGC GTCTTATAGC CGGGGGTGCA  
ATGCGGCCAG 150



TCGGGACCGA GGAACGCGCT TCGGCACGGA CGCTGGCTTA  
 ATGGTCGTAA 200  
 GCGACCCGTC TTGA 214

5

(40) INFORMATION FOR SEQ ID NO 39:

(i) SEQUENCE CHARACTERISTICS:

10

- (A) LENGTH: 213
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

15

(ii) MOLECULE TYPE: Chrysosporium species specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

20

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 39:

AACCAGACTT GCGCGCGGCC GATCATCCGG TGTTCTCACC  
 25 GGTGCACTCG 50  
 GCCGTGCTCA GGCCAGCATC GGTTTTGGCG GCTGGATAAA  
 GGCCCTAGGA 100  
 ATGTGGCTCC TCTCGGGGAG TGTTATAGCC TAGGGTGCAA  
 TGCAGCCTGC 150  
 30 TGGGACCGAG GACCGCGCTT CGGCTAGGAT GCTGGCGTAA  
 TGGTTGTAAG 200  
 CGGCCCGTCT TGA 213

## (41) INFORMATION FOR SEQ ID NO 40:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 207

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Cladosporium species specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 40:

AACCAGACTT GCTCGCGGTG TTCCGCCGGT CTTCTGACCG  
GTCTACTCGC 50  
CGCGTTGCAG GCCAGCATCG TCTGGTGCCG CTGGATAAGA  
CTTGAGGAAT 100  
GTAGCTCCCT CGGGAGTGTT ATAGCCTCTT GTGATGCAGC  
GAGCGCCGGG 150  
CGAGGTCCGC GCTTCGGCTA GGATGCTGGC GTAATGGTCG  
TAATCCGCCC 200  
GTCTTGA 207

## (42) INFORMATION FOR SEQ ID NO 41:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 213

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Curvularia species specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 41:

AGCCAGACTT GCTTGCAGTT GCTCATCCGG GCTTTTGCCC  
 GGTGCACTCT 50  
 TCTGCAGGCA GGCCAGCATC AGTTTGGGCG GTGGGATAAA  
 GGTCTCTGAC 100  
 ACGTTCCTTC CTTCGGGTTG GCCATATAGG GGAGACGTCA  
 TACCACCAGC 150  
 CTGGACTGAG GTCCGCGCAT CTGCTAGGAT GCTGGCGTAA  
 TGGCTGTAAG 200  
 CGGCCCGTCT TGA 213

(43) INFORMATION FOR SEQ ID NO 42:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 213

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Candida albicans specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 42:

GATCAGACTT GGTATTTTGC ATGCTGCTCT CTCGGGGGCG  
GCCGCTGCGG 50  
TTTACCGGGC CAGCATCGGT TTGGAGCGGC AGGATAATGG  
CGGAGGAATG 100  
TGGCACGGCT TCTGCTGTGT GTTATAGCCT CTGACGATGC  
TGCCAGCCTA 150  
GACCGAGGAC TGCGGTTTTT AACCTAGGAT GTTGGCATAA  
TGATCTTAAG 200  
TCGCCCCGTCT TGA 213

(44) INFORMATION FOR SEQ ID NO 43:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 223

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Candida glabrata specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

5 (v) SEQUENCE DESCRIPTION: SEQ ID NO: 43:

GATCAGACAT GGTGTTTTGC GCCCCTTGCC TCTCGTGGGC  
TTGGGACTCT 50  
10 CGCAGCTCAC TGGGCCAGCA TCGGTTTTGG CGGCCGAAA  
AAACCTAGGG 100  
AATGTGGCTC TGC GCCTCGG TGTAGAGTGT TATAGCCCTG  
GGGAATACGG 150  
CCAGCCGGA CCGAGGACTG CGATACTTGT TATCTAGGAT  
15 GCTGGCATAA 200  
TGGTTATATG CCGCCCGTCT TGA 223

(45) INFORMATION FOR SEQ ID NO 44:

20

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 212

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

25

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Candida guilliermondii* specific region of 28S gene.

(iii) HYPOTHETICAL: No

30

(iv) ANTISENSE: No

## (v) SEQUENCE DESCRIPTION: SEQ ID NO: 44:

GATCAGACTC GATATTTTGT GAGCCTTGCC TTCGTGGCGG  
5 GGTGACCCGC 50  
AGCTTATCGG GCCAGCATCG GTTTGGGCGG TAGGATAATG  
GCGTAGGAAT 100  
GTGACTTTRC TTCGGTGAAG TGTTATAGCC TCGGTTGATG  
CTGCCTGCCT 150  
10 AGACCGAGGA CTGCGATTTT ATCAAGGATG CTGGCATAAT  
GATCCCAAAC 200  
CGCCCGTCTT GA 212

## 15 (46) INFORMATION FOR SEQ ID NO 45:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 214  
(B) TYPE: nucleic acid  
20 (C) STRANDEDNESS: single  
(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Coccidioides immitis* specific region of 28S gene.

25 (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 45:

30

AACCAGACTC GGTCGTGGGG GCTCAGCGGG CATGAGTGCC  
CGTGTACTCC 50

CCCATGCTCC GGGCCAGCAT CAGTTCTGGC GGTTGGTTAA  
AGGCCTCTGG 100

5 AATGTATCGT CCTCCGGGAC GTCTTATAGC CAGGGGCGCA  
ATGCGGCCAG 150

CCGGGACTGA GGAACGCGCT TCGGCACGGA TGCTGGCATA  
ATGGTTGTAA 200

GCGGCCCGTC TTGA 214

10

(47) INFORMATION FOR SEQ ID NO 46:

(i) SEQUENCE CHARACTERISTICS:

- 15 (A) LENGTH: 187  
(B) TYPE: nucleic acid  
(C) STRANDEDNESS: single  
(D) TOPOLOGY: linear

20 (ii) MOLECULE TYPE: Candida kefyr specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

25

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 46:

GATCAGACAT GGC GTTTGCT TCGGCTTTCG CTGGGCCAGC  
30 ATCAGTTTAA 50

GCGGTTGGAT AAATCCTCGG GAATGTGGCT CTGCTTCGGT  
AGAGTGTTAT 100

AGCCCGTGGG AATACAGCCA GCTGGGACTG AGGATTGCGA  
CTTTTGTCAA 150

GGATGCTGGC GTAATGGTTA AATGCCGCCC GTCTTGA 187

5

(48) INFORMATION FOR SEQ ID NO 47:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 213

10

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Candida krusei specific region of 28S gene.

15

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

20

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 47:

CGCCCGACAT GGGGATTGCG CACCGCTGCC TCTCGTGGGC  
GGCGCTCTGG 50

25

GCTTTCCTG GGCCAGCATC GGTTCCTGCT GCAGGAGAAG  
GGGTTCTGGA 100

ACGTGGCTCT TCGGAGTGTT ATAGCCAGGG CCAGATGCTG  
CGTGCGGGGA 150

30

CCGAGGACTG CGGCCGTGTA GGTCACGGAT GCTGGCAGAA  
CGGCGCAACA 200

CCGCCCCGTCT TGA

213



## (49) INFORMATION FOR SEQ ID NO 48:

## (i) SEQUENCE CHARACTERISTICS:

5

(A) LENGTH: 236

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

10

(ii) MOLECULE TYPE: *Cryptococcus laurentii* specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

15

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 48:

AGTCAGTCGT GTCTGGGAGG CTCAGCCGGT TCTGCCGGTG  
20 TATTCCTCTC 50  
AGACGGGTCA ACATCAGTTT TGTCCGACGG ATAATGGCGG  
CGGGAAAGTA 100  
GCACCTCCGG GTGTGTTATA GCCCGCTGTC GCATACGCCG  
GATGAGACTG 150  
25 AGGCATGCAG CTCGCCTTTA TGGCAGGGGT TCGCCCACTT  
TCGAGCTTAG 200  
GATGTTGACG TAATGGCTTT AAACGACCCG TCTTGA 236

30

## (50) INFORMATION FOR SEQ ID NO 49:

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 173
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

5

(ii) MOLECULE TYPE: *Candida lusitanae* specific region of 28S gene.

(iii) HYPOTHETICAL: No

10

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 49:

15

AAGCAGACAC GGTTTTACCG GGCCAGCGTC GAAAAGGGGG  
GAGGAACAAG 50

AACTCGAGAA TGTGGCGCGC ACCTTCGGGY GCGCGTGTTA  
TAGCTCGTGT 100

20

TGACGCCTCC ATCCCTTTTC GAGGCCTGCG ATTCTAGGAC  
GCTGGCGTAA 150

TGGTTGCAAG CCGCCCGTCT TGA 173

(51) INFORMATION FOR SEQ ID NO 50:

25

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 238
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

30

(ii) MOLECULE TYPE: *Cryptococcus neoformans* var *gattii* (serotype B)

specific region of 28S gene.

(iii) HYPOTHETICAL: No

5 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 50:

AGTCAGTCGT GTCTATTGGG TTCAGCCAGC TCTGCTGGTG  
 10 TATTCCCTTT 50  
 AGACGGGTCA ACATCAGTTC TGATCGGTGG ATAAGGGCTG  
 GAGGAATGTG 100  
 GCACTCTTCG GGGTGTGTTA TAGCCTCCTG TCGCATACAC  
 TGGTTGGGAC 150  
 15 TGAGGAATGC AGCTCGCCTT TATGGCCGGG GTTCGCCCAC  
 GTTCGAGCTT 200  
 AGGATGTTGA CAAAATGGCT TTAAACGACC CGTCTTGA  
 238

20 (52) INFORMATION FOR SEQ ID NO 51:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 238

(B) TYPE: nucleic acid

25 (C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Cryptococcus neoformans (serotype A) specific  
 region  
 30 of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 51:

5

AGTCAGTCGT GTCTATTGGG TTCAGCCAGT TCTGCTGGTG  
TATTCCTTT 50

AGACGGGTCA ACATCAGTTC TGATCGGTGG ATAAGGGCTG  
10 GGGGAATGTA 100

GCACTCTTCG GAGTGTGTTA TAGCCTCCTG TCGCATACAC  
TGGTTGGGAC 150

TGAGGAATGC AGCTCGCCTT TATGGCCGGG GTTCGCCCAC  
GTTCGAGCTT 200

15 AGGATGTTGA CAAAATGGCT TTAAACGACC CGTCTTGA  
238

(53) INFORMATION FOR SEQ ID NO 52:

20

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 211

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

25

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Candida parapsilosis specific region of 28S gene.

(iii) HYPOTHETICAL: No

30

(iv) ANTISENSE: No

## (v) SEQUENCE DESCRIPTION: SEQ ID NO: 52:

GATCAGACTT GGTATTTTGT ATGTTACTCT CTCGGGGGTG  
 5 GCCTCTACAG 50  
 TTTACCGGGC CAGCATCAGT TTGAGCGGTA GGATAAGTGC  
 AAAGAAATGT 100  
 GGCAGTGCTT CGGTAGTGTG TTATAGTCTT TGTCGATACT  
 GCCAGCTTAG 150  
 10 ACTGAGGACT GCGGCTTCGG CCTAGGATGT TGGCATAATG  
 ATCTTAAGTC 200  
 GCCCGTCTTG A 211

## (54) INFORMATION FOR SEQ ID NO 53:

15

## (i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 238
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

20

(ii) MOLECULE TYPE: *Cryptococcus terreus* specific region of 28S gene.

(iii) HYPOTHETICAL: No

25

(iv) ANTISENSE: No

## (v) SEQUENCE DESCRIPTION: SEQ ID NO: 53:

30

AGTCAGTCAT GTCTATTGGA CTCAGCCGGT TCTGCCGGTG  
 TACTTCCTTT 50

AGATGGGGTC AACATCAGTT TTGATCGCTG GAAAAGGGCA  
GGAGGAATGT 100

AGCACTCTCG GGTGAACTTA TAGCCTTCTG TCGTATACAG  
TGGTTGGGAC 150

5 TGAGGAACGC AGCATGCCTT TATGGCCGGG GTTCGCCCCAC  
GTACATGCTT 200

AGGATGTTGA CATAATGGCT TTAAACGACC CGTCTTGA  
238

10 (55) INFORMATION FOR SEQ ID NO 54:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 211

(B) TYPE: nucleic acid

15 (C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Candida tropicalis specific region of 28S gene.

20 (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 54:

25 GATCAGACTT GGTATTTTGT ATGTTACTTC TTCGGGGGTG  
GCCTCTACAG 50

TTTATCGGGC CAGCATCAGT TTGGGCGGTA GGAGAATTGC  
GTTGGAATGT 100

30 GGCACGGCTT CGGTTGTGTG TTATAGCCTT CGTCGATACT  
GCCAGCCTAG 150

ACTGAGGACT GCGGTTTATA CCTAGGATGT TGGCATAATG  
ATCTTAAGTC 200  
GCCCCGTCTTG A 211

5 (56) INFORMATION FOR SEQ ID NO 55:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 211  
(B) TYPE: nucleic acid  
10 (C) STRANDEDNESS: single  
(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Fusarium species specific region of 28S gene.

15 (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 55:

20

GACCAGACTT GGGCTTGGTT AATCATCTGG GGTTCTCYCC  
AGTGCACTTT 50

TCCAGTCCAG GCCAGCATCA GTTTTCSCCG GGGGATAAAG  
25 RCTTCGGGAA 100

TGTGGCTCYC YYCGGGGAGT GTTATAGCCC GTTGYGTAAT  
ACCCTGGBGG 150

GGACTGAGGT TCGCGCWTCT GCAAGGATGC TGGCGTAATG  
GTCATCAACG 200

30 ACCCGTCTTG A 211

## (57) INFORMATION FOR SEQ ID NO 56:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 238

5 (B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

10 (ii) MOLECULE TYPE: Filobasidium capsuligenum specific region of 28S gene.

(iii) HYPOTHETICAL: No

15 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 56:

AGTCAGTCAT GTCTATTGGA CTCAGCCGGT TCTGCCGGTG  
20 TATTTTCCTTT 50  
AGATGGGGTC AACATCAGTT TTGACCGTTG GATAAAGGCA  
GGAAGAATGT 100  
AGCACTCTCG GGTGAACTTA TAGCTTCTTG TCACATACAA  
TGGTTGGGAC 150  
25 TGAGGAACGC AGCATGCCTT TATGGCCGGG ATTCGTCCAC  
GTACATGCTT 200  
AGGATGTTGA CATAATGGCT TTAAACGACC CGTCTTGA  
238

30

## (58) INFORMATION FOR SEQ ID NO 57:



## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 238

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Filobasidiella neoformans* var *bacillispora* (serotype

C) specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 57:

AGTCAGTCGT GTCTATTGGG TTCAGCCAGC TCTGCTGGTG  
TATTCCCTTT 50

AGACGGGTCA ACATCAGTTC TGATCGGTGG ATAAGGGCTG  
GAGGAATGTG 100

GCACTCTTCG GGGTGTGTTA TAGCCTCCTG TCGCATAAC  
TGGTTGGGAC 150

TGAGGAATGC AGCTCGCCTT TATGGCCGGG GTTCGCCCAC  
GTTCGAGCTT 200

AGGATGTTGA CAAAATGGCT TTAAACGACC CGTCTTGA  
238

(59) INFORMATION FOR SEQ ID NO 58:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 238

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

5 (ii) MOLECULE TYPE: *Filobasidiella neoformans* var *neoformans* (serotype  
D) specific region of 28S gene.

(iii) HYPOTHETICAL: No

10 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 58:

15 AGTCAGTCGT GTCTATTGGG TTCAGCCAGT TCTGCTGGTG  
TATTCCCTTT 50

AGACGGGTCA ACATCAGTTC TGATCGGTGG ATAAGGGCTG  
GAGGAATGTG 100

GCACTCTTCG GGGTGTGTTA TAGCCTCCTG TCGCATAACAC  
20 TGGTTGGGAC 150

TGAGGAATGC AGCTCGCCTT TATGGCCGGG GTTCGCCCAC  
GTTCGAGCTT 200

AGGATGTTGA CAAAATGGCT TTAAACGACC CGTCTTGA  
238

25

(60) INFORMATION FOR SEQ ID NO 59:

(i) SEQUENCE CHARACTERISTICS:

30 (A) LENGTH: 236

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

## (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Filobasidium uniguttulatum* specific region of 28S gene.

5

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

10

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 59:

AGTCAGTCGT GCTCAATGGA CTCAGCCGTT CTGCGGTGTA  
TTTCCATTGG 50

15

GTGGGGTCAA CATCAGTTTT GATCGCTGGA TAAAGGCAGG  
AGGAATGTAG 100

CACCCCCGGG TGAAGTTATA GCCTCTTGTC ACATACAGTG  
GTTGGGACTG 150

20

AGGAACGCAG CATGCCTTTA TGGCCGGGAT TCGTCCACGT  
ACATGCTTAG 200

GATGTTGACA TAATGGCTTT AAACGACCCG TCTTGA 236

(61) INFORMATION FOR SEQ ID NO 60:

25

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 204

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

30

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Geotrichum* species specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 60:

5  
10 AATCAGACTT GGTGCTGTTG TTCAACTRTG TTTCGGCATA  
GTGTACTCAG 50  
CAGTACTAGG CCAAGGTGGG GTGTTTGGGA GTGAAAAAGA  
AGTAGGAACG 100  
TAACTCTTCG GAGTGTTATA GCCTACTTTC ATAGCTCCTC  
AGGCGCCTCA 150  
15 GGACTGCGCT TCGGCAAGGA CCTTGGCATA ATGATTCTAT  
ACCGCCCGTC 200  
TTGA 204

20

(62) INFORMATION FOR SEQ ID NO 61:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 214

25

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Histoplasma capsulatum specific region of 28S gene.

30

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 61:

5

GAYCAGAGTC GGCCGYGGGG GTTCAGCGGG CATTCGTTGC  
CCGTGCAATC 50

CCCCGCGGCC GGGCCAGCGT CGGTTTCGAC GGCCGGTCAA  
AGGCCCCCGG 100

10

AATGTGTCGC CTCTCGGGGC GTCTTAGC CGGGGGTGCA  
ATGCGGCCAG 150

TCGGGACCGA GGAACGCGCT CCGGCACGGA CGCTGGCTTA  
ATGGTCGTCA 200

GCGACCCGTC TTGA

214

15

(63) INFORMATION FOR SEQ ID NO 62:

(i) SEQUENCE CHARACTERISTICS:

20

(A) LENGTH: 215

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

25

(ii) MOLECULE TYPE: Malbranchea species specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

30

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 62:

AGACAGACTC GAGCGCGGGG GCTCAGCGGG TATTGTTATG  
CCCGTGCACT 50

CCCCCGCGCC CGGGCCAGCA TCAGTTTTGG CGGCCGGTCA  
AAGGCCCTTG 100

5 GAATGTATCG TCCTCCGGGA CGTCTTATAG CCAAGGGTGC  
AATGCGGCCA 150

GCCGGGACTG AGGAACGCGC TTCGGCACGG ATGCTGGCGT  
AATGGCTGTA 200

AGCGGCCCGT CTTGA 215

10

(64) INFORMATION FOR SEQ ID NO 63:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 237

15

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Mucor species specific region of 28S gene.

20

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

25

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 63:

AGCCAGACTG GTTTGACTGT AATCAACCTA GAATTCGTTC  
TGGGTGCACT 50

TGCAGTCTAT ACCTGCCAAC AACAGTTTGA TTTGGAGGAA  
30 AAAATTAGTA 100

GGAATGTAGC CTCTCGAGGT GTTATAGCCT ACTATCATAC  
TCTGGATTGG 150

ACTGAGGAAC GCAGCGAATG CCWTTAGGCR AGATTGCTGG  
GTGCTTTCGC 200

TAATAAATGT TAGAATTCT GCTTCGGGTG GTGCTAA 237

5

(65) INFORMATION FOR SEQ ID NO 64:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 209

10

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Paecilomyces species specific region of 28S gene.

15

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

20

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 64:

GACCAGACTT GGGCCCGGTG GATCATCCAG CGTTCTCGCT  
GGTGCACTCC 50

25

GCCGGGTTCA GGCCAGCATC AGTTCGCCGC GGGGGAAAAA  
GGCTTCGGGA 100

ACGTGGCTCC TACGGGAGTG TTATAGCCCG TTGCATAATA  
CCCTGGGGCG 150

GACTGAGGTT CGCGCTCCGC AAGGATGCTG GCGTAATGGT  
CATCAGCGAC 200

30

CCGTCTTGA

209

## (66) INFORMATION FOR SEQ ID NO 65:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 199

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Penicillium species specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 65:

GACCAGACTC GCCCACGGGG TTCAGCCGGC ATTCGTGCCG  
GTGTACTTCC 50

CCGCGGGCGG GCCAGCGTCG GTTTGGKCGG CCGGTCAAAG  
GCCCTCGGAA 100

TRTAACGCCC CCCGGGGCGT CTTATAGCCG AGGGTGCCAT  
GCGGCCAGCM 150

CAGACCGAGG AACGCGCTTC GGCTCGGACG CTGGCATAAT  
GGTCGTAAA 199

## (67) INFORMATION FOR SEQ ID NO 66:

## (i) SEQUENCE CHARACTERISTICS:



- (A) LENGTH: 210
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

5

(ii) MOLECULE TYPE: *Pseudallescheria boydii* region of 28S gene.

(iii) HYPOTHETICAL: No

10

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 66:

15

GACCAGACTT GTGCCCGTCG AATCAGCCGC CGCTCGTCGG  
CGGCGCACTT 50

CGGCGGGCTC AGGCCAGCAT CAGTTCGCTG CAGGGGGAGA  
AAGGCGATGG 100

20

GAATGTGGCT CTTCGGAGTG TTATAGCCCG CCGCGCAATA  
CCCCTCGGCG 150

GACTGAGGAC CGCGCATCTG CAAGGATGCT GGC GTAATGG  
TCGTCAGCGA 200

CCCGTCTTGA

210

25

(68) INFORMATION FOR SEQ ID NO 67:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 244
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

30

(ii) MOLECULE TYPE: Rhizopus species (NO: 1) specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 67:

AGCCAGACTG GCTTGTCTGT AATCAATCTA GGTTCGTGC  
CTGGATGCAC 50  
TTGCAGACTA TTTGCCTGCC AACGACAATT TTTTTGAGT  
GTAAAAACTA 100  
TTGGAAATGT GGCCAATATT TATTTATTGG TGTTATAGTC  
CTTTAGAAAA 150  
TACCTTGAAT TGGATTGAGG AACGCAGCGA ATGCTTCTCT  
TTNGAGGCAA 200  
AGTCTTTTAT TGGGATTAC GGATCAGACT GTGGCATTGT CACA  
244

(69) INFORMATION FOR SEQ ID NO 68:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 215
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: single
- (D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Rhizopus species (NO: 2) specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

5 (v) SEQUENCE DESCRIPTION: SEQ ID NO: 68:

AGCCAGACTG GCTTGTCTGT AATCAATCTA GGCTTCGGCC  
TGGATGCACT 50

TGCAGGCTAT GCCTGCCAAC GACAATTGA CTTGAGGGAA  
10 AAAACTAGGG 100

GAAATGTGGC CCACTTGTGG GTGTTATAGT CCCTTAGAAA  
ATACCTTGGG 150

TTGGATTGAG GAACGCAGCG AATGCTTATT GGCGAGTTTT  
CCAGGAAGGT 200

15 TTTCTGAGGT ACTAC 215

(70) INFORMATION FOR SEQ ID NO 69:

(i) SEQUENCE CHARACTERISTICS:

- 20 (A) LENGTH: 215  
(B) TYPE: nucleic acid  
(C) STRANDEDNESS: single  
(D) TOPOLOGY: linear

25 (ii) MOLECULE TYPE: Rhizopus species (NO: 3) specific region of 28S  
gene.

(iii) HYPOTHETICAL: No

30 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 69:

AGCCAGACTG GCTTGTCTGT AATCAGTCTA AGCTTCGGCT  
 TGGATGCACT 50  
 5 TGCAGGCTAT GCCTGCCAAC GACAATTG GCTTGAGGGAA  
 AAAACTAAGG 100  
 GAAATGTGGC CCATCCGTGG GTGTTATAGT CCCTTAGAAA  
 ATACCTTGGG 150  
 CTGGATTGAG GTACGCAGCG AATGCTATTT GCGGAGTTGG  
 10 CTGGGAATAT 200  
 TTTCTGAGGT GCTTT 215

(71) INFORMATION FOR SEQ ID NO 70:

15

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 210

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

20

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Sporothrix species specific region of 28S gene.

(iii) HYPOTHETICAL: No

25

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 70:

30

GACCAGACTT GCGCCYCGCG GACCACCCGG CGTTCTCGCC  
 GGTGCACTCT 50

GCGKKGCGCA GGCCAGCATC GGTTCCTCCCA GGGGGACAAA  
GGCCGCGGGA 100

ACGTAGCTCC TTCGGGAGTG TTATAGCCCG CGGCGGCATG  
CCCCTGGGGG 150

5 GACCGAGGAC CGCGCTTCGG CAAGGATGCT GGC GTAATGG  
TCACCAGCGA 200  
ACCGTCTTGA 210

10 (72) INFORMATION FOR SEQ ID NO 71:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 208  
(B) TYPE: nucleic acid  
15 (C) STRANDEDNESS: single  
(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Scopulariopsis brevicaulis* specific region of 28S  
gene.

- 20 (iii) HYPOTHETICAL: No  
(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 71:

25 GACCAGACTT GCGCCCGTCG GATCAACCGT CGCTTGCGGC  
GGCGCACTCC 50  
GGCGGGCTCA GGCCAGCATC AGTTCGTCCG GGGGGAGAAA  
GGCGGCGGGA 100  
30 ATGTGGCTCT TCGGAGTGTT ATAGCCCGCC GTGTAATACC  
CTCGGGTGGA 150

CTGAGGACCG CGCGTATGCA AGGATGCTGG CGTAATGGTC  
GTCAGCGACC 200  
CGTCTTGA 208

5 (73) INFORMATION FOR SEQ ID NO 72:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 210

(B) TYPE: nucleic acid

10 (C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: *Scopulariopsis brumptii* specific region of 28S gene.

15 (iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 72:

20

GACCAGACTC GCGCCCGTCG GATCAGCCGT CGCTCGTCGG  
CGGCGCACTC 50

CGGCGGGCTC GGGCCAGCAT CAGTTCGCCT CGGGGGGAGA  
25 AAGGCGGCGG 100

GAATGTGGCT CTACGGAGTG TTATAGCCCG CCGCGTAATA  
CCCCCGGGCG 150

GACTGAGGAC CGCGCGTATG CAAGGATGCT GGC GTAATGG  
TCGTCAGCGA 200

30

CCCGTCTTGA

210

## (74) INFORMATION FOR SEQ ID NO 73:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 214

5 (B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

10 (ii) MOLECULE TYPE: *Saccharomyces cerevisiae* specific region of 28S gene.

(iii) HYPOTHETICAL: No

15 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 73:

20 GATCAGACAT GGTGTTTTGT GCCCTCTGCT CCTTGTGGGT  
AGGGGAATCT 50CGCATTTCAC TGGGCCAGCA TCAGTTTTGG TGGCAGGATA  
AATCCATAGG 100AATGTAGCTT GCCTCGGTAA GTATTATAGC CTGTGGGAAT  
25 ACTGCCAGCT 150GGGACTGAGG ACTGCGACGT AAGTCAAGGA TGCTGGCATA  
ATGGTTATAT 200

GCCGCCCCGTC TTGA

214

30

## (75) INFORMATION FOR SEQ ID NO 74:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 236

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Trichosporon beigeli specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 74:

AGTCAGTCGT GTTCTTTGGA TTCAGCCAGT TCTGCTGGTC  
TACTTCCTTG 50

GAACGGGTCA ACATCAGTTT TGTCCGGTGG ATAAAGGTAG  
TAGGAATGTG 100

ACTTCTCCGG AAGTGTTATA GCCTATTATC ACATACACTG  
GGTGAGACTG 150

AGGACTGCAG CTCGCCTTTA TGGCCGGCCT TCGGGCACGT  
TCGAGCTTAG 200

GATGTTGACA TAATGGCTTT AAACGACCCG TCTTGA 236

(76) INFORMATION FOR SEQ ID NO:75:

## (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear



(ii) MOLECULE TYPE: Nucleic acid probe for *Paracoccidioides brasiliensis*

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO:75:

ACTCCCCCGT GGTC

14

(77) INFORMATION FOR SEQ ID NO:76:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 14

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic acid probe for *Pneumocystis carinii*

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO:76:

GGAAGGGAAA TTGG

14

(78) INFORMATION FOR SEQ ID NO:77:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 223

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Paracoccidioides brasiliensis specific region of 28S

5 gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

10

(v) SEQUENCE DESCRIPTION: SEQ ID NO:77:

GCGCTTGCGA CCAGAGTCGG CCGCGGGGGC TCAGCGGGCA  
15 CTCGTTGCCC 50  
GTGCACTCCC CCGTGGTCGG GCCAGCGTCG GTTTCGACGG  
CCGGTCAAAG 100  
GCCCCCGGAA TGTGTCGCCT CTCGGGGCGT CTTATAGCCG  
GGGGTGCAAT 150  
20 GCGGCCAGTC GGGACCGAGG AACGCGCTCC GGCACGGACG  
CTGGCTTAAT 200  
GGTCGTAAGC GACCCGTCTT GAA 223

(79) INFORMATION FOR SEQ ID NO:78:

25

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 257

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

30

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Pneumocystis carinii specific region of 28S gene.

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO:78:

5  
10 GCGCTTGTGA TCAGACATGC CTTTATAGGA GATGCCATTG  
TTTCGGCATT 50  
GGCAGTATTA TCCGAATTGG CAGGCCAGCA TCGGTTTCAG  
TTACTGGATA 100  
AAACTGGAAG AAGGTAGGCT CTCTTCGGAG GGTTTTTTAG  
CTTCCAGTAG 150  
15 CTGCAGTGAC CGGGACCGGA AGGGAAATTG GGTCTTTGAA  
GACCTTATGA 200  
TGTTGGCAGA AATGGTCCTA AGCGACCCGT CTTGAAACAC  
GGACCAAGGA 250  
GTCTAAT 257

20

(80) INFORMATION FOR SEQ ID NO 79:

(i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 20

25

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic Acid Primer

30

(iii) HYPOTHETICAL: No

(iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 79:

5 ATCAATAAGC GGAGGAAAAG 20

(81) INFORMATION FOR SEQ ID NO 80:

10 (i) SEQUENCE CHARACTERISTICS:

(A) LENGTH: 20

(B) TYPE: nucleic acid

(C) STRANDEDNESS: single

(D) TOPOLOGY: linear

(ii) MOLECULE TYPE: Nucleic Acid Primer

(iii) HYPOTHETICAL: No

20 (iv) ANTISENSE: No

(v) SEQUENCE DESCRIPTION: SEQ ID NO: 80:

CTCTGGCTTC ACCCTATTC 20

CLAIMS

We claim:

1. An oligonucleotide hybridization probe for *Paracoccidioides brasiliensis*, said  
5 probe having the nucleotide residue sequence of (SEQ ID NO:75) or the complement thereof.
2. An oligonucleotide hybridization probe for *Pneumocystis carinii*, said probe having the nucleotide residue sequence of (SEQ ID NO:76) or the complement thereof.
- 10 3. A method of determining whether one or more fungal species selected from the group of fungal species consisting of *Paracoccidioides brasiliensis* and *Pneumocystis carinii* is present in a sample of fungi, said method comprising the following steps:
  - a) extracting nucleic acid material from fungi contained in said sample;
  - 15 b) adding two known oligonucleotide primers, one of said primers being (SEQ ID NO:1) or (SEQ ID NO:2), said primers bracketing a hypervariable region on the 28S rDNA or rRNA present in the fungal species of said group;
  - c) amplifying the sequence between said primers; and
  - 20 d) using one or more labeled probes directed to a portion of the hypervariable region bracketed by said primers, each said labeled probe being specific for one of said fungal species from said group, to determine whether said fungal species identified by each said labeled probe is present in said sample.
- 25 4. The method of claim 3 in which, in said amplifying step, said amplifying procedure is the polymerase chain reaction
5. The method of claim 3 in which said one or more probes is selected from the  
30 group consisting of (SEQ ID NO:75) and (SEQ ID NO:76).

6. The method of claim 3 wherein, in step (d), more than one probe is used, each said probe being connected to (a) a different signal moiety or (b) a moiety which allows separation of said probes.

5 7. A method of determining whether one or more fungal species selected from the group of fungal species consisting of *Paracoccidioides brasiliensis* and *Pneumocystis carinii* is present in a sample of fungi, said method comprising the following steps:

- a) extracting nucleic acid material from fungi contained in said sample;
- b) adding a universal fungal probe selected from the group consisting of  
10 (SEQ ID NO: 1), (SEQ ID NO:2) and the complements thereof;
- c) using one or more second probes, each said second probe being specific for one of said fungal species from said group, wherein said one or more second probes each has a nucleotide residue sequence selected from the group consisting of (SEQ ID NO:75), (SEQ ID NO:76) and the  
15 complements thereof; and
- d) determining whether said fungal species identified by each said second probe is present in said sample,  
wherein at least one of said probes is connected to a signal moiety and at least one of said probes is connected to a moiety that allows separation of  
20 said probes.

8. A method of determining whether one or more fungal species selected from the group of fungal species consisting of *Paracoccidioides brasiliensis* and *Pneumocystis carinii* is present in a sample of fungi, said method comprising the following steps:

- a) extracting nucleic acid material from fungi contained in said sample; and
- b) using one or more labeled probes, each said labeled probe being specific  
25 for one of said fungal species from said group, to determine whether said fungal species identified by each said labeled probe is present in said sample, wherein said one or more labeled probes each has a nucleotide  
30 residue sequence selected from the group consisting of (SEQ ID NO:75), (SEQ ID NO:76) and the complements thereof.

9. A species specific reference oligonucleotide for *Paracoccidioides brasiliensis* having the nucleotide residue sequence of SEQ ID NO:77 or the complement thereof.

10. A species specific reference oligonucleotide for *Pneumocystis carinii* having the nucleotide residue sequence of SEQ ID NO:78 or the complement thereof.

11. A method of determining whether one or more fungal species selected from the group of fungal species consisting of *Paracoccidioides brasiliensis* and *Pneumocystis carinii* is present in a sample of fungi, said method comprising the following steps:

- a) extracting nucleic acid material from fungi contained in said sample;
- b) adding two known oligonucleotide primers, one of said primers being (SEQ ID NO:1) or (SEQ ID NO:2), said primers bracketing a hypervariable region on the 28S rDNA or rRNA present in the fungal species of said group;
- c) amplifying the sequence between said primers; and
- d) using one or more labeled probes directed to a portion of the hypervariable region bracketed by said primers, each said labeled probe being specific for one of said fungal species from said group, wherein each said one or more labeled probes is fully complementary to a species-unique nucleotide sequence in said hypervariable region, to determine whether said fungal species identified by each said labeled probe is present in said sample, wherein, furthermore, said one or more labeled probes each has a nucleotide residue sequence consisting of from 10 to 50 consecutive nucleotide residues from a sequence selected from the group consisting of (SEQ ID NO:77) and (SEQ ID NO:78) and the complements thereof.

12. A method of claim 11 in which said amplifying procedure is the polymerase chain reaction.

13. A method of claim 11 wherein more than one third probe is used, each said third probe connected to a different signal moiety or moiety which allows separation of said third probe.

5 14. A method of determining whether one or more fungal species selected from the group of fungal species consisting of *Paracoccidioides brasiliensis* and *Pneumocystis carinii* is present in a sample of fungi, said method comprising the following steps:

- a) extracting the nucleic acid material from the fungi contained in said sample;
- 10 b) adding a universal fungal probe selected from the group consisting of (SEQ ID NO: 1), (SEQ ID NO:2) and the complements thereof;
- c) using one or more second probes, each said second probe being specific for one of said fungal species from said group, wherein each said one or more second probes is fully complementary to a species-unique  
15 nucleotide sequence in said hypervariable region and wherein, further, said one or more second probes each has a nucleotide residue sequence consisting of from 10 to 50 consecutive nucleotide residues from a sequence selected from the group consisting of (SEQ ID NO:77) and (SEQ ID NO:78) and the complements thereof; and
- 20 d) determining whether said fungal species identified by each said second probe is present in said sample,  
wherein at least one of said probes is connected to a signal moiety and at least one of said probes is connected to a moiety that allows separation of said probes.

25

15. An oligonucleotide probe/primer for fungi, said probe having the nucleotide residue sequence of (SEQ ID NO:79) or the complement thereof.

16. An oligonucleotide probe/primer for fungi, said probe having the nucleotide  
30 residue sequence of (SEQ ID NO:80) or the complement thereof.






FIGURE 1

# INTERNATIONAL SEARCH REPORT

International Application No

PCT/IB 98/00865

## A. CLASSIFICATION OF SUBJECT MATTER

IPC 6 C12Q1/68

According to International Patent Classification (IPC) or to both national classification and IPC

## B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC 6 C12Q

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

## C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category *	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	WO 96 21741 A (CIBA CORNING DIAGNOSTICS CORP ; SANDHU GURPREET S (US); KLINE BRUCE) 18 July 1996	3, 4, 6, 9, 11-14
Y	see the whole document	7, 8
X	SANDHU G S ET AL: "MOLECULAR PROBES FOR DIAGNOSIS OF FUNGAL INFECTIONS" JOURNAL OF CLINICAL MICROBIOLOGY, vol. 33, no. 11, 1 November 1995; pages 2913-2919, XP002002934	15, 16
Y	see the whole document	3, 4, 6-8, 11-14
	--- -/--	

☒ Further documents are listed in the continuation of box C.

☒ Patent family members are listed in annex.

### \* Special categories of cited documents :

"A" document defining the general state of the art which is not considered to be of particular relevance

"E" earlier document but published on or after the international filing date

"L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)

"O" document referring to an oral disclosure, use, exhibition or other means

"P" document published prior to the international filing date but later than the priority date claimed

"T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention

"X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone

"Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.

"&" document member of the same patent family

Date of the actual completion of the international search

1 September 1998

Date of mailing of the international search report

21/09/1998

Name and mailing address of the ISA

European Patent Office, P.B. 5818 Patentlaan 2  
NL - 2280 HV Rijswijk  
Tel. (+31-70) 340-2040, Tx. 31 651 epo nl,  
Fax: (+31-70) 340-3016

Authorized officer

Knehr, M

## INTERNATIONAL SEARCH REPORT

International Application No

PCT/IB 98/00865

## C.(Continuation) DOCUMENTS CONSIDERED TO BE RELEVANT

Category	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
Y	TANG X ET AL.: "A single-tube nested PCR for <i>Pneumocystis carinii</i> f. sp. hominis" JOURNAL OF CLINICAL MICROBIOLOGY, vol. 35, no. 6, 1997, pages 1597-1599, XP002075978 see the whole document ---	3,4,6-8, 11-14
Y	WO 96 29432 A (UNIV BOSTON) 26 September 1996 see the whole document ---	3,4,6-14
Y	LIU Y ET AL.: "Sequence and variability of the 5.8S and 26S rRNA genes of <i>Pneumocystis carinii</i> " NUCLEIC ACIDS RESEARCH, vol. 20, no. 14, 1992, pages 3763-3772, XP002075979 see abstract see page 3763, column 2, paragraph 2 see page 3766, column 1, paragraph 1 - column 2, paragraph 1 see page 3767, column 2, paragraph 2 see page 3769, column 1, paragraph 2 - page 3770, column 2, paragraph 1; figures 1,5,7,8; tables 1,3 ---	3,4,6-14
Y	WO 95 29260 A (CIBA GEIGY AG ;LIGON JAMES M (CH); BECK JAMES J (US)) 2 November 1995 see the whole document ---	3,4,6-8, 11-15
Y	US 5 501 951 A (MILLIMAN CURT L) 26 March 1996 see abstract see column 2, line 57 - column 3, line 2 see column 3, line 29 - line 51 see column 4, line 21 - line 34 see column 6, line 12 - line 25; table 1 ---	3,4,6-8, 11-15
A	WO 94 02636 A (HITACHI CHEMICAL CO LTD ;MITSUHASHI MASATO (US); COOPER ALLAN (US)) 3 February 1994 see the whole document ---	
P,X	SANDHU G S ET AL.: "Molecular detection and identification of <i>Paracoccidioides brasiliensis</i> " JOURNAL OF CLINICAL MICROBIOLOGY, vol. 35, no. 7, 1997, pages 1894-1896, XP002075980 see the whole document -----	1,3-8, 11-14

# INTERNATIONAL SEARCH REPORT

Information on patent family members

International Application No

PCT/IB 98/00865

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
WO 9621741 A	18-07-1996	US 5763169 A	09-06-1998
		US 5707802 A	13-01-1998
		AU 693625 B	02-07-1998
		AU 4313896 A	31-07-1996
		BR 9607497 A	23-12-1997
		CA 2209247 A	18-07-1996
		EP 0804619 A	05-11-1997
		PL 321139 A	24-11-1997
WO 9629432 A	26-09-1996	AU 5366296 A	08-10-1996
WO 9529260 A	02-11-1995	US 5585238 A	17-12-1996
		AU 2424495 A	16-11-1995
		BR 9507513 A	02-09-1997
		CN 1147276 A	09-04-1997
		EP 0758404 A	19-02-1997
		JP 10501965 T	24-02-1998
		PL 316962 A	17-02-1997
US 5501951 A	26-03-1996	US 5284747 A	08-02-1994
WO 9402636 A	03-02-1994	CA 2140763 A	03-02-1994
		EP 0652971 A	17-05-1995
		JP 7509361 T	19-10-1995
		MX 9300494 A	29-07-1994
		US 5639612 A	17-06-1997
		US 5580971 A	03-12-1996